

Product datasheet for **SC314427**

SPTY2D1 (NM_194285) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPTY2D1 (NM_194285) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPTY2D1
Synonyms:	Spt2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314427 representing NM_194285.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGACTTCAGAGAAATTCATGATAGCTTCCAAGGGACAAGGTGTCAACAATGTGCCGAAAAGGTAT
AGTTTGGCAGTGGGGCCTCCAAAAAAGACCCAAAAAGTTAAAGGTGTCCAATCAGCAGCTGTACAAGCT
TTTCTTAAAAGGAAAGAGAGGAGCTGAGACGAAAAGCCTTAGAGGAGAAAAGGAGAAAAGAGGAACTA
GTGAAAAAGCGAATTGAGCTCAACATGACAAGAAAGCAAGAGCTATGGCCAAGAGGACAAAGGATAAT
TTCCATGGTTACAATGGGATTCTATTGAGGAAAAGTCAAAGAAGAGGCAGGCAACAGAAAGCCATACC
AGCCAAGGAACCGACCGAGAGTATGAAATGGAAGAAGAGAATGAATTCCTCGAGTACAATCACGCAGAG
TCAGAGCAGGAGTATGAGGAAGAGCAAGAACCTCCCAAAGTTGAAAGCAAACCAAGGTCCCCCTTAAA
AGTGCCCCACCACCCATGAATTCCTGATTTACTCAGGCTGGCTGAGAAAAAGCAGTTTGAACCAAGT
GAAATCAAGGTAGTGAAGAAATCAGAAGAGCGACCTATGACCGCAGAAGAACTTAGGGAGCGAGAATTC
CTTGAACGAAAGCATAGGAGAAAAAACTTGAGACAGATGAAAACTACCTCCAAGTGTGTCCAAAAAG
GCACCCTCTCAGAAAGAAAGTGTGGGCACAAAACCTTAGCAAAGGTTCTGGAGACAGGCATCCTTCTTCC
AAAGGAATGCCCTTCTCATGTGAGAAGAAATCAAGACCCAGCATGGCCAATGAGAAACACCTCGCT
TTGTCTTCATCCAAATCCATGCCAGGAGAGAGGATCAAGGCAGGATCTGGCAATAGCTCCCAACCTCA
CTTCGTGAGGGCCACGACAAACCTGTTTTAATGGAGCTGGAAAGCCTCATTCCAGCACCTCTTCACCA
AGTGTCCCAAAGACTTCTGCTAGCAGGACTCAGAAATCTGCTGTTGAGCAGCAAGCCAAAAAATCTCTG
TCCCATCTAGCCATTCCAGGCTGGGCCATGGTCACCCACACAATAAGGCTAAGAGTCCAGGTGTC
AGGCAGCCAGGCAGCAGCTCTAGCTCAGCCCCTGGGCAGCCAGCACAGGGGTTGCTCGACCCACAGTT
AGTTCTGGCCCTGTGCCTAGGCGCCAGAATGGCAGCTCCAGCTCAGGACCTGAGCGATCAATCAGTGGG
TCCAAGAAGCCAAACCAATGACTCAAATCCCTCTAGGCGGACAGTCAGTGGTACATGTGGCCCTGGACAA
CCTGCAAGCAGCTCAGGTGGCCCTGGGCGACCCATCAGTGGTTCAAGTATGTTCTGCAAGACCTTGGGC
AGCTCTCGTGGCCCTGGCCGGCCTGTGAGCAGTCCACATGAACCTCGACGACCAAGTGGCTTGGGC
CCCCGGGGCGGTCTGTGAGTGGCCCTGGGAGATCCATAAGTGGCTCAATTCAGCTGGACGGACTGTC
AGTAATTCAGTCCAGGAAGACCAAGTGAAGCAGCTTGGGACCTGGGCAACAGTTAGTAGCTCAGGTCCC
ACTATAAGCCTAAGTGCAGTGTGTCTCCGAAACAATTTCTTCCAAGAATATCATTAGCCGGTCCAGC
AATGGACAGATGAATGGAATGAAGCCTCCCTATCTGGCTACAGAGCTGCCCAAGGTCTCAAAGGCTT
CCCTTCCCTACTGTTACAAAAGGCAGCGAGAATATGAAGAGGAAGATGACGATGATGATGAATACGAC
TCTGAAATGGAAGATTTTATTGAAGATGAAGGAGAGCCTCAGGAAGAAATATCCAAGCACATTAGAGAA
ATCTTTGGTTATGACCGAAAAAATACAAAGATGAAAGTGATTATGCCTTACGTTACATGAAAAGTAGT
TGGAAAGAGCAGCAGAAGGAAGAAGCAAGAGCTTAAGACTGGGTATGCAAGAGGACTTAGAGGAAATG
AGACGTGAAGAAGAAGAAATGAACGTCGAAGGGCCAAGAAGCTGAAGAGGCGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_194285

Insert Size: 2058 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_194285.2</u>
RefSeq Size:	5752 bp
RefSeq ORF:	2058 bp
Locus ID:	144108
UniProt ID:	<u>Q68D10</u>
Cytogenetics:	11p15.1
MW:	75.6 kDa
Gene Summary:	Histone chaperone that stabilizes pre-existing histone tetramers and regulates replication-independent histone exchange on chromatin (PubMed:26109053). Required for normal chromatin refolding in the coding region of transcribed genes, and for the suppression of spurious transcription (PubMed:26109053). Binds DNA and histones and promotes nucleosome assembly (in vitro) (PubMed:23378026, PubMed:26109053). Facilitates formation of tetrameric histone complexes containing histone H3 and H4 (PubMed:26109053). Modulates RNA polymerase 1-mediated transcription (By similarity). Binds DNA, with a preference for branched DNA species, such as Y-form DNA and Holliday junction DNA (PubMed:23378026).[UniProtKB/Swiss-Prot Function]