

Product datasheet for **SC314433**

ZNF569 (NM_152484) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF569 (NM_152484) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF569
Synonyms:	ZAP1; Zfp74; ZNF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314433 representing NM_152484.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGACTGAGTCCCAGGGAACAGTAACATTCAAAGATGTGGCTATCGACTTCACTCAGGAGGAGTGAAG
AGATTGGATCCTGCTCAGAGAAAACCTGTACCGGAATGTGATGCTAGAAAACATAACAACCTTAATCACA
GTAGGCTATCCGTTACCAAACCTGATGTGATTTTCAAATTGGAGCAAGAAGAAGAACCATGGGTGATG
GAGGAAGAAGTATTAAGGAGACACTGGCAAGGAGAAATATGGGGAGTTGATGAGCATCAGAAAAACCAG
GACAGACTTTTGAGACAAGTTGAAGTTAAATTCAGAAAACACTGACTGAAGAAAAGGCAATGAATGT
CAAAAGAAATTTGCAAATGTATTTCTCTGAACCTGATTTTTTCCCTTCCAGACACAATCTCTATGAG
TATGACTTATTTGGAAAGTGTGTTAGAACATAATTTGACTGTCATAAATGTGAATGCCTTATGAGA
AAGGAGCATTGTGAATATAATGAACCTGTGAAATCATATGGTAATAGCTCATCCCATTTTGTCAATTACC
CCCTTTAAGTGAATCATTGTGGAAGGCTTCAATCAGACTTTGGACCTCATCAGACATCTGAGAATT
CATACTGGAGAGAAGCCCTATGAATGTAGTAAGTGTAGAAAAGCCTTCAGTCACAAGGAAAAACTTATT
AAACATTATAAAATTCACAGTAGGGAGCAGTCTTACAAATGTAATGAATGTGGTAAAGCTTTCATTTAAA
ATGTCAAATCTCATTAGACATCAAAGAATTCATACTGGAGAGAAGCCCTATGCATGTAAGGAATGTGAG
AAGTCCTTCAGCCAGAAATCAAATCTTATTGATCATGAAAAATTCATACTGGAGAGAAACCTTATGAA
TGTAATGAGTGTGGAAGCATTGAGCCAGAAAGCAAGCCTCATTGCACATCAGAAAGTTCATACTGGG
GAGAAACCTTATGCATGTAATGAATGTGGTAAAGCCTTCCCTCGAATTGCATCCCTTGCTCTTCATATG
AGAAGTCATACAGGAGAAAAACCTTATAATGTGATAAATGTGGTAAAGCCTTCTCTCAGTTTTCCATG
CTTATTATACATGTTAGAATTCATACAGGTGAAAAACCTATGAATGTAATGAGTGTGGAAGGCTTC
TCTCAAAGCTCAGCCCTTACTGTACATATGAGAAGTCACACTGGTGAGAAACCTATGAATGTAAGGAA
TGCAGAAAAGCCTTCAGCCACAAGAAAAACCTTATTACACACCAGAAAATTCATACTAGAGAGAAACCT
TATGAGTGAATGAATGTGGAAAGCTTTTATACAGATGTCAAATCTTGTAGACACCAGAGAATTCAT
ACTGGGAAAAACCTATATATGTAAGGAATGTGGAAAGCCTTTAGCCAGAAATCAAATCTCATTGCT
CATGAAAAATTCATTCTGGAGAGAAACCTATGAATGCAATGAATGTGGTAAAGCCTTCAGCCAAAAG
CAAACTTCATTACATCAAAGTTCATACTGGAGAGAAACCTTATGATTGTAATGAATGTGGTAAA
GCCTTCTCTCAAATTGCATCCCTTACCCTTCATTTGAGAAAGTCATACAGGGGAAAAGCCTTATGAATGT
GATAAATGTGGTAAAGCCTTCTCTCAGTGCTCACTGCTTAATTTACATATGAGAAGTCACACAGGTGAG
AAGCCCTATGTATGTAATGAATGTGGAAAGCCTTCTCTCAAAGAACTTCCCTTATTGTGCACATGAGA
GGCCATACAGGTGAAAAACCTATGAATGTAATAAATGTGGAAGCCTTCTCCCAAAGCTCATCCCTT
ACTATACATATACGAGGACATACAGGTGAGAAACCTTCGACTGTAGTAATGTGGAAGCCTTCTCT
CAAATCTCATCTTACCTTTCATATGAGAAAACATACAGGTGAGAAGCCTTCACTGTATTGAGTGT
GGCAAGGCTTTCAGCCAAAAGTCGACCTTGTTAGACACCAGAGAATTCATACTCATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_152484

Insert Size: 2061 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152484.2</u>
RefSeq Size:	4071 bp
RefSeq ORF:	2061 bp
Locus ID:	148266
UniProt ID:	<u>Q5MCW4</u>
Cytogenetics:	19q13.12
Protein Families:	Transcription Factors
MW:	79.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]