

Product datasheet for **SC314434**

Calpastatin (CAST) (NM_173060) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calpastatin (CAST) (NM_173060) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAST
Synonyms:	BS-17; PLACK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314434 representing NM_173060.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATCCACAGAAAACCAAGGCTGTAAAAACAGAACCTGAGAAGAAGTCACAGTCAACCAAGCTGTCT
GTGGTTTCATGAGAAAAAATCCCAAGAAGGAAAGCCAAAAAGAACACACAGAGCCAAAAAGCCTACCCAAG
CAGGCATCAGATACAGGAAGTAACGATGCTCACAATAAAAAAGCAGTTTCCAGATCAGCTGAACAGCAG
CCATCAGAGAAATCAACAGAACCAAGACTAAACCACAAGACATGATTTCTGCTGGTGGAGAGAGTGTT
GCTGGTATCACTGCAATATCTGGCAAGCCGGGTGACAAGAAAAAGAAAAGAAATCATTAAACCCAGCT
GTGCCAGTTGAATCTAAACCGGATAAACCATCGGGAAAGTCAGGCATGGATGCTGCTTTGGATGACTTA
ATAGATACTTTAGGAGGACCTGAAGAACTGAAGAAGAAAATACAACGTATACTGGACCAGAAGTTTCA
GATCCAATGAGTTCACCTACATAGAGGAATTGGGTAAAAGAGAAGTCACAATTCCTCCAAAATATAGG
GAACTATTGGCTAAAAGGAAGGGATCACAGGGCCTCCTGCAGACTCTTCCGAAACCCATAGGGCCAGAT
GATGCTATAGACGCCCTTGTCTGACTTCACCTGTGGTGCCTACAGCTGCTGGAAGAAAACCTGAA
AAAGAGGAATCTACAGAAGTTTTAAAAGCTCAGTCAGCAGGGACAGTCAGAAGTGCCTGCCACCCCAA
GAGAAGAAAAGAAAGGTGGAGAAGGATACAATGAGTGATCAAGCACTCGAGGCTCTGTCCGCTTCACTG
GGCACCCGGCAAGCAGAACCTGAGCTCGACCTCCGCTCAATTAAGGAAGTCGATGAGGCAAAAGCTAAA
GAAGAAAACCTAGAGAAGTGTGGTGGAGGATGATGAAACAATCCCATCTGAGTACAGATTAACCAGCC
ACGGATAAAGATGGAAAACCACTATTGCCAGAGCCTGAAGAAAACCCAAAGCCTCGGAGTGAATCAGAA
CTCATTGATGAATTTTCAAGAGATTTTGACCGTCTGAATGTAAGAGAAAACCATCTAAGCCAACCTGAA
AAGACAGAAGAATCTAAGGCCGCTGCTCCAGCTCCTGTGTGGAGGCTGTGTGTGGACCTCCATGTGT
AGTATACAGTCAGCACCCCTGAGCCGGCTACCTTGAAGGGCACAGTGCCAGATGATGCTGTAGAAGCC
TTGGCTGATAGCCTGGGGAAAAAGGAAGCAGATCCAGAAGATGGAAAACCTGTGATGGATAAAGTCAAG
GAGAAGGCCAAAAGAAGACCCTGAAAAGCTTGGTAAAAAGAAAGAAACAATTCCTCCTGATTATAGA
TTAGAAGAGGTCAAGGATAAAGATGAAAAGCCACTCCTGCCAAAAGAGTCTAAGGAACAGCTTCCACCC
ATGAGTGAAGACTTCTTCTGGATGCTTTGTCTGAGGACTTCTCTGGTCCACAAAATGCTTCATCTCTT
AAATTTGAAGATGCTAACTTGCTGCTGCCATCTCTGAAGTGGTTTCCCAAACCCAGCTTCAACGACC
CAAGCTGGAGCCCAACCCGTGATACCTCGCAGAGTGACAAGACCTCGATGATGCCTTGGATAAACTC
TCTGACAGTCTAGGACAAAGGCAGCCTGACCCAGATGAGAACAACCAATGGAAGATAAAGTAAAGGAA
AAAGCTAAAGCTGAACATAGAGACAAGCTTGGAGAAAGAGATGACACTATCCCACCTGAATACAGACAT
CTCCTGGATGATAATGGACAGGACAAACCAAGTGAAGCCACCTACAAGAAATCAGAGGATTCAAAGAAA
CCTGCAGATGACCAAGACCCATTGATGCTCTCTCAGGAGATCTGGACAGCTGTCCCTCCACTACAGAA
ACCTCACAGAACACAGCAAAGGATAAGTGAAGAAGGCTGCTTCCAGCTCCAAAGCACCTAAGAATGGA
GGTAAAGCGAAGGATTGACGAAAGACAACAGAGGAAACTTCCAAGCCAAAAGATGACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_173060

Insert Size: 2061 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173060.3](#)

RefSeq Size: 4304 bp

RefSeq ORF: 2061 bp

Locus ID: 831

UniProt ID: [P20810](#)

Cytogenetics: 5q15

MW: 74.1 kDa

Gene Summary: The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]
Transcript Variant: This variant (2) represents use of an alternate promoter compared to variant 1.