

Product datasheet for **SC314442**

ILF3 (NM_153464) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ILF3 (NM_153464) Human Untagged Clone
Tag:	Tag Free
Symbol:	ILF3
Synonyms:	CBTF; DRBF; DRBP76; MMP4; MPHOSPH4; MPP4; MPP4110; NF-AT-90; NF90; NF90a; NF90b; NF90c; NF90ctv; NF110; NF110b; NFAR; NFAR-1; NFAR-2; NFAR2; NFAR90; NFAR110; TCP80; TCP110
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314442 representing NM_153464.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGTCCAATGCGAATTTTGTGAATGATGACCGCCATGTGATGGCAAAGCATTCTTCGTTTATCCA
ACACAAGAGGAGCTGGAGGCAGTCCAGAACATGGTGTCCACACGAGCGGGCGCTCAAAGCTGTGTCC
GACTGGATAGACGAGCAGGAAAAGGGTAGCAGCGAGCAGGACAGAGTCCGATAACATGGATGTGCCCCCA
GAGGACGACAGTAAAGAAGGGGCTGGGGAACAGAAGACGAGCAGACATGACCAGAACCCTGCGGGGAGTG
ATGCGGGTGGCCTGGTGGCAAAGGGCCTCCTACTCAAGGGGACTTGGATCTGGAGCTGGTGTCTGTG
TGTAAGGAGAAGCCACAACCGCCCTCTGGACAAGGTGGCCGACAACCTGGCCATCCAGCTTGTCTGT
GTAAACAGAAGACAAGTACGAAATACTGCAATCTGTGACGATGCTGCGATTGTGATAAAAAACAAAA
GAGCCTCCATTGTCCCTGACCATCCACCTGACATCCCCTGTTGTCAGAGAAGAAATGGAGAAAGTATTA
GCTGGAGAAACGCTATCAGTCAACGACCCCCGGACGTTCTGGACAGGCAGAAATGCCTTGTGCCTTG
GCGTCCCTCCGACAGCCAAGTGGTCCAGGCCAGAGCCAACGGGCTGAAGTCTTGTGTCATTGTGATC
CGGTCTTGAGGGACCTGTGCACTCGCGTGGCCACCTGGGGTCCCCTCCGAGGCTGGCCTCTCGAGCTC
CTGTGTGAGAAATCCATTGGCACGGCCAACAGACCGATGGGTGCTGGCGAGGCCCTGCGGAGAGTGTG
GAGTGCCTGGCGTCGGGCATCGTGATGCCAGATGGTTCTGGCATTATGACCCCTGTGAAAAAGAAGCC
ACTGATGCTATTGGGCATCTAGACAGACAGCAACGGGAAGATATCACACAGAGTGCAGCAGCAGCAGT
CGGCTCGCTGCCTTCGGCCAGCTCCATAAAGTCTAGGCATGGACCCTCTGCCTTCCAAGATGCCAAG
AAACCAAAGAATGAAAACCCAGTGGACTACACCGTTTCCAGTCCCACCAAGCACCCTATGCCATTACG
CCCATGAAACGCCAATGGAGGAGGACGGGGAGGAGAAGTCCGCCAGCAAAAAGAAGAAGATTGAG
AAGAAAGAGGAGAAGGCAGAGCCCCCAGGCTATGAATGCCCTGATGCGGTTGAACCAGCTGAAGCCA
GGGCTGCAGTACAAGCTGGTGTCCAGACTGGGCCCGTCCATGCCCCATCTTTACCATGTCTGTGGAG
GTTGATGGCAATTCATTCGAGGCCTCTGGGCCCTCCAAAAGACGGCCAAGCTGCACGTGGCCGTTAAG
GTGTTACAGGACATGGGCTTGGCGACGGGTGCTGAAGGCAGGGACTCGAGCAAGGGGGAGGACTCGGCT
GAGGAGACCGAGGCGAAGCCAGCAGTGGTGGCCCTGCCCCAGTGGTAGAAGCTGTCTCCACCCCTAGT
GCGGCCTTCCCTCAGATGCCACTGCCGAGCAGGGCCGATCCTGACAAAGCACGGCAAGAACCAGTC
ATGGAGCTGAACGAGAAGAGGCGTGGGCTCAAGTACGAGCTCATCTCCGAGACCGGGGGCAGCCACGAC
AAGCGCTTCGTCATGGAGTCAAGTGGATGGACAGAAGTCCAAGGTGCTGGTTCCAACAAAAAGGTG
GCGAAGGCCTACGCTGCTTGTGCTGCCCTAGAAAAGCTTTTCCCTGACACCCCTCTCGCCCTTGTGCC
AACAAAAAGAAGAGAGCCCAAGTACCCGTCAGAGGGGGACCGAAATTTGCTGCTAAGCCACATAACCT
GGCTTCGGCATGGGAGGCCCATGCACAACGAAGTGGCCCCACCCCAACCTTCGAGGGCGGGGAAGA
GGCGGGAGCATCCGGGGACGAGGGCGGGGCGAGGATTTGGTGGCGCAACCATGGAGGCTACATGAAT
GCCGGTGTGGGTATGGAAGCTATGGGTACGGAGGCAACTCGGCGACAGCAGGCTACACAGGTTTTGTA
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_153464

Insert Size: 2073 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_153464.2](#)

RefSeq Size: 3302 bp

RefSeq ORF: 2073 bp

Locus ID: 3609

UniProt ID: [Q12906](#)

Cytogenetics: 19p13.2

Protein Families: Druggable Genome, Transcription Factors

MW: 74.6 kDa

Gene Summary: This gene encodes a double-stranded RNA (dsRNA) binding protein that complexes with other proteins, dsRNAs, small noncoding RNAs, and mRNAs to regulate gene expression and stabilize mRNAs. This protein (NF90, ILF3) forms a heterodimer with a 45 kDa transcription factor (NF45, ILF2) required for T-cell expression of interleukin 2. This complex has been shown to affect the redistribution of nuclear mRNA to the cytoplasm. Knockdown of NF45 or NF90 protein retards cell growth, possibly by inhibition of mRNA stabilization. In contrast, an isoform (NF110) of this gene that is predominantly restricted to the nucleus has only minor effects on cell growth when its levels are reduced. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Dec 2014]
Transcript Variant: This variant (3) differs in the 3' UTR and coding region compared to variant 4. The resulting isoform (c) has a distinct C-terminus compared to isoform d. This isoform is also called DRBP76 delta and ILF3-C.