

Product datasheet for **SC314445**

TRAF3IP1 (NM_015650) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRAF3IP1 (NM_015650) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRAF3IP1
Synonyms:	FAP116; IFT54; MIP-T3; MIPT3; SLSN9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC314445 representing NM_015650.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACGCGCGGTGGTGAAGCGGACGACGAGGCGCTGGGAAAGTGATTTCGGAGGCCCGCTGACC
GAGAAGCTGCTGAGCAAGCCCCGTTCCGCTACCTGCACGACATCATCACGAGGTGATTAGAATGACT
GGTTTCATGAAGGCCTCTACACAGACGCCGAGATGAAGTCTGATAATGTGAAGGATAAAGATGCAAAA
ATTAGTTCCTACAAAAGGCCATAGACGTGGTTGTAATGGTGTCTGGGAGAGCCACTGTTGGCCAAACCA
GCCCCAATCGTGGCGGGCATGAGCCTGAAAGAACAACGAGCTGCTCCAGATAATTGAAAAATGCTGT
CTCAACAAGCTCTCTAGTGACGATGCGGTGCGGAGGGTTTTAGCTGGAGAGAAGGAGAAGTGAAAGGC
CGGGCCTCACTGACCTCAAGATCTCAGGAATTGGATAATAAGAATGTGCGAGAAGAAGAGTCCAGAGTT
CACAAAAATACAGAGGATAGAGGAGACGCTGAAATAAAAGAGAGAAGTACAAGCAGAGATCGAAAACAG
AAGGAAGAATTGAAAGAAGACCGCAAGCAAGAGAAAAGGACAAGGACAAGGAGAAGGCCAAGGAGAAT
GGCGGAAACAGACACAGAGAAGGGGAGAGAGAGAGAGCCAAAGCCGGGCCAGGCCAGACAACGAGCGA
CAGAAAGACAGAGGCAACAGGGAGCGGGACAGAGACTCCGAGCGCAAGAAGGAGACAGAGAGAAAGAGT
GAGGGGGGAAAGAGAAGGAGAGACTGAGAGACAGGGACCGAGAGCGCGACCCGGGACAAAGGGAAGGAC
AGGGACAGACGGAGAGTGAAAAACGGGGAGCACTCCTGGGACCTGGACAGGGAGAAGAACAGAGAGCAT
GACAAACCTGAGAAAAAGTCAGCAAGCTCAGGGGAGATGTCTAAAAAGTTATCAGATGGAACCTTTAAA
GACTCCAAGGCTGAAACAGAGACTGAGATTTCCACTAGAGCTTCCAAGTCATTGACAACAAAAACATCA
AAACGGCGATCCAAAAATTCAGTGAAGGAAGGAAGGAGGATAATATTTAGCTAAAAGTTTAGACTCC
ATAGTGTCTGGAATAAATATGAGCCAAATCAGGAAACGACAACATCAGAAATAGGAACAAAAGAAGCT
AATATTAACCTCACTAGTATTTTCAGATGATAATTCAGCTAGTCTGCGGTGTGAGAATATTCAGCCCAAC
CCCACAGAGAAGCAGAAAGGTGACTCCACCAGTGATGCAGAAGGAGATGCTGGACCTGCTGGCCAAGAT
AAGTCTGAGGTGCCAGAGACTCCAGAAATTCCTAATGAGCTTTCATCCAACATCAGAAGAATTCCTCGG
CCTGGGAGTGCAAGACCAGCCCCCTCCCGGGTCAAACGGCAAGACAGCATGGAGGCGCTACAAATGGAT
AGGTCAAGGAGTGGTAAACCGTTTCAAATGTGATTACAGAGTCACACAATTCGACAATGAAGAGGAT
GATCAATTTGGTGGAGCTGCCCTCAGCTCTCTGAAATGTCAGAAATGAAATGGTAACAGCAGTG
GAACTAGAAGAAGAGGAGAAGCATGGTGGACTTGTGAAAAAATTTGGAGACGAAGAAAGATTATGAG
AAATTGCAGCAGTCAACCAACCTGGGAGAAGGAGCGATCTCTTTGAGTCGGCATGGAAGAAGGAG
AAGGACATCGTTTCAAGGAGATAGAGAAGCTCCGCACGTCCATCCAGACCTGTGCAAGAGCGCACTT
CCCTGGGGAAGATCATGGACTACATCCAGGAAGACGTGGATGCCATGCAGAATGAGCTGCAGATGTGG
CACAGCGAGAACAGGCAGCACGCCGAGGCCCTGCAGCAGGAGCAGAGGATCACAGACTGTCCGTGGAG
CCCTTAAGGCTGAGCTCGCGGAGCTGGAGCAGCTGATCAAAGACCAGCAAGACAAGATCTGTGCTGTG
AAGGCCAACATCCTCAAGAAATGAAGAAAAATCCAGAAATGGTATATAGTATCAATTTGACTTCGAGA
AGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_015650
Insert Size: 2076 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015650.3</u>
RefSeq Size:	4183 bp
RefSeq ORF:	2076 bp
Locus ID:	26146
UniProt ID:	<u>Q8TDR0</u>
Cytogenetics:	2q37.3
Protein Families:	Druggable Genome
MW:	78.6 kDa
Gene Summary:	The protein encoded by this gene interacts with TNF receptor-associated factor 3, tethering it to cytoskeletal microtubules. The encoded protein is also an inhibitor of the innate type I IFN response. Defects in this gene are a cause of Senior-Loken syndrome 9. [provided by RefSeq, Mar 2017] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The extent of this RefSeq transcript is supported by transcript alignments.