

Product datasheet for **SC314448**

ZNF180 (NM_013256) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF180 (NM_013256) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF180
Synonyms:	HHZ168
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314448 representing NM_013256.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGCGCATGCCAGGGTCTACGCGGGAAGCTGGAAGCGGTGCTCAGGACCTCAGCACCCCTGCTGTGC
CTGGAGGAGAGCATGGAAGAGCAGGATGAGAAGCCCCCAGAGCCCCCAAGGTCTGTGCACAGATTCT
TTCCTTCTCAAGAGATTATCATCAAAGTCGAGGGAGAAGACACTGGGTCTCTGACCATCCCATCTCAG
GAAGGAGTGAACCTCAAATTTGACTGTGGACTTACACGCGGAGGAACAGGGTACTTGCAACCCTGCT
CAGAGGACCTGGACAGAGATGTGATCCTGGAGAACCACAGGGACCTAGTCTCTGGGACTTGGCAACT
GCAGTTGGAAAAAAGATTCAACTTCAAAGCAGAGGATTTTTGATGAAGAACCAGCTAATGGAGTGAAG
ATAGAAAGGTTTACAAGGGATGATCCTTGGTTATCTTCATGTGAAGAAGTGGATGATTGTAAGACCAG
TTGGAGAAGCAACAGGAAAAACAAGAGATACTTTTGCAGGAAGTGGCATTCACTCAAAGGAAAGCAGTT
ATTCATGAGAGAGTCTGCAAAAGTATGAAACTGGGAGAGAGTGGTCTGAATCCAGTCTATTTTCA
TCCCCAGTTATACCCATAAGAAACCATTTTCATAAACATGTATCACATGCTAAAAATGGCATCTTAAT
GCTGCTGTAAACAGTCACTAGAAGTAAATGAGAATGAGACACTATATGAAAAATGAATGTGGAAAA
CCCCCTCAGAGCATTACCTTATTCAGTTTACAAGAACTCAAACAAAAGATAAATCCTATGGATTTAGT
GACCGTATTCAATCTTTTTGCCATGGTACACCCCTACATATACATGAAAAAATTCATGGAGGAGGAAAA
ACCTTTGATTTTAAAGAAATGTGGGCAGGTTTTGAACCCAAAAATATCCCAATGAACAACAGAGAATT
CCTTTTGAAGAGAGTCAATATAAATGTAGTGAACCTCTCATAGTTCTCCCTTACTCAAACATGAGA
AATAATTCTGAAGAGAAACCTTTTGAATGTAATCAGTGTGGGAAATCCTTCAGCTGGAGCTCGCATCTT
GTTGCACATCAGAGAACTCACACAGGGGAGAAACCTTATGAATGTAGTGAATGTGGAAAAATCCTTCAGC
CGGAAATCCTTTAGCCAGAGTTATGTCCTTGTGTGCATCAAAGAACTCACTACTGGGAGAGAAGCCTTAT
GAATGCAATCAGTGTGGAAGTCATTACGCGAGAGCTATAAACTTATTGCACATCAAAGAACACATACC
GGAGAGAAGCCCTATGAATGTAATCAATGTGGGAAATCATTTATCCAGAGCTATAAACTTATTGCACAT
CAAAGAATTCATACTGGGAAAAACCTATGAATGCAATCAGTGTGGGAAATCCTTTAGTCAAAGTTAT
AACTTGTGCTCATCAGAGAACTCACACAGGAGAAAAACCTTTGATGTAATCAGTGTGGGAAATCC
TTCAGCTGGAGCTCTCAGCTTGTGCACATCAAAGAACTCACACTGGAGAGAAACCGTATGAATGTAGT
GAATGTGGAAAAATCTTTAACCAGCTTCTCACCTTGTATGCATCAGAGAATTCACACTGGGAAAAA
CCGTATGAATGTAATCAGTGTGGGAAATCCTTCAGCCAGAGTTATGTTCTTGTGTACATCAGAGAACT
CACTACTGGAGAAAAGCCCTATGAATGCAGTCAATGTGGGAAAGTCTTCAGACAGAGTTTATGCCTTACT
CAACATCAGAGAACTCATACTGGAGAGAAACATTTGATGTAATCAGTGTGGAAAAACATTTAGCTTG
AGTGCTCGACTTATTGTGCATCAAAGAACTCATACTGGAGAGAAACCTTTACATGTATTAGTGTGGA
AAAGCTTTCATTAATAGCTATAAACTTATTAGGCATCAGGCAACTCATACTGAAGAGAACTCTATGAA
TGTAAC TAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_013256

Insert Size: 2079 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013256.5](#)

RefSeq Size: 4336 bp

RefSeq ORF: 2079 bp

Locus ID: 7733

UniProt ID: [Q9UJW8](#)

Cytogenetics: 19q13.31

Domains: KRAB, zf-C2H2

Protein Families: Transcription Factors

MW: 79 kDa

Gene Summary: Zinc finger proteins have been shown to interact with nucleic acids and to have diverse functions. The zinc finger domain is a conserved amino acid sequence motif containing 2 specifically positioned cysteines and 2 histidines that are involved in coordinating zinc. Kruppel-related proteins form 1 family of zinc finger proteins. See MIM 604749 for additional information on zinc finger proteins.[supplied by OMIM, Jul 2002]
Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.