

Product datasheet for **SC314465**

SRPK2 (NM_182692) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SRPK2 (NM_182692) Human Untagged Clone
Tag:	Tag Free
Symbol:	SRPK2
Synonyms:	SFRSK2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314465 representing NM_182692.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAGCTCCCGAAAGTCTGGCCATTCAGGCCGAAAGCGGAGGCCGAAAAGAGAGAAACATCCGAAA
AAGCCGGAGCCTCAACAGAAAGCTCCTTTAGTTCTCTCCTCCACCGCCACCACCACCACCACCGCCA
CCTTTGCCAGACCCACACCCCGAGCCAGAGGAGAGATCCTGGGATCAGATGATGAGGAGCAAGAG
GACCTGCGGACTACTGCAAAGGTGGATATCATCCAGTAAAAATTGGAGACCTCTTCAATGGCCGGTAT
CATGTTATTAGAAAGCTTGGATGGGGGCACTTCTACTGTCTGGCTGTGCTGGGATATGCAAGGGAAA
AGATTTGTTGCAATGAAAGTTGTAAGTAAAGTGCCAGCATTATACGGAGACAGCCTTGGATGAAATAAAA
TTGCTCAAATGTGTTGAGAAAGTATCCAGTGACCCAAACAAAGACATGGTGGTCCAGCTCATTGAC
GACTTCAAGATTCAGGCATGAATGGGATACATGTCTGCATGGTCTTGAAGTACTTGGCCACCATCTC
CTCAAGTGGATCATCAAATCCAATCAAGGCCTCCAGTACGTTGTGTGAAGAGTATCATTGACAG
GTCCTTCAAGGGTTAGATTACTTACACAGTAAGTCAAGATCATTACTACTGACATAAAGCCGGAAAAT
ATCTTGATGTGTGGTGGATGATGCATATGTGAGAAGAATGGCAGCTGAGGCCACTGAGTGGCAGAAAGCA
GGTGTCTCTCTCTTCCAGGCTGTCAGTGAAGTACGGCTCCACAGCAGAAACCTATAGGAAAAATATCT
AAAAACAAAAAGAAAAACTGAAAAAGAAACAGAAAGAGGCAGGCTGAGTTATTGGAGAAGCGCCTGCAG
GAGATAGAAGAATTGGAGCGAGAAGCTGAAAGGAAAATAATAAGAAGAAACATCACCTCAGCTGCACCT
TCCAATGACCAGGATGGCGAATACTGCCAGAGGTGAAACTAAAAACAACAGGATTAGAGGAGGCGGCT
GAGGCAGAGACTGCAAAGGACAATGGTGAAGCTGAGGACCAGGAAGAGAAAGAAGATGCTGAGAAAGAA
AACATTGAAAAAGATGAAGATGATGTAGATCAGGAACCTGCGAACATAGACCCTACGTGGATAGAATCA
CCTAAAACCAATGGCCATATTGAGAATGGCCATTCTCACTGGAGCAGCAACTGGACGATGAAGATGAT
GATGAAAGAAGACTGCCCAAATCCTGAGGAATATAATCTTGATGAGCCAAATGCAGAAAGTATTACACA
TATAGCAGCTCCTATGAACAATCAATGGTGAATTGCCAAATGGACGACATAAAAATCCCGAGTCACAG
TTCCAGAGTTTTCCACCTCGTTGTTCTCTGGATCCTTAGAACCTGTGGCCTGCGGCTCTGTGCTTTCT
GAGGGATCACCCTTACTGAGCAAGAGGAGAGCAGTCCATCCCATGACAGAAGCAGAACGGTTTCAGCC
TCCAGTACTGGGATTTGCCAAAAGCAAAAACCCGGGCAGCTGACTTGTTGGTGAATCCCCTGGATCCG
CGGAATGCAGATAAAATTAGAGTAAAAATTGCTGACCTGGGAAATGCTTGTGGTGCATAAACACTTC
ACGGAAGACATCCAGACGCGTCAGTACCGCTCCATAGAGGTTTTAATAGGAGCGGGTACAGCACCCCT
GCGGACATCTGGAGCACGCGTGTATGGCATTGAGCTGGCAACGGGAGATTATTTGTTGAACCACAT
TCTGGGAAGACTATTCCAGAGACGAAGACCACATAGCCACATCATAGAGCTGCTAGGCAGTATCCA
AGGCACTTTGCTCTATCTGGAAAATATTCTCGGAATTCTTCAATCGCAGAGGAGAACTGCGACACATC
ACCAAGCTGAAGCCCTGGAGCCTCTTTGATGACTTGTGGAAAAGTATGGCTGGCCCCATGAAGATGCT
GCACAGTTTACAGATTTCTGATCCCGATGTTAGAAAATGGTTCAGAAAAACGAGCCTCAGCTGGCGAA
TGCCTTCGGCATCCTTGGTTGAATTCTTAG
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

Plasmid Map:


ACCN: NM_182692

Insert Size: 2100 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_182692.2](#)

RefSeq Size: 3703 bp

RefSeq ORF: 2100 bp

Locus ID: 6733

UniProt ID: [P78362](#)

Cytogenetics: 7q22.3

Protein Families:	Druggable Genome, Protein Kinase
MW:	79 kDa
Gene Summary:	<p>Serine/arginine-rich protein-specific kinase which specifically phosphorylates its substrates at serine residues located in regions rich in arginine/serine dipeptides, known as RS domains and is involved in the phosphorylation of SR splicing factors and the regulation of splicing. Promotes neuronal apoptosis by up-regulating cyclin-D1 (CCND1) expression. This is done by the phosphorylation of SRSF2, leading to the suppression of p53/TP53 phosphorylation thereby relieving the repressive effect of p53/TP53 on cyclin-D1 (CCND1) expression. Phosphorylates ACIN1, and redistributes it from the nuclear speckles to the nucleoplasm, resulting in cyclin A1 but not cyclin A2 up-regulation. Plays an essential role in spliceosomal B complex formation via the phosphorylation of DDX23/PRP28. Can mediate hepatitis B virus (HBV) core protein phosphorylation. Plays a negative role in the regulation of HBV replication through a mechanism not involving the phosphorylation of the core protein but by reducing the packaging efficiency of the pregenomic RNA (pgRNA) without affecting the formation of the viral core particles.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes isoform a.</p>