

## Product datasheet for **SC314509**

### ZDHC5 (NM\_015457) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZDHC5 (NM_015457) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZDHC5
Synonyms:	DHC5; ZNF375
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC314509 representing NM\_015457.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCCGCAGAGTCTGGAAAGAGATTCAAACCCAGCAAGTATGTCCCGGTCTCTGCAGCCGCCATCTTC
CTAGTGGGAGCTACGACACTCTTCTTTGCCTTTACGTGTCCAGGACTAAGCCTGTATGTGTACCTGCA
GTGCCATCTACAATGCAATTATGTTTCTTTGTGTTGGCCAACCTCAGCATGGCCACCTTCATGGAC
CCAGGGATTTTCCCTCGAGCTGAGGAGGATGAGGACAAGGAAGATGATTTCCGAGCTCCCTTTACAAA
ACAGTGGAGATAAAGGGCATCCAGGTGCGCATGAAATGGTGTGCCACCTGCCGCTTTTACCGTCCCCT
CGATGTTCCCACTGCAGTGTCTGTGACAACTGTGTGGAGGAATTTGATCATCACTGCCCTGGGTGAAT
AACTGTATTGGTCGCCGGAACCTACCGTTATTTTTCTTTTCTCTTTCCCTTTCCCTGACAGCCACATTATG
GGTGTGTTGGCTTTGGCTCCTTATGTCCTCTACCACATAGAGGAACCTCAGGGGTCCGCACGGCT
GTCACAATGGCAGTAATGTGTGGCTGGCTTATTTCTTCCCTGTAGCTGGCCTCACGGGATTTAC
GTGGTTCTGGTGGCCAGGGGACGCACAACCAATGAACAGGTTACGGGTAATTCCGGGGAGGTGTGAAC
CCCTTACCAATGGCTGCTGTAACAATGTACGCGTGTCTCTGAGTTCCTCAGCACCCAGGATTTG
GGGAGACCAAAGAAAGAGAAGACAATTGTAATCAGACCTCCCTTCCCTTCGACCAGAAGTTTCAGATGGG
CAGATAACTGTGAAGATCATGGATAATGGCATCCAGGGAGAGCTGAGGAGAACAAAGTCTAAGGGAAGC
CTGGAGATAACAGAGAGCCAGTCTGCAGATGCTGAACCTCCACCTCCTCCTAAGCCAGACCTGAGCCGT
TACACAGGGTTGCGAACACACCTCGGCTTGCTACTAATGAGGATAGTAGCTTATTGGCCAAGGACAGC
CCCCGACACCTACCATGTACAAGTATCGGCCGGGTTACAGTAGCAGCAGTACGTGAGCTGCCATGCCG
CATTCTCCAGCCCAAGTTGAGTCGTGGGACAGCTTGAAGGAGCAACCTCAATTGCAGAGAGCAGC
CGTACCCCAAGCTACCGCTCAGAGCCAGCTTGAAGGAGCAACCTCAATTGCAGAGAGCAGC
AGTTTTCACTTCGATCCACTATCCAGTGGCTCAGCTCCTCCAGCCTCAAGTCAGCCAGGGCACAGGC
TTTGAGCTGGCCAGTTGCAATCCATTCGTTTCCAGAGGGCACACCTCCACCTCCTATAAGAGCCTGGCC
AACCAGACACGCAATGGAAGCCTATCTTATGACAGCTTGTCTCACACCTTCAGACAGCCCTGATTTTGG
TCAGTGCAGGCAGGGCTGAGCCAGACCCACCTTAGGCTATACCTCCTCCTTCTGTCAGCCAGGCTG
GCCAGCAACGGGAAGCTGAGAGGCACCCACGTTTGGTGCCAACCTGGCCCAACACACCGAGAGCCCTCA
CCAGTCCGTTACGACAATCTGTCGCGCCACATTGTGGCTCTCTCCAGGAACGAGAGAAGTTGTGCGC
CAGTACCCCCACTCCCGGGCCGTGAGGAAGAACCAGGCTTGGGGGACTCAGGCATTGAGTCAACACCA
GGCTCGGGCCATGCCCTCGTACTAGTTCTCCTCAGATGATTCAAAGAGATCACCTTTGGGCAAGACT
CCACTGGGACGCCAGCTGTCCCGGTTTTGGCAAGCCAGATGGGCTAAGGGGCCGGGAGTAGGGTCC
CCTGAACCAGGCCAACAGCCCCATACCTGGGCGGATCGATGTCTTACAGCAGCCAAAAAGCCAACT
GGTGTCTCTGAGACAGAAGAAGTGGCCTTGACGCCATTACTGACACCCAAAGATGAAGTACAGCTGAAG
ACCACCTACAGCAAAATCCAACGGGCAGCCCAAGAGCTTAGGCTCAGCCTCCCTGGCCAGGCCAGCCA
CCTCTCAGTAGCCCCACGAGGGGAGGAGTCAAGAAGGTGTAGGGGTTGGTGGTACCACCTATGAGATT
TCGGTGTGA
ACGGTACGGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_015457

**Insert Size:** 2148 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_015457.2](#)

**RefSeq Size:** 4582 bp

**RefSeq ORF:** 2148 bp

**Locus ID:** 25921

**UniProt ID:** [Q9C0B5](#)

<b>Cytogenetics:</b>	11q12.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	77.5 kDa
<b>Gene Summary:</b>	Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5. Also palmitoylates FLOT2 (By similarity).[UniProtKB/Swiss-Prot Function]