

## Product datasheet for **SC314534**

### ELMO2 (NM\_182764) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ELMO2 (NM_182764) Human Untagged Clone
Tag:	Tag Free
Symbol:	ELMO2
Synonyms:	CED-12; Ced-12A; CED12; ELMO-2; VMPI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314534 representing NM\_182764.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCACCACCGTCAGACATTGTCAAAGTGGCCATTGAGTGGCCAGGTGCTAACGCCAGCTCCTTGAA
ATCGACCAGAAACGGCCCTGGCATCCATTATCAAGGAAGTTTGATGGGTGGTCGTTGCCAAACCCA
GAGTATTATACCCTCCGTTATGCAGATGGTCTCAGCTGTACATACCGAACAGACTCGCAGTGACATT
AAGAATGGGACAATCTTACAACCTGGCTATCTCCCGTCCCGGGCTGCACGCCAGCTGATGGAGAGGACC
CAGTCATCCAACATGGAGACCCGGCTGGATGCCATGAAGGAGCTGGCCAAGCTCTCTGCCAGCTGACT
TTCGCTACTGAGTTCATCAACATGGATGGCATCATTGTGCTGACAAGGCTCGTGGAAAGTGAACCAAG
CTCTTGCCCACTACAGTGAGATGCTGGCATTACCCCTGACTGCCTTCTAGAGCTCATGGACCATGGC
ATTGTCTCCTGGGACATGGTTTCAATCACCTTTATTAAGCAGATTGCAGGGTATGTGAGCCAGCCCATG
GTGGACGTGTCAATCCTCAGAGGTCCCTGGCCATCCTGGAGAGCATGGTCTTGAACAGCCAGAGTCTG
TACCAGAAGATAGCCGAGGAAATCACCGTGGGACAGCTCATCTCACACCTCCAGGTCTCCAACCAGGAG
ATTCAGACCTACGCCATTGCACTGATTAATGCACTTTTTCTGAAGGCTCCTGAGGACAAAACGACAGGAT
ATGGCAAATGCATTTGCACAGAAGCATCTCCGGTCTATAATCCTGAATCATGTGATCCGAGGGAACCCG
CCCATCAAACACTGAGATGGCCCATCAGCTATATGTCTTCAAGTCTAACCTTTAACCTTCTGGAAGAA
AGGATGATGACCAAGATGGACCCCAATGACCAGGCTCAAAGGGACATCATATTTGAACTGAGGAGGATT
GCATTTGACGCAGAGTCTGATCCTAGCAATGCCCTGGGAGTGGGACCGAAAAACGAAAGCCATGTAC
ACAAAGGACTACAAAATGCTGGGATTTACCAACCACATCAATCCAGCCATGGACTTTACCCAGACTCCT
CCTTGGAAATGCTGGCCTTGGACAACATGCTGTACTTGGCTAAAGTCCACCAGGACACCTACATCCGGATT
GTCTTGGAGAACAGTAGCCGGGAAGACAACATGAATGCCCTTTGGCCGAGTGCCATTGAGCTCAC
AAAATGCTCTGTGAAATCCTGCAGGTTGGGGAACACCAATGAAGGACGCAATGACTACCACCCGATG
TCTTTACCCATGACCGAGCCTTTGAAGAGCTCTTTGGAATCTGCATCCAGCTGTTGAACAAGACCTGG
AAGGAGATGAGGGCAACAGCAGAGGACTTCAACAAGGTTATGCAAGTCGTCGAGAGCAAATCACTCGA
GCTTTGCCCTCAAACCCAACCTTTGGATCAGTTCAGAGCAAATTCGTAGCCTGAGTACTCTGAG
ATTCTACGACTGCGCCAGTCTGAGAGGATGAGTCAGGATGACTTCCAGTCCCGCCAATTGTGGAGCTG
AGGGAGAAGATCCAGCCGAGATCCTTGAGCTGATCAAGCAGCAGCGCTGAACCGGCTCTGTGAGGGC
AGCAGCTTCCGAAAGATTGGGAACCGCCGAAGGCAAGAACGGTTCTGGTACTGCCGTTGGCACTGAAC
CACAAGGTCTTCACTATGGTACTTGGATGACAACCCACAAGGGGAGGTGACATTTGAATCCCTGCAG
GAGAAAATTCCTGTTGCAGACATTAAGGCCATTGTCACTGGGAAAGATTGCCCCACATGAAAGAGAAA
AGTGCTCTGAAACAGAACAAAGGAGGTGTTGGAATTGGCCTTCTCCATCCTGTATGACCCTGATGAGACC
TTAAACTTCATCGCACCTAATAAATATGAGTACTGCATCTGGATTGATGGCCTCAGTGCCCTTCTGGGG
AAGGACATGTCCAGTGAGCTGACCAAGAGTGACCTGGACACCCTGCTGAGCATGGAGATGAAGCTGCGG
CTCCTGGACCTGGAGAACATCCAGATTCCTCGAAGCCCCACCCCCATCCCCAAGGAGCCAGCAGCTAT
GACTTTGTCTATCACTATGGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_182764

**Insert Size:** 2163 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_182764.2](#)

**RefSeq Size:** 4021 bp

**RefSeq ORF:** 2163 bp

**Locus ID:** 63916

**UniProt ID:** [Q96JJ3](#)

**Cytogenetics:** 20q13.12

**MW:** 82.6 kDa

**Gene Summary:** The protein encoded by this gene interacts with the dedicator of cyto-kinesis 1 protein. Similarity to a *C. elegans* protein suggests that this protein may function in phagocytosis of apoptotic cells and in cell migration. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015]

Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.