

## Product datasheet for **SC314549**

### VWA2 (NM\_198496) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	VWA2 (NM_198496) Human Untagged Clone
Tag:	Tag Free
Symbol:	VWA2
Synonyms:	AMACO; CCSP-2; NET42
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC314549 representing NM\_198496.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCCCCTTTCCCTGTTGCTGGAAGCCGCTGTGTTTTCCCTGTTTTCCAGAGTGCCCCCATCTCTCCCT
CTCCAGGAAGTCCATGTAAGCAAAGAAACCATCGGGAAGATTTTCAGCTGCCAGCAAAATGATGTGGTGC
TCGGTGCAGTGGACATCATGTTTCTGTTAGATGGGTCTAACAGCGTCGGGAAAGGGAGCTTTGAAAGG
TCCAAGCACTTTGCCATCACAGTCTGTGACGGTCTGGACATCAGCCCCGAGAGGGTCAGAGTGGGAGCA
TTCCAGTTCAGTTCCTACTCCTCATCTGGAATTCCTTGGATTCATTTTCAACCAACAGGAAGTGAAG
GCAAGAATCAAGAGGATGGTTTTCAAAGGAGGGCGCACGGAGACGGAACCTGCTCTGAAATACCTTCTG
CACAGAGGGTTGCCCTGGAGGCAGAAATGCTTCTGTGCCCCAGATCCTCATCATCGTCACTGATGGGAAG
TCCCAGGGGGATGTGGCACTGCCATCCAAGCAGCTGAAGGAAAGGGGTGCTACTGTGTTTGTGTGGGG
GTCAGGTTTCCCAGGTGGGAGGAGCTGCATGCACTGGCCAGCGAGCCTAGAGGGCAGCACGTGCTGTTG
GCTGAGCAGGTGGAGGATGCCACCAACGGCCTTTCAGCACCTCAGCAGCTCGGCCATCTGCTCCAGC
GCCACGCCAGACTGCAGGGTCGAGGCTCACCCCTGTGAGCACAGGACGCTGGAGATGGTCCGGGAGTTC
GCTGGCAATGCCCATGCTGGAGAGGATCGCGGGACCCTTGGCGTGTGGCTGCACACTGTCCCTTC
TACAGCTGGAAGAGAGTGTTCCTAACCCACCCTGCCACCTGCTACAGGACCACCTGCCAGGCCCTGT
GACTCGCAGCCCTGCCAGAATGGAGGCACATGTGTTCCAGAAGGACTGGACGGCTACCAAGTGCCTCTGC
CCGCTGGCCTTTGGAGGGGAGGCTAACTGTGCCCTGAAGCTGAGCCTGGAATGCAGGGTCGACCTCCTC
TTCCTGCTGGACAGCTCTGCGGGCACCCTCTGGACGGTTCCTGCGGGCCAAAGTCTTCGTGAAGCGG
TTTGTGCGGGCCGTGCTGAGCGAGGACTCTCGGGCCGAGTGGGTGTGGCCACATACAGCAGGGAGCTG
CTGGTGGCGGTGCCTGTGGGGAGTACCAGGATGTGCTGACCTGGTCTGGAGCCTCGATGGCATTCCC
TTCCGTGGTGGCCCCACCCTGACGGGCAGTGCCTTGGCGCAGGGCGCAGAGCGTGGCTTCGGGAGCGCC
ACCAGGACAGGCCAGGACCGGCCACGTAGAGTGGTGGTTTTGCTCACTGAGTCACACTCCGAGGATGAG
GTTGCGGGCCAGCGCTCACGCAAGGGCGGAGAGCTGCTCCTGCTGGGTGTAGGCAGTGGGCGGTG
CGGGCAGAGCTGGAGGAGATCACAGGCAGCCAAAGCATGTGATGGTCTACTCGGATCCTCAGGATCTG
TTCAACCAAAATCCCTGAGCTGCAGGGGAAGCTGTGCAGCCGGCAGCGGCCAGGGTCCGGACACAAGCC
CTGGACCTCGTCTTCATGTTGGACACCTCTGCCTCAGTAGGGCCCGAGAATTTTGTCTCAGATGCAGAGC
TTTGTGAGAAGCTGTGCCCTCCAGTTTGGAGTGAACCCTGACGTGACACAGGTCGGCCTGGTGGTGTAT
GGCAGCCAGGTGCAGACTGCCTTCGGGCTGGACACCAAAACCACCGGGCTGCGATGCTGCGGGCCATT
AGCCAGGCCCCCTACCTAGGTGGGGTGGGCTCAGCCGGCACCGCCCTGCTGCACATCTATGACAAAGTG
ATGACCGTCCAGAGGGGTGCCCGCCCTGGTGTCCCAAAGCTGTGGTGGTGTCTACAGGGCGGAGAGGC
GCAGAGGATGCAGCCGTTCTGCCAGAAGCTGAGGAACAATGGCATCTCTGTCTTGGTCTGGGCGTG
GGGCTGTCTAAGTGAAGGCTGCGGAGGCTTCAGGTCCCCGGGATTCCCTGATCCACGTGGCAGCT
TACGCCGACCTGCGGTACCACAGGACGTGCTCATTGAGTGGCTGTGTGGAGGTGAGTGGGGGAATCCA
CACCTCAGGGCTGCCCCATGGCAGGCCCTCAGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-Mlul

**Plasmid Map:**


**ACCN:** NM\_198496

**Insert Size:** 2178 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_198496.1](#)

**RefSeq Size:** 2946 bp

**RefSeq ORF:** 2178 bp

**Locus ID:** 340706

**Cytogenetics:** 10q25.3

**MW:** 78.4 kDa

**Gene Summary:** This gene encodes a member of the von Willebrand factor A-like domain protein superfamily. The encoded protein is localized to the extracellular matrix and may serve as a structural component in basement membranes or in anchoring structures on scaffolds of collagen VII or fibrillin. This gene has been linked to type 1A diabetes and is a candidate serological marker for colon cancer. [provided by RefSeq, Jan 2013]