

## Product datasheet for **SC314565**

### **RAP1GAP2 (NM\_015085) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	RAP1GAP2 (NM_015085) Human Untagged Clone
Tag:	Tag Free
Symbol:	RAP1GAP2
Synonyms:	GARNL4; RAP1GA3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC314565 representing NM\_015085.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGTTTGGCCGGAAGCGCAGTGTCTCCTTTGGGGCTTCGGATGGATCGACAAGACCATGCTGGCAAGT
CTGAAGTCAAGAAGCAGGAGCTGGCCAACAGCTCGGATGCGACCCTCCAGACCGGCCGCTCTCCCT
CCTCTCACGGCACCTCCACCATGAAGTCGTCGGAGTTCTTTGAGATGCTGGAGAAAATGCAGGGGATC
AAGCTTGAAGAGCAGAAGCCGGGACCCAGAAGAACAAGGACGACTATATCCCATACCCAGCATCGAC
GAGTTGTGGAGAAGGGAGGCCGTACCCTCAGGTCATCCTGCCACAGTTTGGGGCTATTGGATCGAG
GACCCGGAGAACGTGGGCACCCCAACATCGCTGGGGAGCAGCATCTGTGAGGAGGAGGAAGAGGACAAC
CTCAGCCCAACACATTTGGCTACAAGCTCGAGTGAAGGGTGAAGCCAGGGCCTACCGGAGGCACTTC
CTGGGGAAGGATCATCTAACTTTACTGTACCGGCAGCAGCCTGGGAACTTGATCCTGTCCGTCAAG
TGCGAGGAAGCAGAGGGGATCGAGTACCTCCGGGTCATCCTTAGTCCAACTGAAGACGGTACATGAG
CGGATCCCCTTGCTGGACTGAGCAAGCTTCCAGTGTCCCTCAGATTGCAAAGGCTTTCTGTGATGAT
GCAGTGGGACTGAGATTCATCCTGTCCGTACCCCAAGGCCCTCCAAATGATTGTGCCTATGATGAG
CATGAAGTCAACAACACATTCAAATTCGGAGTCATTTATCAAAAAGCCAGGCAGACCCTGGAGGAGGAG
CTATTTGGGAACAATGAGGAGAGCCAGCTTTAAGGAGTTCTTGGACCTGCTGGGGACACGATCACA
CTGCAGGATTTCAAAGGTTTCCGAGGAGGCTGGACGTGACCCACGGACAGACAGGGGTGGAATCAGTG
TACACAACATTCGGGACAGGGAGATCATGTTTACGTTTCCACAAGCTGCCATTTACCGACGGAGAC
GCCAGCAGCTCCAGAGAAAGAGACACATTGAAATGACATCGTGCCATCATCTTCAAAGGAAAAC
ACGCCGTTTGTCCAGACATGATAGCTCCAATTTCTTACATGCCTACATCGTCGTGCAGGTCGAGACC
CCAGGCACAGAGACCCCATCCTACAAGTCTCTGTCACTGCGCGGGAAGATGTGCCACCTTTGGTCCA
CCTCTGCCAGTCCCCCGTTTTCCAGAAGGGCCCGGAATTCAGGGAGTTTCTGTCTACCAAGCTCACC
AATGCCGAGAACGCTGTGCAAGTCGGACAAGTTTGAAGCTGGAGGACCGGACCAGGCTGCCCTC
CTGGACAACCTTCACGATGAGCTCCACGCCACACAGGCCATGCTGGGACTGGGCCAGAGGAGGAC
AAGTTTGAGAAATGGAGGCCACGGGGGTTCTGGAGTCTTTAAGAGGGCCATCCGCTACGCAGCCAC
TCCATGGAGACCATGGTGGCGGCCAGAAGAAGTCGCACAGTGGGGCATCCCTGGCAGCCTCAGCGGG
GGCATCTCCACAACAGCATGGAGGTCACCAAGACCACCTTCTCGCTCCAGTGGTGGCGGAACGGTG
AAGAACCAGTACGGAGTCCCATCAAGCGACGCTCGGGGCTTCCCCCGCTGCACACGGGCTCAGAA
GGCCAGGGCGACAGCCGGGACGATGTGACAGCACATCCAGCACCCCAAGACCCAGATGGTGGACAC
TCCTCTCAGGAGATAAAGTCTGAGACCTCATCCAATCCCAGCTCTCCGGAATCTGCCCAACAAGGAG
AAGCCCTTCATGAAGTTGAAGGAAAACGGCGTGCCATCTCCCGTCTCTCCTCCAGCACCAGCAGGTC
AGCAGCACTGCAGGGGAGGGCGAGGCCATGGAGGAGGGGACAGTGGGGGAGCCAGCCGTCCACGACC
TCACCTTCAAGCAGGAGGTGTTTGTCTACAGCCCGTCCCGAGCAGCGAGGCCAGCCTGGGGCA
GCTGCCACCCGATCATCATGAGCCGAGTCCCACAGATGCCAAAAGCAGAAAACCTCCCGAGATCGAAC
CTGAAATCCGCTTTGACAAGCTCAGCCATGCCAGCTCTGGTGGGGTCACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_015085

**Insert Size:** 2193 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_015085.4](#)

**RefSeq Size:** 6680 bp

**RefSeq ORF:** 2193 bp

**Locus ID:** 23108

**UniProt ID:** [Q684P5](#)

**Cytogenetics:** 17p13.3

**MW:** 80.1 kDa

**Gene Summary:** This gene encodes a GTPase-activating protein that activates the small guanine-nucleotide-binding protein Rap1 in platelets. The protein interacts with synaptotagmin-like protein 1 and Rab27 and regulates secretion of dense granules from platelets at sites of endothelial damage. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]  
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). This variant is also known as variant b, as described in PMID: 15632203. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.