

## Product datasheet for **SC314645**

### **ZNF509 (ZBTB49) (NM\_145291) Human Untagged Clone**

#### **Product data:**

|                           |                                                  |
|---------------------------|--------------------------------------------------|
| Product Type:             | Expression Plasmids                              |
| Product Name:             | ZNF509 (ZBTB49) (NM_145291) Human Untagged Clone |
| Tag:                      | Tag Free                                         |
| Symbol:                   | ZNF509                                           |
| Synonyms:                 | ZNF509                                           |
| Mammalian Cell Selection: | Neomycin                                         |
| Vector:                   | pCMV6-Entry (PS100001)                           |
| E. coli Selection:        | Kanamycin (25 ug/mL)                             |



[View online »](#)

**Fully Sequenced ORF:** >SC314645 representing NM\_145291.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACCCCTGTTGCTACCCACAGCTGCCATCTGCTCCAGCAACTGCATGAGCAGCGAATCCAAGGCCTG
CTTTGTGACTGTATGTTGGTGGTAAAAGGAGTCTGCTTTAAAGCGCATAAGAATGTCCTGGCAGCATTC
AGCCAGTATTTTAGGAGCCTCTTTCAGAATTCCTCAAGCCAGAAGAATGATGTTTTTCACTTGGATGTT
AAAAATGTCAGTGGCATAGGGCAGATCCTGGACTTCATGTACACTTCTCATCTAGATCTTAACCAGGAC
AATATACAAGTAATGCTGGACACAGCACAGTGTGCAAGTTCAAAATGTTCTGAGTCTGTGTACACACA
TTTTTAAATCAGCCACTGTAGTACAGCCACCTGGCATGCCTTGAATAGTACATTGTCTCTACAAAGC
ACCCTGACCCAGATGCCACTTGTGTTATCAGTGAAAACCTACCCCTCATTTACTGCAGGAATGTTCA
GCAGATGCACAGCAGAACAAAACGTTGGATGAATCGCATCCGCATGCTTACCATCAGTTAATCGTCAT
CACTCCGAGGTGAAATCTCAAAACAAGCTCCTGATACTTCAGATGGCAGCTGCACAGAAGTGCCTTTC
AAACAGCCAAATTACTATTACAAACTCAGAACTTTTACAGTAAGCAGTACCATAAACACGCAGCTGGT
CCCAGTCAGGAGAGAGTTGTTGAGCAGCCTTTTGCCTTTCAGCACCTCTACAGACCTTACCACGGTAGAG
AGCCAGCCTTGTGCCGTGAGTCATTCTGAATGCATCCTGGAGTCTCCCGAGCACTTACCTTCCAACCTTC
CTGGCCAGCCTGTGAATGACTCTGCCCCACACCTGAGTCAGACGCCACATGCCAACAACTGTCAAG
CAGATGAGGCTCAAAAAGGCCATTATCTGAAGAAGCTCAATTTCTGAAGTACAGAAATACGCAGAG
CAAGTATCTGAACCAAGTCAGATGATGGTTTGACAAAGAGGTTGGAATCTGCTAGTAAAAATACCTTA
GAGAAAGCTAGCAGCCAAAGTGTGAAGAAAAAGAAAGTGAAGAAGTCGTGAGTTGTGAGATTTTAAAT
TGCATTAGTGAGACGGAGAGGCTGAAGACCCGGCTGCCCTGGAAGACCAGTCCAGACACTTCAGTCC
CAGAGACAATACGCGTGTGAATTATGCGGGAACCTTTTAAACACCCAAGCAACTTGGAGCTTCAAAA
CGGTCTCATACAGGTGAGAAACCTTTTGAATGTAAACATTTGTGGGAAACATTTCTCAGGCAGGTAAC
TTGCAGACTCACTTACGACGGCATTCTGGTGAAAAACCATACATCTGCGAGATCTGTGGAAGAGGTTT
GCAGCCTCTGGCGACGTCCAGCGTCACATTATTATTCACTCAGGAGAAAAACCACTTGTGTGACATC
TGTGGTCGAGGTTTAGTAACCTCAGTAATTTGAAGGAGCAGAAAAAGACACACACGGCTGATAAAGTC
TTCACCTGTGATGAGTGTGGAAGTCTTTTAAATATGCAAAGGAAGTTAGTAAAGCACAGAATTCCGCAC
ACGGGGGAGCGGCTTACAGCTGCTCTGCCTGCGGGAATGTTTTGGGGGATCAGGTGACCTCCGCAGG
CATGTCCGCACTCACACTGGGAGAAAGCGTACACATGTGAGATCTGTAACAAGTGCTTTACCCGCTCT
GCGGTGCTCCGGCGGCACAAGAAGATGCACTGCAAAGCTGGTGACGAGAGCCAGATGTGCTGGAGGAG
CTCAGCCAAGCCATCGAGACCTCCGACCTCGAGAAATCTCAGAGCTCAGACTCTTTCTCCCAAGACACG
TCTGTGACGCTGATGCCAGTGTCGGTTAAACTCCCTGTCCACCCAGTGGAAAAATTCTGTGGCAGAATTT
GATAGCCACTCTGGCGGCTCCTATTGTAAGTTACGGTCCATGATCCAACCTCATGGAGTTAGTGACCAG
GAGAAGCTGAGTTTGGATCCTGGTAAACTTGCCAGCCCGAGATGCAGCAGACACAGCCTCAGGCCTAT
GCTTACTCGGATGTGGACACCCAGCCGGTGGCGAACCCTGCAGGCCGATGGCATGGCCATGATCCGT
TCCTCTCTGGTGCTTTTGGACAACACGGCGGTGACCCCTGGGCAGTCGAGCATCTTCCACCATTAT
AGGAACCTCAGAGGGTCAGTTTTTCTCCAGCATGACTCTCTGGGGGCTAGCGATGAAGACGCTGCAGAAT
GAAAACGAGTTAGACCACTGA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:** SgfI-RsrII

**ACCN:** NM\_145291

**Insert Size:** 2298 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <a href="#">NM_145291.3</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>           | 2960 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2298 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 166793                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <a href="#">Q6ZSB9</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 4p16.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>MW:</b>                    | 85.1 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Gene Summary:</b>          | <p>Transcription factor. Inhibits cell proliferation by activating either CDKN1A/p21 transcription or RB1 transcription.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Both variants 1 and 2 encode the same protein.</p>                                                                                                                                                                                                                               |