

Product datasheet for **SC314680**

ACAD11 (NM_032169) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACAD11 (NM_032169) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACAD11
Synonyms:	ACAD-11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314680 representing NM_032169.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAGCCAGGTGCTACTGGCGAGTCCGATTTGGCCGAAGTCTGCCCCAGCACAAAGTTCGACAGCAAG
TCCCTGGAGGCCCTACCTAAACCAGCACTTGTCTGGCTTTGGGGCCGAACGTGAGGCTACGCTGACCATT
GCCAGTACAGAGCAGGAAAGTCCAATCCAACCTTTTATCTCCAGAAGGGCTTTCAAACATATGTGCTC
AGGAAAAAACCACCAGGTTCACTTCTTCTAAAGCACATCAGATTGATAGAGAATTTAAAGTCCAGAAA
GCCTTGTTTTCAATTGGATTCCCGTTCCTCAAGCCTATACTGACTGCAGTGATACTTCTGTCATTGGA
ACAGAATTTTACGTAATGGAACATGTGCAGGGTCAATCTTCCGTGATTTAAACAATTCCTGGACTTAGC
CCAGCAGAACGTTCCAGCCATATATGTGGCCACGGTAGAAAATTGGCTCAGTTACATTCTTGAATATA
CAGTCACTGCAGCTGGAAGGATATGGTATAGTGCTGGTACTGCAAAAGACAGGTATCAACCTGGACA
AAGCAATATCAAGCTGCAGCTCATCAGGACATCCCTGCCATGCAACAGCTATCGGAGTGGCTAATGAAG
AACTTGCCCGATAATGACAATGAAGAGAATTTGATTCATGGAGATTTAGACTAGATAACATAGTTTTC
CACCCTAAAGAGTGTGCGATTATAGCAGTGTGGATTGGGAGCTGTCAACCATTGGTCATCTTTGTCA
GACTTAGCTCATTTTTCCCTGTTCTACTTTGGCCAAGGACAGTTCCAATGATAAATCAAGTTCCTTAT
AGTGAAAACCTCAGGGATACCATCAATGGAAGAAGTGAATTTCAATATATTGCCGCTGCAGGGGAATTAAT
TCTATTCTTCTAACTGGAATTTCTTCTTCCCTTTTCAATTTTAAGATGGCTGGAATAGCACAGGGA
GTATATAGCAGATATCTTCTGGGAAATAATTCACTGAGGATAGCTTTTTATTGGCAATATTGTGCAA
CCTCTGGCAGAACTGGACTACAACCTCCAAACGAACTTTCAGTACTGTACTACCACAGATTGATACT
ACTGGACAGTTGTTTTGTACAGACTCGGAAAGGTCAGGAAGTCTTATTAAGGTGAAGCATTTCATGAAA
CAACACATTTCTCCAGCTGAAAAGGAGGTAAGTGAATTTCTATGTTCAAAATGAAAATTCAGTGGACAAG
TGGGAAAAACCTTTAGTGATTGATAAACTCAAGGAAATGGCCAAAGTCGAGGGTCTCGGAACCTGTTT
TTGCCAGCTGTCAGCGGACTCAGCCACGTGGACTATGCCTTGATTGCTGAAGAAACAGGAAAATGCTTT
TTTGCTCCAGATGCTTTAACTGCCAAGCACCAGACACAGGGAATATGGAGTTCTGCACCTGTATGGA
AGTGAGGAACAGAAAGAACAGTGGCTTGAGCCTTCTTCAAGGGAACATTACCTCTTGCTTCTGTATG
ACAGAACCTGATGTAGCTTCAAGTGATGCCACGAATTTGAATGCAGCATCCAACGAGATGAAGATAGC
TATGTAATTAACGGCAAAAAATGGTGGAGCAGTGGAGCTGGGAATCCAAGTGCAAAATGCAATTGTT
TTGGGAAGAACTCAAAATACTTCTCTCCAGACACAAACAGCACAGCATGATTCTTGTCCCATGAAC
ACACCTGGAGTAAAAATAATAAGCCTTTGTCAAGTTTTTGGCTACACAGATAATTTTCATGGAGGACAT
TTTGAGATCCATTTAATCAAGTGCAGATTCTGCCACAATCTAATACTAGGTGAAGGTAGGGGATTT
GAAATTTCCCAAGGCCGCTTGGACCTGGCAGAATCCACCACTGTATGAGAACAGTAGGTTTGGCGGAA
CGCGCTTTGCAGATCATGTGTGAGCGGGCAACACAAAGGATAGCTTTCAAGAAAGAGTTGTATGCACAT
GAGGTTGTGGCTCACTGGATTGCTGAAAGCCGCTTGGCATTGAGAAGATCCGCTTGTGACTCTGAAA
GCTGCTCACAGCATGGACTCTGGGCAGTGTGGCGCTAAGAAAGAGATTGCAATGATCAAAGTGGCT
GCCCCACGGGCTGTCAGCAAAATCGTTGACTGGCCATCCAGGTGTGCGGAGGTGCTGGTGTTCACAG
GATTACCTCTGGCTAACATGTATGCTATAACCCGAGTTTTGCGTTTAGCAGATGGACCTGACGAAGTT
CATCTTTACGAATCGCAACAATGGAGCTGCGGGACCAAGCCAAAAGACTGACAGCCAAGATATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_032169

Insert Size: 2343 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032169.4](#)

RefSeq Size: 3556 bp

RefSeq ORF: 2343 bp

Locus ID: 84129

UniProt ID: [Q709F0](#)

Cytogenetics: 3q22.1

Domains: Acyl-CoA_dh, Acyl-CoA_dh_M

Protein Families: Druggable Genome

MW: 87.3 kDa

Gene Summary: This gene encodes an acyl-CoA dehydrogenase enzyme with a preference for carbon chain lengths between 20 and 26. Naturally occurring read-through transcription occurs between the upstream gene NPHP3 (nephronophthisis 3 (adolescent)) and this gene. [provided by RefSeq, Aug 2015]

Transcript Variant: This variant (1) encodes the protein.