

Product datasheet for **SC314687**

GCSF Receptor (CSF3R) (NM_172313) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GCSF Receptor (CSF3R) (NM_172313) Human Untagged Clone
Tag:	Tag Free
Symbol:	GCSF Receptor
Synonyms:	CD114; GCSFR; SCN7
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC314687 representing NM_172313.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGGCAAGGCTGGGAACTGCAGCCTGACTTGGGCTGCCCTGATCATCTGCTGCTCCCCGGAAGTCTG
GAGGAGTGCAGGACATCAGTGTCTCAGCCCCATCGTCCACCTGGGGGATCCCATCACAGCCTCCTGC
ATCATCAAGCAGAAGTGCAGCCATCTGGACCCGGAGCCACAGATTCTGTGGAGACTGGGAGCAGAGCTT
CAGCCCCGGGGCAGGCAGCAGCCTGTCTGTGATGGGACCCAGGAATCTATCATCACCTGCCACCTC
AACACACTCAGGCCTTCTCTCCTGCTGCCTGAACTGGGGCAACAGCCTGCAGATCCTGGACCAGGTT
GAGCTGCGCGCAGGCTACCCTCCAGCCATACCCACAACCTCTCCTGCCTCATGAACCTCACAACCAGC
AGCCTCATCTGCCAGTGGGAGCCAGGACCTGAGACCCACCTACCCACCAGCTTCACTCTGAAGAGTTTC
AAGAGCCGGGCAACTGTGAGACCAAGGGGACTCCATCCTGGACTGCGTGCCCAAGGACGGGCAGAGC
CACTGCTGCATCCACGCAACACCTGTGTGTACCAGAATATGGGCATCTGGTGCAGGCAGAGAAT
GCGCTGGGACCAGCATGCCCCACAACCTGTGTCTTGATCCCATGGATGTTGTGAACTGGAGCCCCC
ATGCTGCGGACCATGGACCCAGCCTGAAGCGGCCCTCCCCAGGCAGGCTGCCTACAGCTGTGCTGG
GAGCCATGGCAGCCAGGCTGCACATAAATCAGAAGTGTGAGCTGCGCCACAAGCCGAGCGTGGAGAA
GCCAGCTGGGCACTGGTGGGCCCTCCCTTGGAGGCCCTTCACTATGAGCTCTGCGGGCTCCTCCA
GCCACGGCTACACCCTGCAGATACGCTGCATCCGCTGGCCCCTGCCTGGCCACTGGAGCGACTGGAGC
CCCAGCCTGGAGCTGAGAATAACGAAACGGGCCCCACTGTGAGACTGGACACATGGTGGCGGCAGAGG
CAGCTGGACCCAGGACAGTGCAGCTGTTCTGGAAGCCAGTGCCTGGAGGAAGACAGCGGACGGATC
CAAGGTTATGTGGTTTCTTGGAGACCCTCAGGCCAGGCTGGGGCCATCCTGCCCTCTGCAACACCACA
GAGCTCAGCTGCACCTCCACCTGCCTTCAAGAGCCAGGAGTGGCCCTGTGGCCTATAACTCAGCC
GGGACCTCTCGTCCACTCCGGTGGTCTTCTCAGAAAGCAGAGGCCAGCTCTGACCAGACTCCATGCC
ATGGCCCCGAGACCCTCACAGCCTCTGGGTAGGCTGGGAGCCCCCAATCCATGGCCTCAGGGCTATGTG
ATTGAGTGGGGCTGGGCCCCCCAGCGGAGCAATAGCAACAAGACCTGGAGGATGGAACAGAATGGG
AGAGCCACGGGGTTTCTGCTGAAGGAGAACATCAGGCCCTTTCAGCTCTATGAGATCATCGTACTCCC
TTGTACCAGGACACCATGGGACCCTCCCAGCATGTCTATGCCTACTCTCAAGAAATGGCTCCCTCCCAT
GCCCCAGAGCTGCATCTAAGCACATTGGCAAGACCTGGGCACAGCTGGAGTGGTGCCTGAGCCCCCT
GAGCTGGGGAAGAGCCCCCTTACCCACTACACCATCTTCTGGACCAACGCTCAGAACCAGTCTTCTCC
GCCATCCTGAATGCCTCCTCCCGTGGCTTTGCTCCTCCATGGCCTGGAGCCCGCCAGTCTGTATCACATC
CACCTCATGGCTGCCAGCCAGGCTGGGGCCACCAACAGTACAGTCTCACCTGATGACCTTGACCCCA
GAGGGTTCGGAGCTACACATCATCTGGGCTGTTGCGCCTCCTGCTGTTGCTCACCTGCCTCTGTGGA
ACTGCCTGGCTCTGTTGCAGCCCAACAGGAAGAATCCCCTCTGGCCAAGTGTCCAGACCCAGCTCAC
AGCAGCCTGGGCTCCTGGGTGCCACAATCATGGAGGAGGATGCCTTCCAGCTGCCCGGCCCTTGGCAGC
CCACCCATCACCAAGCTCACAGTGTGGAGGAGGATGAAAAGAAGCCGGTGCCTGGGAGTCCCATAAC
AGCTCAGAGACCTGTGGCCTCCCCACTCTGGTCCAGACCTATGTGCTCCAGGGGGACCCAAGAGCAGTT
TCCACCCAGCCCCAATCCCAGTCTGGCACCAGCGATCAGGCTGGGCTCCAGGGGATCTGCATACTTT
AAGGACCAGATCATGCTCCATCCAGCCCCACCAATGGCCTTTGTGCTTGTTCCTATAACTTCAGTA
TTGTAA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_172313
Insert Size: 2352 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172313.1
RefSeq Size:	2583 bp
RefSeq ORF:	2352 bp
Locus ID:	1441
UniProt ID:	Q99062
Cytogenetics:	1p34.3
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway, Pathways in cancer
MW:	86.6 kDa
Gene Summary:	<p>The protein encoded by this gene is the receptor for colony stimulating factor 3, a cytokine that controls the production, differentiation, and function of granulocytes. The encoded protein, which is a member of the family of cytokine receptors, may also function in some cell surface adhesion or recognition processes. Alternatively spliced transcript variants have been described. Mutations in this gene are a cause of Kostmann syndrome, also known as severe congenital neutropenia. [provided by RefSeq, Aug 2010]</p> <p>Transcript Variant: This variant (4) lacks an alternate segment compared to variant 1, that causes a frameshift. The resulting isoform (d) is shorter and has a distinct C-terminus compared to isoform a.</p>

Product images:



Circular map for SC314687