

Product datasheet for **SC314731**

PCDHB8 (NM_019120) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCDHB8 (NM_019120) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCDHB8
Synonyms:	PCDH-BETA8; PCDH3I
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC314731 representing NM_019120.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGCCAGCGGGAAGCTCATTTCAGACAAAGGCAAGTCCTTTTTCTTTCTCTTTTGGGCTTA
TCTCTGGCGGGCGCGCGGAACCTAGAAGCTATTCTGTGGTGGAGGAACTGAGGGCAGCTCCTTTGTC
ACCAATTTAGCAAAGGACCTGGGTCTGGAGCAGAGGGAATTCTCCAGCGGGGGTTAGGGTTGTTTCC
AGAGGGAACAAACTACATTTGCAGCTCAATCAGGAGACCGCGATTGTTGCTAAATGAGAAATGGAC
CGTGAGGATCTGTGCGGTCACACAGAGCCCTGTGTGCTACGTTTCCAAGTGTGCTAGAGAGTCCCTTC
GAGTTTTTTCAAGCTGAGCTGCAAGTAATAGACATAAACGACCACTCTCCAGTATTTCTGGACAAACAA
ATGTTGGTGAAGTATCAGAGAGCAGTCTCTGGGACTGCGTTTCTCTGAAGAATGCTGAAGACTTA
GATATAGGCCAAAACAATATTGAGAACTATAAATCAGCCCCAACTCTATTTTCGGGTCTCACCCGC
AAACGCAGTGATGCGAGGAAATACCCAGAGCTGGTGTGACAAAGCGCTGGACCGAGAGGAAGAAGCT
GAGCTCAGGTTAACTCACAGCACTGGATGGTGGCTCTCGGCCAGATCTGGCACTGCTCAGGTCTAC
ATTGAAGTTGTCGATGTCAATGATAATGCCCTGAATTTGAGCAGCCTTCTATAGGGTGCAGATCTCT
GAGGACAGTCCAATAAGCTTCTGGTGTGAAGGTCTCTGCCACGGATGTAGACACAGGAGTCAACGGA
GAGATTTCTATTCACTTTTCCAAGCTTCAGATGAGATAAGCAAAACTTTAAGGTGATTTCTTGACA
GGAGAAATTCGACTAAAGAAACAACCTGATTTGAAAAATTTTCAGTCTATGAAGTCAATATCGAGGCG
AGAGATGCTGGAGGCTTTTCTGAAAAATGCACCGTTCTGATTCAAGTGATAGATGTGAACGACCATGCC
CCAGAAGTTACCATGTCTGCATTTACCAGCCCAATACCTGAGAATGCGCCTGAACTGTGGTTGCACTT
TTCAGTGTTCAGACCTGATTCAGGAGAAAAATGGGAAAAAAGTTGCTCCATTCAGGAGGATCTACCC
TTCTCTGAAATCTTCTGTGGGAACTTTACACCCTACTAACAGAGACCACTAGACAGAGAAAAAGC
AGAGCCGAGTACAACGTCACTATCACCGTCACTGACTTAGGGACACCCAGGCTGACAACACATCTCAAT
ATGACCGTGTGGTGTGCGGACGTCATGACAACGCCCCCGCTTACCCAAACCTCTACACCCTGTTC
GTCCGCGAGAAACAACAGCCCCGCTGCACATCGGACGCGTACGCGCCACAGACAGAGACTCGGGCACC
AACGCCAGGTCACTACTCGTGTGCTGCCGCCAGGATCCGCACCTGCCCTCGCTCCCTGGTCTCC
ATCAACACAGACAACGGCCACCTGTTGCCCCTCAGGTGCTGACTACGAGGCCCTGCAGGCGTTCGAG
TTCCGGGTGGGCGCTTACAGACCGGCTCCCGGCTTTGAGCAGCGAGGCGCTGGTGCAGGCTGGTG
CTGGACGCCAACGACAACCTCGCCCTCGTGTGTACCCGCTGCAGAATGGCTCCGCGCCCTGCACCGAG
CTGGTCCCCGGGCGGCCGAGCCGGCTACCTGGTGACCAAGGTGGTGGCGGTGGACGGCGACTCGGGC
CAGAACGCTGGCTGTCTGTAACGAGTGTCAAGGCCACGGAGCCCGGGCTGTTCCGGTGTGTGGGCGCAC
AATGGCGAGGTGCGCACCGCCAGGCTGTGAGCGAGCGGACGCGGCCAAGCAGAGGCTGGTGGTGTG
GTCAAGGACAATGGCGAGCCTCCGTGCTCGGCCACCGCCACGCTGCACGTGCTCCTGGTGGACGGCTTC
TCCCAGCCCTACCTGCCGCTTCCGGAGGCTGCCAGCCAGGGCCAGGGCCGACTCTCTACCCGTCTAC
CTGGTGGTGGCGTTGGCCTCGGTGTCTTCGCTCTTCTCTTCTCGGTGCTCCTGTTCTGTGGCGGTGCTG
CTGTGTAGGAGGAGCAGGGCGGCTCGGTGGTGTGCTGCTCAGTGCCTGAGGGCCCTTTCCAGGGCAT
CTGGTGGACGTGAGGGGACCGGGAGCCTGTCTCAGAATATCAGTACGAGGTGTGCTGGCAGGAGGC
TCAGGGACGAATGAGTTCCAGCTCCTGAAACAGTATTACCTAATATTCAGGGCCATTCTTTGGGCCA
GAAATGGAACAAAACCTAACTTTAGGAATGGCTTTGGTTTCAGCCTTCAGTTAAAGTAA
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_019120

Insert Size: 2406 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019120.4](#)

RefSeq Size: 2711 bp

RefSeq ORF: 2406 bp

Locus ID: 56128

UniProt ID: [Q9UN66](#)

Cytogenetics: 5q31.3

Domains: CA

Protein Families: Transmembrane

MW: 87.6 kDa

Gene Summary: This gene is a member of the protocadherin beta gene cluster, one of three related gene clusters tandemly linked on chromosome five. The gene clusters demonstrate an unusual genomic organization similar to that of B-cell and T-cell receptor gene clusters. The beta cluster contains 16 genes and 3 pseudogenes, each encoding 6 extracellular cadherin domains and a cytoplasmic tail that deviates from others in the cadherin superfamily. The extracellular domains interact in a homophilic manner to specify differential cell-cell connections. Unlike the alpha and gamma clusters, the transcripts from these genes are made up of only one large exon, not sharing common 3' exons as expected. These neural cadherin-like cell adhesion proteins are integral plasma membrane proteins. Their specific functions are unknown but they most likely play a critical role in the establishment and function of specific cell-cell neural connections. [provided by RefSeq, Jul 2008]