

Product datasheet for **SC314740**

ZFX (NM_003410) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZFX (NM_003410) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZFX
Synonyms:	ZNF926
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC314740 representing NM_003410.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATGAAGATGGGCTTGAATTACAACAAGAGCCAAACTCATTTTTTGTGCAACAGGAGCTGATGGT
ACACACATGGATGGTGATCAAATTTGTTGTGGAAGTACAAGAAACTGTTTTGTTTCAGATGTTGTGGAT
TCAGACATAACTGTGCATAACTTTGTTCTGATGACCCAGATTCAAGTTGTAAATCCAAGATGTTATTGAG
GACGTTGTTATAGAAGATGTTCACTGCCCAGATATCATGGAAGAAGCAGATGTGTCTGAAACGGTCAATC
ATTCTGAGCAAGTGCTGGACTCAGATGTAAGTGAAGAAGTTCTTTAGCACATTGCACAGTCCCAGAT
GATGTTTTAGCTTCTGACATTACTTCAGCCTCAATGTCTATGCCAGAACACGTCTTGACGGGTGATTCT
ATACATGTGTCTGACGTTGGACATGTTGGACATGTTGGACATGTTGAACATGTGGTTCATGATAGTGA
GTGGAAGCAGAAATTGCTACTGATCCTCTGACTACCGACGTAGTTTCAGAAGAAGTATTGGTAGCAGAC
TGTGCCTCTGAAGCAGTCATAGATGCCAATGGGATCCCTGTGGACCAGCAGGATGATGACAAAGGCAAC
TGTGAGGACTACCTTATGATTTCTTGGATGATGCTGGCAAAATAGAACACGATGGTTCCTCTGGAATG
ACCATGGACACAGAGTCGGAAATGATCCTTGTAAAGTGGATGGCACTTGCCTGAGGTCATCAAGGTG
TACATTTTTAAAGCTGACCCCTGGAGAAGATGACTTAGGTGGAAGTGTAGACATTGTGGAGAGTGAGCCT
GAGAATGATCATGGAGTTGAACTGCTTGATCAGAACAGCAGTATTCGTGTTCCAGGGAAAAAGATGGTT
TATATGACTGTCAATGACTCTCAGCCAGAAGATGAAGATTTAAATGTTGCTGAAATCGCTGACGAAGTT
TATATGGAAGTGATCGTAGGAGAGGAGGATGCTGCAGCAGCAGCGCCAGCCGCCCGTGCACGAGCAG
CAAATGGATGACAATGAAATCAAACCTTCATGCCGATTGCATGGGCAGCAGCTTATGGTAATAATTCT
GATGGAAATGAAAACCGAATGGCACTGCAAGTGCCCTCTTGCACATAGATGAGTCTGCTGGCCTCGGC
AGACTGGCTAAACAAAAACCAAGAAAAGGAGAAGACTGATTCAGGCAGTACCAACAGCAATAATT
ATTGGCCCTGATGGACATCCTTTGACTGTCTATCCTTGCATGATTTGTGGGAAGAAGTTTAAAGTCGAGA
GGTTTTTTGAAAAGGCACATGAAAAACCATCCCGAACCTTGCCAAGAAGAAATACCGCTGTACTGAC
TGTGATTACACTACCAACAAGAAGATAAGTTTACACAACCACCTGGAGAGCCACAAGCTGACCAGCAAG
GCAGAGAAGGCCATTGAATGCGATGAGTGTGGGAAGCATTTCTCTCATGCAGGGGCTTTGTTTACTCAC
AAAATGGTGCATAAGGAAAAAGGAGCCAACAAAATGCACAAGTGTAAATCTGTGAATACGAGACAGCT
GAACAAGGGTTATTGAATCGCCACCTCTTGGCAGTCCACAGCAAGAACTTCTCTCATATTTGTGTGGAG
TGTGGTAAGGGTTTTCTCACCCGTCAGAGCTCAAAAAGCACATGAGAATCCATACTGGGAGAAGCCG
TACCAATGCCAGTACTGCGAATATAGGTCTGCAGACTCTTCTAACTTGAAAACGCATGTCAAACCTAAG
CATAGTAAAGAGATGCCATTCAAGTGTGACATTTGTCTTCTGACTTTCTCGGATACCAAAGAGGTGCAG
CAACATGCTCTTATCCACCAAGAAAGCAAAACACACCAGTGTTCATTGCGACCACAAGAGTTCGAAC
TCAAGTGATTTGAAACGACACATAATTTCAAGTTCACACGAAAGACTACCCCATAAAGTGTGACATGTGT
GATAAAGGCTTTTACAGGCCTTCAGAACTCAAGAAACACGTGGCTGCCCAAGGGCAAAAAATGCAC
CAGTGTAGACATTGTGACTTTAAGATTGCAGATCCATTTGTTCTAAGTCGCCATATTCTCTCAGTTCAC
ACAAAGGATCTTCCATTTAGGTGCAAGAGATGTAGAAAGGGATTTAGGCAACAGAGTGAGCTTAAAAAG
CATATGAAGACACACAGTGGCAGGAAAGTGTATCAGTGTGAGTACTGTGAGTATAGCACTACAGATGCC
TCAGGCTTTAAACGGCAGGTTATTTCCATTCACACGAAAGACTATCCTCACCGGTGTGAGTACTGCAAG
AAAGGCTTCCGAAGACCTTCAGAAAAGAACCAGCACATAATGCGACATCATAAAGAAGTTGGCCTGCC
TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_003410

Insert Size: 2418 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_003410.3
RefSeq Size:	7604 bp
RefSeq ORF:	2418 bp
Locus ID:	7543
UniProt ID:	P17010
Cytogenetics:	Xp22.11
Domains:	Zfx_Zfy_act, zf-C2H2
Protein Families:	Transcription Factors
MW:	90.5 kDa
Gene Summary:	<p>This gene on the X chromosome is structurally similar to a related gene on the Y chromosome. It encodes a member of the krueppel C2H2-type zinc-finger protein family. The full-length protein contains an acidic transcriptional activation domain (AD), a nuclear localization sequence (NLS) and a DNA binding domain (DBD) consisting of 13 C2H2-type zinc fingers. Studies in mouse embryonic and adult hematopoietic stem cells showed that this gene was required as a transcriptional regulator for self-renewal of both stem cell types, but it was dispensable for growth and differentiation of their progeny. Multiple alternatively spliced transcript variants encoding different isoforms have been identified, but the full-length nature of some variants has not been determined. [provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (1) and variants 2 and 3 encode the same isoform (1).</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>