

## Product datasheet for **SC314763**

### **PAXBP1 (NM\_013329) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PAXBP1 (NM_013329) Human Untagged Clone
Tag:	Tag Free
Symbol:	PAXBP1
Synonyms:	BM020; C21orf66; FSAP105; GCFC; GCFC1
Vector:	<u><a href="#">pCMV6 series</a></u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_013329, the custom clone sequence may differ by one or more nucleotides

```

ATGTTCCGAAAGGCCCGGGTGAACGTGCGCAAGCGGAACGACTCCGAAGAGGAAGAG
CGGGAACGCGATGAGGAGCAGGAGCCGCCCGTTGTTGCCGCCCGGGCAGGGCGAA
GAGGCGGGCCCCGGTGGCGGCGACAGGGCCCTGGCGGGAGTCGCTGCTGGCCCCGGG
CCGTCGCCGCCTCCGCGCTGACCCCGGCCTCGGGCTGAGGCCGGGGCGGCTTCCCC
GGCGGCGCGAGCCCGCAACGGGCTGAAGCCGCAAGAGGCCCTCGCGAGAACAAGAG
GTGCCCCGGGCCAGCCTGCTCAGCTTCCAGGACGAGGAGGAAGAAAATGAAGAAGTTTC
AAAGTGAAGAAATCAAGTTATAGCAAAAAGATAGTAAAATTGCTCAAGAAGGAATATAAA
GAAGATCTTGAAAAATCGAAGATTAAGACAGAACTCAACTCATCAGTGAAGTGAACAA
CCTTTGGACAAAACAGGACATGTTAAGGATACAAATCAAGAAGATGGAGTTATCATCAGT
GAACATGGTGAAGATGAAATGGATATGGAAGTGAAGAAAGAGGAAGAAAAGCCAAAGACT
GGTGGAGCTTTTTCAAATGCTTTATCTTATTGAATGTTCTTCGTCAGGAGAAATCCA
GATGCAGCTTTTATACATGCTGCAAGGAAAAGCGCCAAATGGCCCGAGAATTGGGAGAT
TTCCTCTCATGATAATGAGCCTGGTAAAGGCCCTTGTAGAGAAGATGAGAATGAT
GCCAGTGATGATGAAGATGACGATGAGAAAACGCCGATAGTTTTTCTGTGAAAAGAAAAG
TCACAAAGACAAAAATTGCTGAGGAAATAGGAATTGAGGGGAGTGATGATGATGCTTTA
GTAAGTGGAGAACAGGATGAAGAGCTCAGCCGATGGGAACAGGAGCAGATAAGGAAAAGGA
ATTAATATCCCTCAGGTTCAAGCCAGTCAACCCGAGAAGTGAATATGTACTACCAGAAC
ACTTACCAGACAATGCCTTACGGCTCATCTATGGCATTCTTATAGTTATACGGCTAT
GGATCATCAGATGCCAAATCTCAAAAAACAGATAATACAGTCCCTTTCAAACTCCCAGT
AATGAGATGACTCCCGTTACTATTGATTTGGTAAAGAAACAGCTTAAAGACAGGTTGGAC
TCCATGAAAGAATTGCACAAAACAAATCGACAGCAGCATGAGAAACATCTGCAAAGCCGA
GTGGACTCTACCAGGGCTATTGAAAGATTAGAAGGGTCTTCTGGGGGTATTGGTGAACGG
TATAAATTTTTGCAAGAAATGCGAGGGTATGTCGAAGACTTGCTTGAGTGTTCAGTGAA
AAGGTGCCACTGATTAATGAAGTGAATCAGCAATACATCAGCTGTACAAAACAGCGAGCT
TCCCGCTTGTCCAAAGACGACAAGATGATATTAAGATGAATCTTCGGAGTTTTCAAGC
CATTCAAACAAAGCTCTGATGGCACCAATCTTGACTCTTTGGACGCGATCGGGCACTG
TATCAAGAGCATGCAAAACGTCGATTGCAGAGCGGGAGGCCAGGAGACTCGTCGTAGA
CAAGCCAGAGAACAACCGGTAAGATGGCAGATCACCTTGAAGGCCTTTCAGTGATGAT
GAAGAACTTCTACAGATATTACTAATTTCAATCTGAAAAAGATCGAATTTCAAAGAA
TCCGGCAAAGTTTTGAAGATGTCTTGAAGTTTCTATTCAATTGACTGTATTAATCA
CAGTTTGAAGCATGGCGTTCAAAATACTACACATCCTACAAAGATGCTTACATTGGCCTT
TGTTTTGCCAAAATTATTCAACCCCTCATACGACTTCAGCTCCTCACTTGGACTCCTCTT
GAGGCAAAATGTCGTGACTTTGAGAATATGCTGTGGTTTGAATCTTTGCTGTTTTATGGT
TGTGAAGAACGAGAGCAAGAAAAGATGATGTAGATGTTGCCCTACTACCTACCATTGTG
GAAAAGGTGATTCTTCTAACTAACAGTGATAGCTGAAAATATGTGGGACCCCTTTTCT
ACAACACAGACTTCAAGAATGGTGGGAATTACACTAAAATTAATCAATGGATATCCTTCA
GTAGTGAATGCAGAAAATAAAAATACACAGGTATACCTAAAGGCACTTTTATTGAGAAT
AGAAGAACTTTAGATGATGATGATTTTATGCCCTTATATCCCAAAAATGTCTTAGAAAAT
AAAAATCTGGGCCTTACTTGTTTTTTCAACGACAGTTTTGGTCTTCAAGTAAAGTCATA
AAACCCCATTCAGAGAGGGTCTGCCCATACCCAGAAGGAAAGAAATGCTGCTCAGAG
AGGCCAAGAAGAACTGGACGGACAGGCCTTGTGTGTCTTCTCT

```

**Restriction Sites:** Please inquire

**ACCN:** NM\_013329

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_013329.2</a></u> , <u><a href="#">NP_037461.2</a></u>
<b>RefSeq Size:</b>	2467 bp
<b>RefSeq ORF:</b>	2448 bp
<b>Locus ID:</b>	94104
<b>UniProt ID:</b>	<u><a href="#">Q9Y5B6</a></u>
<b>Cytogenetics:</b>	21q22.11
<b>Protein Families:</b>	Transcription Factors
<b>Gene Summary:</b>	<p>This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009]</p> <p>Transcript Variant: This variant (2) differs in the 3' coding region and 3' UTR, compared to variant 1. The encoded isoform (2) has a distinct C-terminus and is shorter than isoform 1.</p>