

## Product datasheet for **SC314765**

### Cadherin 24 (CDH24) (NM\_022478) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cadherin 24 (CDH24) (NM_022478) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDH24
Synonyms:	CDH11L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC314765 representing NM\_022478.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTGGGGCTGGTGGGCTCCTGCTGGCTGGCTGGTGGTGGGCTGCATGGGGCTCTGGCAGCC
CCAGCCCGGGCTGGGCAGGGTCCCGGGAACACCCAGGGCTGCTCTGCTGCGGACTCGAAGGAGCTGG
GTCTGGAACCAAGTCTTTTGTATTGAGGAATATGCTGGTCCAGAGCCTGTTCTCATTGGCAAGCTGCAC
TCGGATGTTGACCGGGGAGAGGGCCGACCAAGTACCTGTTGACCGGGGAGGGGGCAGGCACCGTATTT
GTGATTGATGAGGCCACAGGCAATATTCATGTTACCAAGAGCCTTGACCGGGAGGAAAAGGCGCAATAT
GTGCTACTGGCCCAAGCCGTGGACCGAGCCTCAAACCGGCCCTGGAGCCCCATCAGAGTTCATCATC
AAAGTGAAGACATCAACGACAATCCACCCATTTTCCCCTTGGGCCCTACCATGCCACCGTGCCCGAG
ATGTCCAATGTCGGGACATCAGTGATCCAGGTGACTGCTCAGGATGCTGATGACCCAGCTATGGGAAC
AGTGCCAAGCTGGTGTACTGTCTGGATGGACTGCCTTCTTCTCTGTGGACCCCGAGCTGGAGTG
GTGCGTACAGCCATCCCCAATGGACCGGGAGACACAGGAGGAGTTCTTGGTGGTGTCCAGGCCAAG
GACATGGCGGCCACATGGGGGGCTGTACGGCAGCACTACGGTACTGTACGCTCAGCGATGTCAAC
GACAACCCCCCAAGTCCCACAGAGCCTATACCAAGTTCCTCGTGGTGGAGACAGCTGGACCTGGCACA
CTGGTGGGCCGGCTCCGGGCCAGGACCCAGACCTGGGGGACAACGCCCTGATGGCATAACAGCATCTG
GATGGGGAGGGTCTGAGGCCTTCAGCATCAGCACAGACTTGACGGGTCGAGACGGGCTCCTCACTGTC
CGCAAGCCCTAGACTTTGAGAGCCAGCGCTCCTACTCCTCCGTGTCGAGGCCACCAACACGCTCATT
GACCCAGCCTATCTGCGGCGAGGGCCCTTCAAGGATGTGGCTCTGTGCGTGTGGCAGTGAAGATGCC
CCAGAGCCACCTGCCTTACCCAGGCTGCCTACCACCTGACAGTGCCTGAGAACAAGGCCCGGGGACC
CTGGTAGGCCAGATCTCCGCGGTGACCTGGACTCCCCTGCCAGCCCAATCAGATACTCCATCCTCCCC
CACTCAGATCCGGAGCGTTGCTTCTCTATCCAGCCCGAGGAAGGACCATCCATACAGCAGCACCCCTG
GATCGCGAGGCTCGCGCCTGGCACAACCTCACTGTGCTGGCTACAGAGCTCGGCTGGAGCTGGGGCCA
GAAAGGGGCTGGGTACCTCTTCTGGTTGCTGAGTGGAGTGCACCAGCAGCCCCACCCAGAGAAGCCCT
GTTGGAAGCGCTGTGGGAATCCCCAAGACAGTTCTGCACAGGCCTCGCGCGTGAAGTGGCCATCCAG
ACCCTGGATGAGAATGACAATGCTCCCCAGCTGGCTGAGCCCTACGATACTTTTGTGTGACTCTGCA
GCTCCTGGCCAGCTGATTAGGTGATCCGGGCCCTGGACAGAGATGAAGTTGGCAACAGTAGCCATGTC
TCCTTTCAAGGTCTCTGGGCCCTGATGCCAATTTACTGTCCAGGACAACCGAGATGGCTCCGCCAGC
CTGCTGCTGCCCTCCCGCCTGCTCCACCCCGCATGCCCTACTTGGTTCCCATAGAAGTGTGGGAC
TGGGGGACGCGCGCTGAGCAGCACTGCCACAGTACTGTTAGTGTGTGCCGCTGCCAGCTGACGGC
TCTGTGGCATCCTGCTGGCCTGAGGCTCACCTCTCAGTGTGGGCTCAGCACCGGCGCCCTGCTTGCC
ATCATCACCTGTGTGGGTGCCCTGCTTGCCCTGGTGGTGTCTTTCGTGGCCCTGCGGCGGAGAAGCAA
GAAGCACTGATGGTACTGGAGGAGGAGGACGTCCGAGAGAACAATCATCACCTACGACGACGAGGGCGGC
GGCGAGGAGGACACCGAGGCCTTCGACATCACGGCCTTGCAAGCCCGGACGGGGCGGCCCGGGCGG
CCCGGCCCTCCCGCGCCGAGACGTGTTGCCCGGGCCCGGGTGTGCGCCAGCCAGACCCCGGGC
CCCGCCGACGTGGCGCAGCTCCTGGCGCTGCGGCTCCGCGAGGCGGACGAGGACCCCGGCGTACCCCG
TACGACTCGGTGCAGGTGTACGGCTACGAGGGCCGCGGCTCCTTTCGGGCTCCCTCAGCTCCCTGGGC
TCCGGCAGCGAAGCCGGCGGCGCCCGGCCCGGAGCCGCTGGACGACTGGGTCCGCTCTTCCGC
ACCCTGGCCGAGCTGTATGGGGCAAGGAGCCCCGCGCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_022478

**Insert Size:** 2460 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_022478.3](#)

**RefSeq Size:** 3567 bp

**RefSeq ORF:** 2460 bp

**Locus ID:** 64403

**UniProt ID:** [Q86UP0](#)

**Cytogenetics:** 14q11.2

<b>Domains:</b>	Cadherin_C_term
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	87.8 kDa
<b>Gene Summary:</b>	<p>Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediate strong cell-cell adhesion.</p> <p>[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).</p>