

Product datasheet for **SC314907**

anoctamin 5 (ANO5) (NM_213599) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	anoctamin 5 (ANO5) (NM_213599) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANO5
Synonyms:	GDD1; LGMD2L; LGMDR12; TMEM16E
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC314907 representing NM_213599.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGCGACCCGGATCTCCTGGAAGTGTGGCGGAGGAAGGGGAAAAAGTCAATAAGCATATAGACTAC
CTTTCCAAATGAGTGAGCAGAGCCTGAGCAGCAGAGAGACCAGCTTTTCATCAATGAAGAAACAATG
CCTGCAAAGCGATTCAATTTGTTCTGAGGCGCGGCTTATGTTTCAAAAAATCAGCAAAGCAAAGAT
TCTATCTTCTCCGAGATGGGATTAGGCAAATTGATTTTGTGCTTTCCTACGTTGATGATGTAAGAAA
GACGCAGAGTTAAAGCGGAAAGAAGAAAAGAGTTGAACTAATCTCAGAAAAACAGGCTTGGAGTTG
GAAATAGAAGACAAAAGGGACTCGGAAGATGGAAGAACTATTTTGTCAAGATCCATGCCCTTGGGAG
GTATTAGTTACCTATGCTGAAGTCTTGGGAATCAAATGCCTATTAAGGAGAGTGATATCCCGCCCT
AAGCACACTCCTATAAGCTATGTGCTTGGACCTGAAGACTCCCACTGAGTGTGAAGTATCCCATCCT
GAATATTTACTGCACAATTCAGCAGACATCGGCAGGAGCTCTTCTCATCGAAGATCAGGCAACCTTC
TTCCATCCTCATCAAGAAACAGAAATTGTACTATATTTCTCAAGATGCCTTTTGGCATAGAAGAT
GGGAAGAAAAGGTTTGGGATTGAAAGACTGCTAAACTCTAACACTTACTCATCTGCTATCCACTCCAT
GATGGCCAAATATTGGAAGCCATCAGAACCTCCCAATCCTACCAATGAAAGATACACACTTACCAGAAT
TGGGCTCGATTTTCTATTTCTACAAGGAGCAGCCTTTAGACTTGATTAAGAATTATTATGGAGAAAA
ATTGGTATCTATTTTGTCTTTCTTGGATTTTACACAGAAATGCTATTCTTGCAGCTGTAGTTGGCTTA
GCTTGTATTTATTTATGGCTTATTATCAATGGAACATAACACAAGCAGCACTGAAATCTGTGACCCTGAG
ATTGGTGGTCAGATGATCATGTGCCACTCTGTGATCAAGTGTGTGATTATTGGAGACTAAATAGTACG
TGTTTGGCTTCAAAGTTCTCCATTTGTTTGATAATGAGTCAACAGTGTCTTTGCAATTCATGGGA
ATTTGGTCCACTTATTTTGGAGTTTGGAAACAACGACAAGCCAGACTGGAATATGAATGGGACCTG
GTGGACTTTGAAGAGGAACAGCAGCAGCTTCAGCTGAGACCAGAATTTGAAGCTATGTGTAACACAGG
AAATTGAATGCAGTGAATAAGGAGATGGAACCTTACATGCCTCTATACACGCGTATTCCATGGTACTTT
CTTTCAGGAGCCACAGTGACATTATGGATGTCTTGTGCTCACCAGTATGGTAGCTGTAATTGTGTAC
CGCCTGTCAGTCTTGTACATTTGCTAGTTTTCATGGAAGTGTGATGCATCCTTAAAGCAGGTCAAAGC
TTCCTTACTCCTCAGATAACCACATCACTCACAGGATCATGCTTGAACCTTATTGTCATCTTGATCTTG
AATTTCTTTTATGAAAAGATATCTGCCTGGATCACAAAAATGGAATTCCTCGAACATACCAGGAGTAT
GAGAGCAGTCTTACCTTGAAGTGTTCCTGTTTTCAGTTTGTAAATTTTACTCATCTGCTTCTACGTA
GCTTTCTTTAAAGGGAAGTTCGTAGGCTATCCTGGAAAAACACATATTTATTTAATGAGTGGAGAAGT
GAAGAGTGTGATCCTGGAGGCTGTCTTATAGAATTGACAACCCAATTGACCATTATAATGACCGGAAA
CAGATTTTGGAAACATTAAGAAGCCATTTATCCCTTGCTTTGAATTGGTGGAGACGCCGAAAAGCT
CGGACAACTCTGAGAAGCTGTATAGTCGATGGGAGCAGGATCATGACCTTGAAGTTTTGGACCCCTT
GGGCTTTTCTATGAGTACTTAGAAAACAGTTACTCAATTTGGATTTGTTACTACTATTTGTGGCCTCTTT
CCTTTGGCTCCTCTTCTGCTCTCATAAATAATATTGTAGAGATTCGAGTGGATGCCTGGAACTTACC
ACTCAATACAGGAGAAGTGTAGCTTCTAAAGCTCATAGCATAGGTGTTTGGCAAGACATTCTTTATGGA
ATGGCTGTCTTTCTGTGCAACTAATGCCTTTATTGTTGCATTTACGTACAGACATCATTCCCGTCTA
GTTTACTACTATGCTTACTCAACAAATGCCACACAGCCTATGACAGGATATGTGAATAATAGCCTGTCA
GTATTCTGATAGCTGATTTTCCAAACCACACTGCACCTTCGGAAAAACGAGACTTCATCACTTGCAGG
TACAGAGATTACAGATATCCTCCTGATGACGAGAATAAATATTTTCATAATATGCAATTCTGGCATGTC
CTTGCTGCCAAGATGACCTTTCATCATTGTTATGGAACATGTTGTGTTTTTGTAAATTTTGTGGCC
TGGATGATACCTGATGTTCCAAAAGATGTTGTGGAGAGAATCAAGAGAGAAAAGTTAATGACTATCAAG
ATTCTCCATGATTTTGTGCTCAACAAATTAAGAAGAACTTGGGAATTAATTCTAATGAATTTGCCAAG
CATGTCATGATTGAGGAAAACAAAGCACAGCTGGCTAAATCAACACTTAA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

Plasmid Map:



ACCN: NM_213599

Insert Size: 2742 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_213599.2](#)

RefSeq Size: 6661 bp

RefSeq ORF: 2742 bp

Locus ID: 203859

UniProt ID: [Q75V66](#)

Cytogenetics: 11p14.3

Protein Families: Transmembrane

MW: 107.2 kDa

Gene Summary: This gene encodes a member of the anoctamin family of transmembrane proteins. The encoded protein is likely a calcium activated chloride channel. Mutations in this gene have been associated with gnathodiaphyseal dysplasia. Alternatively spliced transcript variants have been described. [provided by RefSeq, Nov 2009]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).