

## Product datasheet for **SC314925**

### ZMIZ2 (NM\_031449) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZMIZ2 (NM_031449) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZMIZ2
Synonyms:	hZIMP7; NET27; TRAFIP20; ZIMP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC314925 representing NM\_031449.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAACTCCATGAACCCCATGAAACCTGCCCTGCCCCCTGCGCCACACGGTGATGGTTCATTGCGCATAT
GAGTCTGTGCCCTTGGCAACAAAGCGCCACTCAGCCGGCTGGATCGCTGTCTGTGGTCACTACTGTGTGG
GGAGTTGGCAACGCGACACAGAGCCAGGTTTTGGGGAACCCCATGGGCCCTGCAGGGAGTCCCTCTGGC
AGCTCCATGATGCCTGGTGTGGCAGGGGCGAGCTCCGCCTTGACCTCCCCACAGTGCCTGGGACAGCAG
GCGTTTGTGAAGGCGCGCAACAAGGGCTACGTGCAGCAAGGCGTGTACAGCCGCGGGGGCTACCCT
GGGGCCCCCGGCTTACCACCGGTATGCAGGCGGCCCGGGGGCTGGGCTCCCTCACATGTGCA
AGACCCTCACTGACTTCACGCAAGCGGCAGTGTGCAGCTGTGGTGTGCGGCAGCCACTGCCACC
GCCACAGCCACAGCCACCGTGGTGTCTCCTCAGGAGAAGCAGAGCCAGGAGCTGAGCCAGTATGGAGCG
ATGGGGGCGGACAGTCTTTAACAGCCAGTTTCTGCAGCATGGAGGTCCCCGGGGGCTAGTGTCCCC
GCTGGCATGAACCTACTGGCATAGGAGGGTAATGGGCCCTCTGGCCTCTCCCTTGGCTATGAAC
CCCACCCGGGACAGGAATGACACCCTTGTATGCAGGGCAGCGTCTGCCCAACATGGGTATCCTGGG
CCTCCCCAGGCCAGCCACTGCCCGACAGGGGTCAAGAGAACCTACTCTGAGGTGTATCCAGGGCAG
CAGTATCTGCAAGGAGGCCAGTATGACCCAGCACCGCCAGTTTGCGCCAGCCCTGGGACGCCCCCT
GCCCCCTCCCCTTCTACCCTGGGCACAGGCTGCCCTGCAGCAGGGCATGACCCAGTCCCTGTCCGTG
CCTGGCCCCACGGGACTGCATTAAAGCCACAGAGCAGTTCAACGGGCAGGGCGCCAGTTCACAGGG
GGCAGCGTCAGCTACAGCCAACCTGGCTGAGTGGGCTACCCGTTCCATCCCGGGCTATCCAGTTCC
CCACTGCCAGGGAACCCACGCCACCCATGACCCCAAGCAGCAGCGTCCCTTACATGTACCAAACCAA
GAGGTCAAGTCTCCCTTCTTGCCTGATCTCAAGCCCAACCTCAACTCCTGCACATCGCCCTTGGA
AGCGGGCTTGTGACGAGTTGCGGCTGACCTTCCCTGTGCGCGATGGGGTGGTCCCTGGAGCCCTCCGC
CTGCAGCAACAACCTGGTGTAAAGCAACCATGTCTTCCAGCTGCGAGACTCAGTCTACAAGACCCTGATA
ATGAGGCTGACCTGGAGCTGCAATTCAAGTGCTACCACCAGGAGACCGGCAGATGAACACCAACTGG
CCAGCCTCGGTGCAGGTCAAGTCAATGCCACGCGCTCACCATCGAGCGTGGGACAACAAGACCTCG
CACAAGCCACTCTACCTGAAGCATGTGTGCCAGCCAGGCCGAACACCATCCAGATCACCGTACCAGCC
TGCTGTGCTCCCACCTTCTGTGCTGCAGCTAGTGCACCGCCATCCGTCGGTGTGCAAGGGC
CTCCTCAAAAAGCGCCTCTGCCTGCTGAGCACTGCATACCAAAGATAAAGCGGAACTTACAGCAGGGC
ACCATCCCTGGCACCCCTGGGCCAACGGAGAGGACGGGTGGAGCAGACAGCTATCAAGGTGTCCCTG
AAGTGCCCATACCTTCCGAGGATCCAGTCCCTGCCCGAGGTGATGACTGTCCGCACATACAGTGC
TTTGACCTGGAGTTCGACCTGCAGCTCAACTGTGAGCGGGGACTTGGAGGTGTCCTGTGTGCAACAAG
ACAGCTTTGCTGGAGGGCTGGAGGTGGACCAGTACATGCTGGGCATCCTGATTTACATTAGAAGTCT
GACTATGAGGAGATCACCATCGACCCACGTGCAGCTGGAAGCCAGTGGCCGTGAAGCCTGACATGCAC
ATCAAGGAGGAGCCGGATGGGCCAGCACTGAAGCGCTGCCGCACCGTGCAGCCCGCCACGTGCTCATG
CCCAGCGTGATGGAGATGATCGCCGCCCTGGGCCCGGGCTGCCCTTGGCCCCCTGCAGCCCCC
TCAGTCCCTGCCCCAGGCACTACCCTGGCCAGGGTCCAGCTTCTGGGGCTGGAATTTCCCTGAG
TCCTTCCCACCCACCAGCCAGCACCCCAACCCTTGTGAGTTCACCCGGGACCACCCCATCTCC
TACCAGTCTGACATTCAGCAGCCTCCTGACTTCAGAGAAGTCTACCGCTGCCTCCCAAGCCAGATG
GCACCAGCAGGTACCTGGACCCCACTACAATCCTGGGACACCAGGACTACACACCTCCAACTTGGG
GCCCTCCAGGTCCCAGCTGCACCAATTCAAACCTCCCCAGCGTCCCGGCAGTCTTGGGCCAAGCG
AGCTTAGGACCTACGGGTGAAGTGGCTTCACTCCTGCCACAGGCGTGTGGGGCCCCCAGCATGTCT
GGAGCCGGGGAGGCCCAAGAACAGCTCTGGACCTGCTCCCGAACTGACCAACCCTGATGAGCTACTG
TCCTACTTGGGCCACCCGACCTCCCTACGAACAACAATGACGACCTGCTTCTGTGTTGAGAACAAC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_031449

**Insert Size:** 2763 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_031449.3](#)

**RefSeq Size:** 5159 bp

**RefSeq ORF:** 2763 bp

**Locus ID:** 83637

**UniProt ID:** [Q8NF64](#)

**Cytogenetics:** 7p13

**Domains:** zf-MIZ

**MW:** 96.5 kDa

**Gene Summary:** ZMIZ2 and ZMIZ1 (MIM 607159) are members of a PIAS (see MIM 603566)-like family of proteins that interact with nuclear hormone receptors. ZMIZ2 interacts with androgen receptor (AR; MIM 313700) and enhances AR-mediated transcription (Huang et al., 2005 [PubMed 16051670]).[supplied by OMIM, May 2010]  
Transcript Variant: This variant (1) represents the longest transcript and it encodes the longest isoform (1).