

Product datasheet for **SC315088**

ZNF197 (NM_006991) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF197 (NM_006991) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF197
Synonyms: D3S1363E; P18; VHLaK; ZKSCAN9; ZNF20; ZNF166; ZSCAN41
Vector: pCMV6 series
Fully Sequenced ORF: >NCBI ORF sequence for NM_006991, the custom clone sequence may differ by one or more nucleotides

```

ATGACAAGAGAAAATGTAGCCCAATGCTCTGAGACAAGAGGGCCTTGTGAAGGGGAAG
GATGATACCTGGAAATGGGGAACCACTTCCAAGGAAGTAGCTCCTCTGTTGGGAGACC
TCCACCTACTTTAGACAATTACGTTACCATGAGACATCTGGACCCAGGAAGCCCTG
AGCCGGCTCAGGGAAGCTGTGCGCGGTGGCTGAGACCAGAAGCACGCACCAAGGCACAG
ATCCTGGAGCTGCTGGTCTGGAGCAGTTTCTGAGCATCCTGCCTGGGGAGATTCGGACC
TGGGTACAGCTCCATCACCTGGAAGTGGCGAGGAGGCTGTGGCCCTGGTAGAGGAGCTG
CAGAAAGACCTTGATGGACCAGCAATACAAGTTCCAGTCTTGCAAGGATCAGGACACT
CTCCAGAAGGTGGTGTGAGTGCCCAAGGAACAACACTTCTCCTGTACTTCTGGCAGCCAC
ATAGCAGCTGAAATTTGCCCGCATCCTCCTACTGACCTAGTGGCATTCAACCTCCAGGAT
CCTCAGCATGATTCTCCTGCCCTGAAGCTTCTGCCCTTTCCAGGAAGAGAACCCAAGA
AATCAATTAATGGCACTTATGCTCCTAACAGCCAGCCCAAGGAGTTGGTGTGTTTCGAG
GAGGTGTGAGTATGCTTCACTTCAAGGAATGGGCATGTCTGGGCCCAATCCAGAGGGCC
TTGTACTGGGATGTGATGCTGGAGAATTATGGAAATGTGACCTCCCTAGAATGGGAGACC
ATGACCGAGAATGAGGAGGTGACATCAAAGCCAAGTAGTTCTCAAAGAGCAGACTCTCAT
AAAGGAACATCAAAAAGACTTCAAGGAAGTGTCCCAAGTCCCTGATTTTGAAGAAGAG
TGTGAATGGCAAGTTTTGGCAAGTCAAGTGGGAAATGAAACAGATGAAAGGGCAGATACA
GTGAAGAAAGTTTCCCTTTGTGAACGAGACAAGAAGAAAAGGACTCCACCAGAGAAAACA
GGCCAAAAGTGGAAAGGAATTAGGAGACAGCTTGACTTTTCGGTTACAGTATTTCTGAAAGT
TTAATAGTACTGAAGGAAAGAAGTTTTATAAATGTGATATGTGTTGTAACATTTTAAAT
AAAATCTCCATCTTATAAACCATCGGAGAATCCACACTGGTGAGAAACCTCATAAATGT
AAGGAATGTGGAAAAGGCTTTATTCAGCGTTTCGAGCCTTCTAATGCATTTACGGAACCAT
TCAGGGGAGAAAACCTTATAAATGTAATGAATGTGGAAAGCATTCTCAAAGTGTCTTAC
CTTCTAAACCATCAGAGGATCCACACTGGGGAGAAAACCTTATAAGTGAAGGAGTGTGGA
AAGGGCTTCTATAGGCACTCAGGCCTAATTATACATCTAAGGCGCCATTCAGGGGAGAGA
CCTTATAAGTGAATGAATGTGGAAAAGTCTTCTCAGAAATGCTTACCTCATTGACCAT
CAGAGGCTCCACAAAGGGGAAGAACCTTATAAATGTAATAAGTGTGAGAAAGCTTTCATT
CTGAAGAAGAGCCTCATTCTGCACCAGAGAATCCACTCTGGGGAAAAACCTTATAAATGT
GATGAATGTGGAAAAGACCTTTGCTCAGACCCTTATCTTATTGACCATCAGCGACTCCAC
AGTGCAGAGAACCCTTACAAGTGTAAAGAATGTGGAAAAGTTTTTCATTGAAAGCAAAAGC

```



[View online »](#)

```

CTCCTCTTACATCAGAGAGTCCACACAGAAAAGAAAACCTTTGGTTGTA AAAAGTGTGGG
AAGATTTTCAGTTCTAAGTCAAACCTTCATTGACCATAAGAGGATGCACAGCAGAGAGAAA
CCTTACAAATGCACTGAATGTGGGAAAGCCTTTACTCAAAGTGCTTACCTTTTTGACCAC
CAGAGACTCCACAATGGGGAGAAGCCCTATGAATGTAATGAATGTGGGAAAGTTTTTATT
CTGAAGAAGAGCCTCATTTCATCAAAGGTTCCACACTGGAGAGAATCTCTATGAATGT
AAAGATTGTGGTAAGGTCTTCGGTTCAAACAGAAAACCTCATTGACCATGAGAGACTCCAC
AATGGGGAGAAGCCATATGAATGTCGAGAGTGTGGGAAAACCTTTATTATGAGCAAAAAGT
TTTATGGTCCATCAGAACTCCATACACAAGAGAAAGCCTACAAATGTGAGGATTGTGGG
AAGGCTTTCAGTTACAATTCAAGCCTGCTTGTACATCGGAGAATCCACACCGGAGAAAAA
CCCTTTGAATGCAGTGAGTGTGGAAGAGCTTTCAGTTCAAACAGAAAACCTCATTGAGCAT
AAGAGAAATCCACAGTGGTGAGAAACCTATGAGTGTGATGAGTGTGGCAAATGCTTCATT
CTGAAGAAAAGCCTCATTGGACATCAGAGAATTCACACGAGGGAAAAATCTTATAAATGC
AATGACTGTGGGAAGGTCTTCAGTTACCGCTCAAACCTTATAGCCCATCAGAGGATCCAT
ACTGGTGAGAAGCCCTATGCGTGTAGTGAGTGTGAAAAAGTTTTACGTACAACAGAAAAC
CTGATTGAACATCAAAGAATTCACAGTGGAGAAAAACCTACGAATGTCATGTATGTAGG
AAAGTCCTTACCTCTAGTAGAAATCTTATGGTACATCAAAGAATCCATACTGGAGAGAAA
CCTTATAAATGTAATGAGTGTGGAAAAGACTTTAGTCAGAATAAAAAACCTTGTGTACAT
CAGAGAATGCACACTGGGGAAAAACCTTATGAGTGTGACAAGTGTAGGAAATCCTTTACT
TCTAAGAGGAATTTAGTTGGCCACCAGAGAATTCACACAGGGGAGAAACCTATGGGTGT
AATGATTGTAGTAAAGTTTTAGGCAAAGAAAAACCTTACTGTACATCAGAAAAATCCAC
ACAGATGAAAAACCTTGTGAATGTGATGTGTCTGAAAAAGAATTCTCTCAGACTTCCAAC
CTTCATCTTCAACAGAAAATCCATACCATTGAGGAATTCTCTTGGCTACAAAACACCAAT
GAGTCCAAGATTGAGATTCAGAAAATC

```

Restriction Sites:

Please inquire

ACCN:

NM_006991

OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:

[NM_006991.3](#), [NP_008922.1](#)

RefSeq Size:

4338 bp

RefSeq ORF:

3090 bp

Locus ID:	10168
UniProt ID:	O14709
Cytogenetics:	3p21.31
Domains:	KRAB, LER, zf-C2H2
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	<p>This gene product belongs to the zinc finger protein superfamily, members of which are regulatory proteins characterized by nucleic acid-binding zinc finger domains. The encoded protein contains 20 tandemly arrayed C2H2-type zinc fingers, a Kruppel-associated box (KRAB) domain, and a SCAN box. This transcript turns over rapidly and contains 3' UTR AUUUA motifs, which are often a hallmark of rapid turnover. It is overexpressed in some thyroid papillary carcinomas. This gene is located in a cluster of zinc finger genes at 3p21. Naturally-occurring readthrough transcription is observed between this gene and the upstream zinc finger protein 660 gene and is represented by GeneID:110354863. [provided by RefSeq, May 2017]</p> <p>Transcript Variant: Variants 1 and 3 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>