

Product datasheet for **SC315090**

MMS19 (NM_022362) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: MMS19 (NM_022362) Human Untagged Clone
Tag: Tag Free
Symbol: MMS19
Synonyms: CIAO4; hMMS19; MET18; MMS19L
Mammalian Cell Selection: None
Vector: pCMV6-XL4
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_022362 edited
GAGGCGCATGCGCGGAGAGCCCGTCTCGAGCCACCGCCCATATCCCTCCCACGGTCTCT
AGTTCGCGTTATGGCCGCTGCCCGGCTGTGGAGCGCGCGCCTATGGGTGCCCTATG
GGGCTCGTGCACGACTTCGTCTGGGTGAGCAAGAGGGCCCGCTGACCAGGTGGCTGC
AGATGTGAAATCTGGCACTATACAGTGTACAAGTTGTGGAAGCCCTTGGTCTCTCT
AGAGAATCCAGAACCCGAACTCGGCACGAGGAATCCAGCTTTTGTACAGGTGCTACT
CCACTGTCACACCTTGCTCCTGGAGAAGGAAGTGTACACCTGATACTGTTCTATGAGAA
CCGGCTGAAGGACCATCATCTTGTGATCCCATCTGCTCCTGCAGGGTTTGAAGGCACTTAG
CCTGTGTGGCCCTGCCCCAGGGCTGGCTGTTTCTGTGCTTAAAGCCATCTTCCAGGA
AGTGCATGTACAGTCCCTGCCACAGGTGGACCGACACACAGTCTACAATATCATACCAA
TTTTATGCGAACC CGGAAGAAGAGCTAAAGAGCCTAGGAGCTGACTTACCTTTGGCTT
CATCCAGGTGATGGATGGGAAAAGGATCCCCGTAACTTCTGGTGGCCTTCCGCATCGT
CCATGACCTCATCTCCAGGGACTATAGCCTGGGACCCCTTGTGGAGGAGTTGTTTGAAGT
GACATCCTGTTATTTCCCTATCGATTTTACCCCTCCACCTAATGATCCCCATGGTATCCA
GAGAGAAGACCTCATCTGAGTCTTCGCGCTGTGCTGGCTTCTACACCACGATTTGCTGA
GTTTCTGCTGCCCTGTTGATTGAGAAAGTGGATTCTGAGGTTCTGAGTGCCAAGTTGGA
TTCTCTACAGACTCTGAATGCTTGTGCTGTGTATGGACAGAAGGAAGTGAAGGACTT
CCTCCCCAGCCTTTGGGCTTCTATCCGCAGAGAGGTGTTCCAGACGGCAAGTGAGCGGT
GGAGGAGAGGGCCTGGCGCCCTCCACTCCCTGACTGCGTGTGTTGCTCCTGCTGTGCT
GAGGGCTGATGCTGAGGACCTCCTTGACTCCTTCTTAGCAACATTCTACAGGACTGCAG
GCACCACCTGTGTGAACCGGACATGAAACTGGTGTGGCCTAGTGCCAAGCTGTTGCAGGC
AGCTGCAGGTGCATCTGCCCGGCTGTGACTCTGTACCCAGCAATGACTGCCTTTACT
GCTGGAACAGTTCCACAAGCACAGTCAGAGCAGCCAGCGGACAATCCTTGAATGCT
CCTGGGTTTCTTGAAGCTGCAGCAGAAATGGAGCTATGAAGACAAAGATCAAAGGCCTCT
GAATGGCTTCAAGGACCAGCTGTGCTCACTGGTATTATGGCTTAACAGACCCAGCAC
CCAGCTTCAAGCTTGTGGCATCCGTACACTGACAGTCTTGGGTGCCAGCCAGATCTCCT
ATCTTATGAGGACTTGGAGCTGGCAGTGGGTACCTGTACAGACTGAGCTTCTGAAGGA



[View online »](#)

```

GGATCCCAGAGTTGCAGGGTGGCAGCACTGGAAGCATCAGGAACCCTGGCTGCTCTCTA
CCCTGTGGCCTTCAGCAGCCACCTCGTACCCAAGCTCGCTGAGGAGCTGCGTGTAGGGGA
GTCAAATTTGACTAACGGAGATGAGCCCACCAATGCTCCCGCATCTGTGCTGTCTGCA
AGCCTTGTGAGCTGTATCAACACATCCCAGCATCGTCAAGGAGACACTGCCTCTGCTGCT
GCAGCATCTCTGGCAAGTGAACAGAGGGAATATGGTTGCACAATCCAGTGACGTTATTGC
TGTCTGTGAGAGCCTCAGACAGATGGCAGAAAAATGTGAGCAGGACCCTGAGAGTTGCTG
GTATTTCCACCAGACAGCTATACCTTGCCTGCTTGCCTTGGCTGTGAGGCTCTATGCC
AGAGAAGGAGCCCTCAGTTCTGAGAAAAGTACTATTGGAGGATGAGGTGTTGGCTGCCAT
GGTGTCTGTCTATTGGCACTGCTACAACCACCTGAGCCCTGAGTTAGCTGCCAGAGTGT
GACACACATTGTGCCCTCTTCTTGGATGGCAACGTGCTCTTTCTGCCTGAAAACAGCTT
CCCAGCAGATTCAGCCATTCCAGGATGGCTCCTCAGGGCAGAGGCGGTGATTGCACT
GCTTATGGCCTTTGTCTGCTCCCTGCCTCGAAATGTGAAATCCCTCAGCTGAACCACT
CATGCGGGAGCTTTTGGAACTGAGCTGCTGCCACAGCTGCCCTTTTCTCCACCCTGC
TGCCAAGTGCTTTCAGGACTCCTCAACAAGCACCTGCAGGGCAGCAGCTGGATGAATT
CCTACAGCTAGCTGTGGACAAAGTGGAGGCTGGCCTGGACTCTGGGCCCTGTCGTAGTCA
GGCCTTCACTCTTCTTCTGGTAACAAAGGCCCTAGTCTCAGATACCATCCTCTCAG
CTCCTGCCTTACAGCCCGGCTCATGGGCCTCCTGAGTGACCCAGAATTAGGTCCAGCAGC
AGCTGATGGCTTCTCTGCTCATGTCTGACTGCACTGATGTGCTGACTCGTGTGGCCA
TGCCGAAGTGCGGATCATGTTCCGCCAGCGGTTCTTACAGATAATGTGCTGCTTTGGT
CCAGGGCTTCCATGCTGCTCCCCAAGATGTGAAGCCAACTACTGAAGGGTCTTTCTCA
TGTACTTAACAGGCTGCCAAGCCTGTACTCTTGCCAGAGCTGCCACGCTTCTTCTCCT
GCTGCTGGAGGCCCTGTCTGCCCTGACTGTGTGGTGCAGCTCTCCACCCTCAGCTGCCT
TCAGCTCTTCTACTGGAAGCACCCCAAGTCATGAGTCTTACGTGGACACCCCTCGTCAC
CAAGTTTCTGAACCTCAGCTCTAGCCCTTCCATGGCTGTCCGGATCGCCGCACTGCACTG
CATGCATGCTCTCACTCGCTGCCACCCTGTGCTGCTGCCGTACAAACCACAGGTGAT
TCGGGCCTTAGCCAAACCCTGGATGACAAGAAGAGACTGGTGCACAAGGAAGCAGTGTG
AGCCAGAGGGGAGTGGTTTCTGTTGGGGAGCCCTGGCAGCTGAGCCCTCAGTCTGGCCT
AGACTGTTCTGACAATCTAACCTGGGATTACTAACTGTTGAGCCATCTTCCCAAAGCAG
GGAAACCCTGGTCTCTGACTGCCTTTCCACAGACACAGCACAAATGCTAGGCCTCTGT
TGCATGGCTGTACAAAGAACATAAGAGTCCATATTTCTAGTGGATTTGAAAATAAGTGT
GTGTGAGACACTTGCCTTTGAAGAAAGTCTAGGGTCTGGGCTCTTGCATTTATATGT
CAGAAAAGGGGCGATATGCTGCTGAGGGGTGAGTGCATATGAGTGTGGCCCTGAGGACCA
GGGCTGGCAGATGTTGTCTACCTGCTGAAGAATAAAGATTTCTTTTGGTAAAAAAAAAA
AAAAAAAAAA
    
```

- Restriction Sites:** Please inquire
- ACCN:** NM_022362
- Insert Size:** 3550 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** The ORF of this clone has been fully sequenced and found to contain one SNP compared with NM_022362.2. The SNP is a valid as published.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022362.2](#), [NP_071757.2](#)

RefSeq Size: 3504 bp

RefSeq ORF: 3093 bp

Locus ID: 64210

UniProt ID: [Q96T76](#)

Cytogenetics: 10q24.1

Protein Families: Druggable Genome, Transcription Factors

Gene Summary: Key component of the cytosolic iron-sulfur protein assembly (CIA) complex, a multiprotein complex that mediates the incorporation of iron-sulfur cluster into apoproteins specifically involved in DNA metabolism and genomic integrity. In the CIA complex, MMS19 acts as an adapter between early-acting CIA components and a subset of cellular target iron-sulfur proteins such as ERCC2/XPD, FANCI and RTEL1, thereby playing a key role in nucleotide excision repair (NER), homologous recombination-mediated double-strand break DNA repair, DNA replication and RNA polymerase II (POL II) transcription (PubMed:22678362, PubMed:22678361, PubMed:29225034, PubMed:23585563). As part of the mitotic spindle-associated MMXD complex, plays a role in chromosome segregation, probably by facilitating iron-sulfur cluster assembly into ERCC2/XPD (PubMed:20797633). Indirectly acts as a transcriptional coactivator of estrogen receptor (ER), via its role in iron-sulfur insertion into some component of the TFIIH-machinery (PubMed:11279242).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) encodes isoform 1. Variants 1 and 4 encode the same protein. Sequence Note: A downstream start codon is selected for this RefSeq based on better conservation with homologous proteins. The use of an alternative upstream start codon, which is present in some species, including human, chimp and macaque (with a subsequent CDS frameshift in the latter species), would increase the protein length from 1030 aa to 1051 aa. The shorter protein is referred to in the literature, including PMIDs 11071939, 11279242 and 11328871.