

## Product datasheet for SC315423

### ZNF391 (NM\_001076781) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF391 (NM_001076781) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF391
Synonyms:	dj153G14.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315423 representing NM_001076781. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAAGCCTCAGAGGGAATACTGCTCAGGGTCTACAAATGAAGAAGACTATAAAACGAAGGCCAA
TTATCAAGGCAAAACAAAATGTCTGCACAGAAGAAATCCTCTTTGAGAACACAGTGGTCAGAAAAGTG
TCAGTGACACTCAAAGAAATTTTCACAGGGGAGGAAGGCCCTGAATCCAGTGAATTTAGTCTAAGCCCA
AACCTTGACGCACAACAGAAAATTCCAAAGGGACATGGATCCCAATATCTAGGAAAACTCCAAGAT
AATTCAGACTTAATTAACACCAAAGACTTTTCTCACAAGAAAACCTTGTAATGCAATGAATGTGAA
AAAGCCTTTAGTTACCAATCAGACCTTCTGTACACAGTAGAATTCATGGTGGAGAAAAGCCTTTTGAA
TGCAACAAATGTGGGAAATCTTTAGCCGAAAGTACACACCTTATTGAACATCAAAGAACTCACACTGGA
GAGAAACCTTATGAATGCAATGAATGTGGAAAAGCTTTTAGCCGGAGCACACATCTTAGTCTACATCAG
AGAATCCATACTGGAGAAAAACCATATGAATGTAGTGAATGTGGAAAAGCCTTTAGCCGAAAGCACTAAC
CTTAGTCAGCATCAGCGAACTATACTCAAGAAAGGCCTTACAAATGTAATGAATGTGGGAAAGCCTTC
GGTGACCGTTCAACCATAATTCAGCATCAACGAATACACACTGGAGAGAATCCCTATGAATGCAGTAAA
TGTGGAAAAGCTTTCAGTTGGATCTCATCGCTTACTGAACATCAGAGAACACACACTGGGGAGAACCCC
TATGAGTGCAGTGAATGTGGGAAAGTGTTCAGTGAAGCTCGTCTTTACAGAACATCAGAGAATCCAC
AGTGGAGAAAAGCCTCACGAGTGTAGAGTGTGTGGAAAGGGCTTCAGTCGAAGCTCATCCCTTATTATT
CATCAGAGAACTCATACCGGGGAGAAGCCGTACAATGTAATGACTGTGGAAAAGCCTTCTGTGAGAGT
TCAACTCTGATCAGACATCAGCACCTTCTACTAAAGAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001076781



[View online >](#)

<b>Insert Size:</b>	1077 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001076781.2</a></u>
<b>RefSeq Size:</b>	4082 bp
<b>RefSeq ORF:</b>	1077 bp
<b>Locus ID:</b>	346157
<b>UniProt ID:</b>	<u><a href="#">Q9UJN7</a></u>
<b>Cytogenetics:</b>	6p22.1
<b>MW:</b>	40.7 kDa
<b>Gene Summary:</b>	<p>May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 all encode the same protein.</p>