

Product datasheet for **SC315433**

RCC1 (NM_001048194) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RCC1 (NM_001048194) Human Untagged Clone
Tag:	Tag Free
Symbol:	RCC1
Synonyms:	CHC1; RCC1-I; SNHG3-RCC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315433 representing NM_001048194. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCACCCAAGCGCATAGCTAAAAGAAGGTCCCCCAGCAGATGCCATCCCCAAAAGCAAGAAGGTG
AAGGACACGAGGGCCGCTGCCTCCCGCCGCTTCTGGCGCCGCTCCTGCCAAGGTGCCTGCGGGCCG
AGCCCTCCTGACCAGAAAACCCGACCAGTCTCACACAGGTCCCACAGCACAGAACCCGGCTTGGTGTG
ACACTAGGCCAGGGCGACGTGGGCCAGCTGGGGCTGGGTGAGAATGTGATGGAGAGGAAGAAGCCGGCC
CTGGTATCCATTCGGAGGATGTTGTGCAGGCTGAGGCTGGGGCATGCACACCGTGTGTCTAAGCAAA
AGTGGCCAGGTCTATTCCTTCGGCTGCAATGATGAGGGTGCCCTGGGAAGGGACACATCAGTGGAGGGC
TCGGAGATGGTCCCTGGGAAAGTGGAGCTGCAAGAGAAGGTGGTACAGGTGTGAGCAGGAGACAGTCAC
ACAGCAGCCCTCACCGATGATGGCCGTGTCTTCTCTGGGGCTCCTCCGGGACAATAACCGTGTGATT
GGACTGTTGGAGCCCATGAAGAAGAGCATGGTGCCTGTGCAGGTGCAGCTGGATGTGCCTGTGGTAAAG
GTGGCCTCAGGAAACGACCCTTGGTGTGCTGACAGCTGATGGTGCCTTACACCTTGGGCTGCGGG
GAACAGGGCCAGCTAGGCCGTGTGCCTGAGTATTTGCCAACCGTGGTGGCCGGCAAGGCCCTCGAACGA
CTCCTGGTCCCCAAGTGTGTGATGCTGAAATCCAGGGGAAGCCGGGGCCACGTGATCCAGGATGCC
TTTTGTGGTGCCTATTTACCTTTGCCATCTCCATGAGGGCCACGTGTACGGCTTCGGCCTCTCCAAC
TACCATCAGCTTGAACCTCGGGCACAGAATCTTGCTTCATACCCAGAACCTAACATCCTTCAAGAA
TCCACCAAGTCTGGGTGGCTTCTCTGGTGGCCAGCACCATACAGTCTGCATGGATTTCGGAAGGAAAA
GCATACAGCCTGGGCCGGCTGAGTATGGCGGGCTGGCCCTGGAGAGGGTGTGAGGAGAAGAGCATA
CCCACCCTCATCTCCAGGCTGCCTGCTGTCTCCTCGGTGGCTTGTGGGCCTCTGTGGGGTATGCTGTG
ACCAAGGATGGTCTGTTTTCCGCTGGGGCATGGGCACCAACTACCAGCTGGGCACAGGGCAGGATGAG
GACGCTGGAGCCCTGTGGAGATGATGGGCAAACAGCTGGAGAACCGTGGTCTTATCTGTGTCCAGC
GGGGGCCAGCATACAGTCTTATTAGTCAAGGACAAAGAACAGAGCTGA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001048194
Insert Size:	1359 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001048194.2
RefSeq Size:	2633 bp
RefSeq ORF:	1359 bp
Locus ID:	1104
UniProt ID:	P18754
Cytogenetics:	1p35.3
MW:	48.1 kDa
Gene Summary:	<p>Guanine-nucleotide releasing factor that promotes the exchange of Ran-bound GDP by GTP, and thereby plays an important role in RAN-mediated functions in nuclear import and mitosis (PubMed:1944575, PubMed:17435751, PubMed:20668449, PubMed:22215983, PubMed:11336674). Contributes to the generation of high levels of chromosome-associated, GTP-bound RAN, which is important for mitotic spindle assembly and normal progress through mitosis (PubMed:12194828, PubMed:17435751, PubMed:22215983). Via its role in maintaining high levels of GTP-bound RAN in the nucleus, contributes to the release of cargo proteins from importins after nuclear import (PubMed:22215983). Involved in the regulation of onset of chromosome condensation in the S phase (PubMed:3678831). Binds both to the nucleosomes and double-stranded DNA (PubMed:17435751, PubMed:18762580). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (a).</p>