

Product datasheet for **SC315443**

CLCC1 (NM_001048210) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLCC1 (NM_001048210) Human Untagged Clone
Tag:	Tag Free
Symbol:	CLCC1
Synonyms:	MCLC; RP32
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315443 representing NM_001048210.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGTGTTCTTTGCTCCTTTGTAATGTCTGTTGCTGGTAGCTGGTTATGCTCATGATGACTGG
ATTGACCCACAGACATGCTTAACTATGATGCTGCTTCAGGAACAATGAGAAAATCTCAGGCAAAATAT
GGTATTTCAAGGGAAAAGGATGTCAGTCCTGACTTGTCATGTGCTGATGAAATATCAGAATGTTATCAC
AAACTTGATTCTTTAACTTATAAGATTGATGAGTGTGAAAAGAAAAGAGGGAAGACTATGAAAGTCAA
AGCAATCTGTTTTAGGAGATACTTAAATAAGATTTTAAATTGAAGCTGGAAAGCTTGGACTTCCTGAT
GAAAACAAAGGCGATATGCATTATGATGCTGAGATTATCCTTAAAAGAGAAAATTTGTTAGAAAATACAG
AAGTTTCTCAATGGAGAAGACTGAAAACCGGTCCTGGATGATGACTAAGTGATATTTTAAATTAAT
TTTAAAGTTTCATGATTTTGAACATGGAAGTGGCGATTGGAAGATTCCTTTGGAGTGGATCCATATAAT
GTGTTAATGGTACTTCTTTGCTGCTCTGCATCGTGGTTTTAGTGGCTACTGAGCTGTGGACATATGTA
CGTTGGTACACTCAGTTGAGACGTGTTTTAATCATCAGCTTTCTGTTGAGTTGGATGGAATGGATG
TATTTATATAAGCTAGCTTTTGCACAGCATCAGGCTGAAGTCGCCAAGATGGAGCCATTAACAATGTG
TGTGCCAAAAAGATGGACTGGACTGGAAGTATCTGGGAATGGTTAGAAGTTCATGGACCTATAAGGAT
GACCCATGCCAAAAATACTATGAGCTCTTACTAGTCAACCCTATTTGGTTGGTCCCAACAAAGGCA
CTTGCAAGTTACATTCACCACATTTGTAACGGAGCCATTGAAGCATATTGGAAAAGGAACTGGGGAATTT
ATTAAGCACTCATGAAGGAAATTCAGCGCTGCTTCATCTTCCAGTCTGATAATTATGGCATTAGCC
ATCCTGAGTTTCTGCTATGGTGTGAAAATCAGTTTATGTGCTGAGACATATAGGCGGCTCTGAGAGC
GAACCTCCCAGGCACTTCGGCCACGGGATAGAAGACGGCAGGAGGAAATGATTATAGACCTGATGGT
GGAGCAGGTGATGCCGATTTCCATTATAGGGGCCAAATGGGCCCACTGAGCAAGGCCCTTATGCCAAA
ACGTATGAGGGTAGAAGAGAGATTTTGAAGAGAGAGATGTTGACTTGAGATTTTCAACTGGCAACAAG
AGCCCTGAAGTGCTCCGGGCATTTGATGTACCAGACGCAGAGGCACGAGAGCATCCACGGTGGTACCC
AGTCATAAATCACCTGTTTTGGATACAAAAGCCCAAGGAGACAGGTGGAATCCTGGGGGAAGGCACACCG
AAAGAAAGCAGTACTGAAAGCAGCCAGTCGGCCAAAGCCTGTCTCTGGCCAAGACACATCAGGGAATACA
GAAGGTTACCCGCAGCGGAAAAGGCCAGCTCAAGTCTGAAGCCGCAGGCAGCCAGACCAAGGCAGC
ACATACAGCCCCGAAGAGGTGTGGCTGGACCAGTGGACAGGATCCGGTCAGCAGCCCCCTGTGGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001048210

Insert Size: 1656 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001048210.2](#)

RefSeq Size: 4835 bp

RefSeq ORF: 1656 bp

Locus ID: 23155

UniProt ID: [Q96S66](#)

Cytogenetics: 1p13.3

Protein Families: Ion Channels: Other, Transmembrane

MW: 62 kDa

Gene Summary: Seems to act as a chloride ion channel.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).