

Product datasheet for SC315545

G protein alpha S (GNAS) (NM_001077488) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	G protein alpha S (GNAS) (NM_001077488) Human Untagged Clone
Tag:	Tag Free
Symbol:	GNAS
Synonyms:	AHO; C20orf45; GNAS1; GPSA; GSA; GSP; NESP; PITA3; POH; SCG6; SgVI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315545 representing NM_001077488. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGGCTGCCTCGGGAACAGTAAGACCGAGGACCAGCGCAACGAGGAGAAGGCGCAGCGTGAGGCCAAC
AAAAAGATCGAGAAGCAGCTGCAGAAGGACAAGCAGGTCTACCGGGCCACGCACCGCCTGCTGCTGCTG
GGTGTGGAGAATCTGGTAAAAGCACCATTGTGAAGCAGATGAGGATCCTGCATGTTAATGGGTTTAAAT
GGAGAGGGCGGCAAGAGGACCCGAGGCTGCAAGGAGCAACAGCGATGGCAGTGAGAAGGCAACCAAAA
GTGCAGGACATCAAAAACAACCTGAAAGAGGCGATTGAAACCATTGTGGCCGCCATGAGCAACCTGGTG
CCCCCGTGGAGCTGGCCAACCCGAGAACCAGTTCAGAGTGGACTACATCCTGAGTGTGATGAACGTG
CCTGACTTTGACTTCCCTCCCGAATTCTATGAGCATGCCAAGGCTCTGTGGGAGGATGAAGGAGTGCCT
GCCTGCTACGAACGCTCCAACGAGTACCAGCTGATTGACTGTGCCAGTACTTCTCTGGACAAGATCGAC
GTGATCAAGCAGGCTGACTATGTGCCGAGCGATCAGGACCTGCTTCGCTGCCGTGTCTGACTTCTGGA
ATCTTTGAGACCAAGTTCAGGTGGACAAAGTCAACTCCACATGTTTGACGTGGGTGGCCAGCGCGAT
GAACGCCGCAAGTGGATCCAGTGCTTCAACGATGTGACTGCCATCATCTTCGTGGTGGCCAGCAGCAGC
TACAACATGGTCATCCGGGAGGACAACCAGACCAACCCGCTGCAGGAGGCTCTGAACCTCTCAAGAGC
ATCTGGAACAACAGATGGCTGCGCACCATCTCTGTGATCCTGTTCCAGAATTTGCTCGTACACTACTCCT
GAGGATGCTACTCCCGAGCCGGAGAGGCCACGCGTGACCCGGCCAAGTACTTCATTGAGATGAG
TTTCTGAGGATCAGCACTGCCAGTGGAGATGGGCGTCACTACTGTACCCTCATTTCACCTGCGCTGTG
GACACTGAGAACATCCGCGTGTGTTCAACGACTGCCGTGACATCATTCAGCGCATGCACCTTCGTCAG
TACGAGCTGCT TAA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
```

Restriction Sites: SgfI-RsrII



[View online >](#)

Plasmid Map:



ACCN: NM_001077488

Insert Size: 1188 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_001077488.3
RefSeq Size:	1997 bp
RefSeq ORF:	1188 bp
Locus ID:	2778
UniProt ID:	O95467
Cytogenetics:	20q13.32
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Calcium signaling pathway, Dilated cardiomyopathy, Gap junction, GnRH signaling pathway, Long-term depression, Melanogenesis, Taste transduction, Vascular smooth muscle contraction, Vibrio cholerae infection
MW:	45.8 kDa
Gene Summary:	<p>This locus has a highly complex imprinted expression pattern. It gives rise to maternally, paternally, and biallelically expressed transcripts that are derived from four alternative promoters and 5' exons. Some transcripts contain a differentially methylated region (DMR) at their 5' exons, and this DMR is commonly found in imprinted genes and correlates with transcript expression. An antisense transcript is produced from an overlapping locus on the opposite strand. One of the transcripts produced from this locus, and the antisense transcript, are paternally expressed noncoding RNAs, and may regulate imprinting in this region. In addition, one of the transcripts contains a second overlapping ORF, which encodes a structurally unrelated protein - Alex. Alternative splicing of downstream exons is also observed, which results in different forms of the stimulatory G-protein alpha subunit, a key element of the classical signal transduction pathway linking receptor-ligand interactions with the activation of adenylyl cyclase and a variety of cellular responses. Multiple transcript variants encoding different isoforms have been found for this gene. Mutations in this gene result in pseudohypoparathyroidism type 1a, pseudohypoparathyroidism type 1b, Albright hereditary osteodystrophy, pseudopseudohypoparathyroidism, McCune-Albright syndrome, progressive osseous heteroplasia, polyostotic fibrous dysplasia of bone, and some pituitary tumors. [provided by RefSeq, Aug 2012]</p> <p>Transcript Variant: This variant (6) is biallelically expressed and uses an alternate in frame splice site, compared to variant 1. It encodes isoform f which has an additional amino acid in the internal region compared to isoform GNASL.</p>