

Product datasheet for SC315555

GRIA4 (NM_001077244) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GRIA4 (NM_001077244) Human Untagged Clone
Tag:	Tag Free
Symbol:	GRIA4
Synonyms:	GluA4; GluA4-ATD; GLUR4; GLUR4C; GLURD; NEDSGA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315555 representing NM_001077244. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGATTATTTCCAGACAGATTGTCTTGTTATTTTCTGGATTTGGGGACTCGCCATGGGAGCCTTT
CCGAGCAGCGTGCAAATAGGTGGTCTTTCATCCGAAACACAGATCAGGAATACACTGCTTTTCGATTA
GCAATTTTCTTCATAACACCAGCCCAATGCGTCGGAAGCTCCTTTAATTTGGTACCTCATGTGGAC
AACATTGAGACAGCCAACAGTTTGTCTGTAACAAACGCCTTCTGTTCCAGTATTCTAGAGGAGTATTT
GCCATTTTGGACTCTATGATAAGAGGTCGGTACATACCTTGACCTCATTCTGCAGCGCCTTACATATC
TCCCTCATCACACCAAGTTTCCCTACTGAGGGGAGAGCCAGTTTGTGCTGCAACTAAGACCTTCGTTA
CGAGGAGCACTCTTGAGTTTGTGGATCACTACGAATGGAAGTGTTCCTGTATGACACAGAC
AGGGGATACTCGATACTCCAAGCTATTATGAAAAAGCAGGACAAAATGTTGGCATGTCAGCGCTATA
TGTGTGGAAAATTTAATGATGTGAGCTATAGGCAACTTCTAGAAGAACTTGACAGAAGACAAGAGAAG
AAGTTTGTAAAGACTGTGAGATAGAGAGACTTCAAACATATTAGAACAGATTGTAAGTGTGGAAAG
CATGTTAAAGGCTACCATTATATCATTGCAAACCTGGGATTCAAGGATATTTCTCTTGAGAGGTTTATA
CATGGTGGAGCCAATGTTACTGGATTCCAGTTGGTGGATTTTAAACACCTATGGTAATCAAACCTAATG
GATCGTGAAGAACTAGATCAGAGAGAGTATCCAGGATCTGAGACTCCTCAAAGTACACCTCTGCT
CTGACTTATGATGGAGTCTTGTGATGGCTGAAACTTTCCGAAGTCTTAGGAGGCAGAAAATGATATC
TCAAGGAGAGGAAATGCTGGGGATTGTCTGGCAAATCCTGCTGCCATGGGGCCAGGGAATTGACATG
GAGAGGACACTCAAACAGGTTCCAATTCAGGGCTGACAGGGAATGTTCAAGTTTGACCACTATGGACGT
AGAGTCAATTACACAATGGATGTGTTTGTGCTGAAAAGCACAGGACCTAGAAAAGTTGGTACTGGAAT
GATATGGATAAGTTAGTCTTGATTCAAGATGTACCAACTCTGGCAATGACACAGCTGCTATTGAGAAC
AGAACAGTGGTTGAACCAATATGCCTCTGATGAAGAATCCTATTTAAGAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001077244

Insert Size: 1302 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001077244.1](#)

RefSeq Size: 3345 bp

RefSeq ORF: 1302 bp

Locus ID: 2893

UniProt ID:	P48058
Cytogenetics:	11q22.3
Protein Families:	Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction
MW:	49.1 kDa
Gene Summary:	<p>Glutamate receptors are the predominant excitatory neurotransmitter receptors in the mammalian brain and are activated in a variety of normal neurophysiologic processes. These receptors are heteromeric protein complexes composed of multiple subunits, arranged to form ligand-gated ion channels. The classification of glutamate receptors is based on their activation by different pharmacologic agonists. The subunit encoded by this gene belongs to a family of AMPA (alpha-amino-3-hydroxy-5-methyl-4-isoxazole propionate)-sensitive glutamate receptors, and is subject to RNA editing (AGA->GGA; R->G). Alternative splicing of this gene results in transcript variants encoding different isoforms, which may vary in their signal transduction properties. Some haplotypes of this gene show a positive association with schizophrenia. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) is missing several coding exons at the 3' end, and contains a novel 3' terminal exon compared to transcript variant 1. This results in a shorter isoform (3) with a different C-terminus compared to isoform 1. Variants 3 and 4 encode the same isoform.</p>