

Product datasheet for **SC315560**

ENT1 (SLC29A1) (NM_001078174) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ENT1 (SLC29A1) (NM_001078174) Human Untagged Clone
Tag:	Tag Free
Symbol:	ENT1
Synonyms:	ENT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315560 representing NM_001078174. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACAACCAAGTCAACAGCCTCAGGACAGATACAAAGCTGTCTGGCTTATCTTCTCATGCTGGGTCTG
GGAACCGTGTCCCGTGGAAATTTTTCATGACGGCCACTCAGTATTTACAAAACCGCCTGGACATGTCC
CAGAATGTGTCTTGGTCACTGCTGAACAGCAAGGACGCCAGGCGTCAGCCGCCCTGCAGCACCC
TTGCCTGAGCGGAACCTCTCAGTGCCATCTTCAACAATGTCATGACCCTATGTGCCATGCTGCCCTG
CTGTTATTCACCTACCTCAACTCCTTCTGCATCAGAGGATCCCCAGTCCGTACGGATCCTGGGCAGC
CTGGTGGCCATCCTGCTGGTGTCTGATCACTGCCATCCTGGTGAAGGTCAGCTGGATGCTCTGCCC
TTCTTTGTCATCACCATGATCAAGATCGTGCTCATTAAATTCATTTGGTGGCATCCTGCAGGGCAGCCTG
TTTGGTCTGGCTGGCCTTCTGCCTGCCAGCTACACGGCCCCATCATGAGTGGCCAGGGCCTAGCAGGC
TTCTTTGCCTCCGTGGCCATGATCTGCGCTATTGCCAGTGGCTCGGAGCTATCAGAAAAGTGCCTTCGGC
TACTTTATCACAGCCTGTGCTGTATCATTTTGACCATCATCTGTTACCTGGGCCTGCCCGCCTGGAA
TTCTACCGCTACTACCAGCAGCTCAAGCTTGAAGGACCCGGGGAGCAGGAGACCAAGTTGGACCTCATT
AGCAAAGGAGAGGAGCCAAGAGCAGGCAAGAGGAATCTGGAGTTTCAGTCTCCAACCTCAGCCCCACC
AATGAAAGCCACTCTATCAAAGCCATCCTGAAAAATATCTCAGTCTGGCTTTCTCTGTCTGCTTCATC
TTCATATCACCATTTGGGATGTTTCCAGCCGTGACTGTTGAGGTCAGTCCAGCATCGCAGGCAGCAGC
ACCTGGGAACGTTACTTCATTCCTGTGTCTGTTTCTTGACTTTCAATATCTTTGACTGGTTGGGCCGG
AGCCTCACAGCTGATTATCATGTGGCCTGGAAAGGACAGCCGCTGGCTGCCAAGCCTGGTGTGGCCCGG
CTGGTGTTTGTGCCACTGCTGCTGTGCAACATTAAGCCCCGCGCTACCTGACTGTGGTCTTCGAG
CACGATGCCTGGTTCATCTTCTCATGGCTGCCTTTGCCCTTCTCAACGGCTACCTGCCAGCCTCTGC
ATGTGCTTCGGGCCAAGAAAGTGAAGCCAGCTGAGGCAGAGACCGCAGGAGCCATCATGGCCTTCTC
CTGTGTCTGGGTCTGGCACTGGGGCTGTTTTCTCCTTCTGTTCCGGGCAATTGTGTA
ACGGTACGGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



View online »

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001078174
Insert Size:	1371 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001078174.1</u>
RefSeq Size:	2344 bp
RefSeq ORF:	1371 bp
Locus ID:	2030
Cytogenetics:	6p21.1
Protein Families:	Transmembrane
MW:	50.2 kDa

Gene Summary:

This gene is a member of the equilibrative nucleoside transporter family. The gene encodes a transmembrane glycoprotein that localizes to the plasma and mitochondrial membranes and mediates the cellular uptake of nucleosides from the surrounding medium. The protein is categorized as an equilibrative (as opposed to concentrative) transporter that is sensitive to inhibition by nitrobenzylthioinosine (NBMPR). Nucleoside transporters are required for nucleotide synthesis in cells that lack de novo nucleoside synthesis pathways, and are also necessary for the uptake of cytotoxic nucleosides used for cancer and viral chemotherapies. Multiple alternatively spliced variants, encoding the same protein, have been found for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. There is an upstream in-frame AUG at nt 255-257. This AUG is not annotated because it has a weak Kozak signal and N-terminal sequencing (PMID 8986748) identified the downstream annotated AUG as the translation start. Variants 1, 2, 3, 4, and 5 encode the same protein.