

Product datasheet for **SC315579**

DEPDC7 (NM_001077242) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DEPDC7 (NM_001077242) Human Untagged Clone
Tag:	Tag Free
Symbol:	DEPDC7
Synonyms:	dj85M6.4; TR2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315579 representing NM_001077242.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCACCGTGCAGGAGAAGGCTGCTGCGCTGAACCTCTCGGCTCTCCACAGCCCCGCGCACAGGCCT
CCGGGTTTCAGTGTAGCTCAGAAGCCATTTGGAGCCACGTATGTATGGAGCAGCATATAACACTCTT
CAAACACAAGTGGAAGTGAAAAACGAAGGCACCGTTTAAACGACATAATGACTGCTTTGTTGGTTCA
GAAGCTGTGGATGTCATTTTTCTCACCTAATTCAGAATAAGTATTTGGTGATGTAGATATTCCTCGA
GCCAAAGTGGTGAGAGTGTGTCAAGCGCTTATGGACTACAAAGTATTTGAAGCAGTTCCAACCAAAGTC
TTTGGAAAAGACAAAAACCTACATTTGAAGATAGTAGTTGCAGCCTTTATAGATTCACCACAATACCT
AACCAAGACAGTCAGTTAGGCAAAGAGAACAACTATATTCACCTGCCAGGTATGCAGATGCATTATTT
AAGTCATCCGATATCAGATCAGCCAGTTTAGAGGACCTGTGGGAAATCTGAGTTAAAGCCTGCCAAC
TCCCTCATGTAAATATCTCTGCAACCTTGTCTCCACAAGTTATTAATGAAGTGTGGCAAGAAGAAACA
ATTGGGCGTCTACTACAACCTGTAGACCTTCCACTTCTTGACTCCTTACTGAAACAGCAAGAGGCTGT
CCTAAAATTCCTCAACCTAAGAGGCAGTCCACCATGGTCAACAGCAGTAACTATCTGGATCGAGGGATT
CTCAAGGCTTATAGTGACTCTCAGGAAGATGAGTGGCTCTCGGCAGCAATTGACTGTTTAGAATACCTT
CCAGACCAAAATGGTGGTGGAAATAAGCAGAAGCTTTCTGAGCAACCAGACCGAACAGACTTAGTGAAA
GAACCTTCTGTTTGATGCCATTGGCAGATATTACAGTAGTAGGGAACCTCTGTTAAATCACTTATCTGAC
GTTTATAATGGAATTGCAGAACCTTAGTGAATGGGAAGACGGAATAGCTTTAGAAGCTACCCAGCTC
CTTCTAAAGCTTTTAGATTTCCAAAATAGAGAAGATTTAGAAGACTACTGTATTTATGCTGTTGCA
GCAATCCTTCTGAGTTTAAATTACAGAAAGAAAGTGACAACCGAATGGTGTGAAAAGGATATTCTCA
AAAGCTATTGTTGACAATAAAAAATTTATCCAAAGGCAAAACAGATCTTCTGGTACTCTTTTAAATGGAT
CATCAGAAAGATGTTTTTAAGATTCTGGAACCTACATAAAATTGTAAGTGTTAAGCTTATGGCCATA
CAGAACGGAAGAGATCCAAATAGAGATGCAGGATATATTTATTGCCAGAGAATTGATCAACGTGACTAT
TCCAACAATACAGAGAAGACAACCAAGATGAGCTGTTGAATTTACTAAAACTCTTGATGAGGATTCA
AACTTTCTGCCAAAGAGAAGAAAAAATTGCTAGGTCAATTCTATAAGTGTCAACCCAGACATCTTTATT
GAGCATTTTGGAGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001077242

Insert Size: 1536 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001077242.1](#)

RefSeq Size: 1774 bp

RefSeq ORF: 1536 bp

Locus ID: 91614

UniProt ID: [Q96QD5](#)

Cytogenetics: 11p13

MW: 58.3 kDa