

Product datasheet for **SC315583**

LAP3 (NM_015907) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LAP3 (NM_015907) Human Untagged Clone
Tag:	Tag Free
Symbol:	LAP3
Synonyms:	HEL-S-106; LAP; LAPEP; PEPS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC315583 representing NM_015907.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTCTTGCTGCCTCTCCGGCTGCGGGGCGAGTAGTCGTCCGACGTCTGGCCGTGAGACGTTTCGGG
AGCCGGAGTCTCCACCGCAGACATGACGAAGGCCTTGTTTTAGGAATCTATTCCAAAGAAAAAGAA
GATGATGTGCCACAGTTCACAAGTGCAGGAGAGAATTTTGATAAATTGTTAGCTGGAAAGCTGAGAGAG
ACTTTGAACATATCTGGACCACCTCTGAAGGCAGGGAAGACTCGAACCTTTTATGGTCTGCATCAGGAC
TTCCCCAGCGTGGTGCTAGTTGGCCTCGGCAAAAAGGCAGCTGGAATCGACGAACAGGAAAACTGGCAT
GAAGGCAAAAGAAAACATCAGAGCTGCTGTTGCAGCGGGTGCAGGCAGATTCAAGACCTGGAGCTCTCG
TCTGTGGAGTGGATCCCTGTGGAGACGCTCAGGCTGCTGCGGAGGGAGCGGTGCTTGGTCTCTATGAA
TACGATGACCTAAAGCAAAAAAGAAGATGGCTGTGTCGGCAAAGCTCTATGGAAGTGGGGATCAGGAG
GCCTGGCAGAAAGGAGTCTGTTTCTTCTGGCAGAACTTGGCAGCCAATTGATGGAGACGCCAGCC
AATGAGATGACGCCAACCAGATTTGCTGAAATTATTGAGAAGAATCTCAAAGTGTAGTAGTAAAACC
GAGGTCCATATCAGACCCAAGTCTTGGATTGAGGAACAGGCAATGGGATCATTCCCTCAGTGTGCCAAA
GGATCTGACGAGCCCCAGTCTTCTTGGAAATCACTACAAGGCAGCCCCAATGCAAACGAACCCCC
CTGGTGTTTGTTGGGAAAGGAATTACCTTTGACAGTGGTGGTATCTCCATCAAGGCTTCTGCAAATAG
GACCTCATGAGGGCTGACATGGGAGGAGCTGCAACTATATGCTCAGCCATCGTGTCTGCTGCAAAGCTT
AATTTGCCATTAATATTATAGGTCTGGCCCTCTTTGTGAAAATATGCCAGCGGCAAGGCCAACAAAG
CCGGGGATGTTGTTAGAGCCAAAAACGGGAAGACCATCCAGTTGATAACACTGATGCTGAGGGGAGG
CTCATACTGGCTGATGCGCTCTGTTACGCACACACGTTTAAACCGAAGGTCATCCTCAATGCCGCCACC
TTAACAGTGGCATGGATGTAGCTTGGGATCAGTGCCACTGGGGTCTTTACCAATTCATCCTGGCTC
TGGAAACAACTCTTCGAGGCCAGCATTGAAACAGGGGACCGTGTCTGGAGGATGCCTCTCTTCGAACAT
TATAAAGACAGTTGTAGATTGCCAGTTGCTGATGTTAACAACATTGGAAAATACAGATCTGCAGGA
GCATGTACAGCTGCAGCATTCTGAAAGAATTCGTAACCTCATCCTAAGTGGGCACATTTAGACATAGCA
GGCGTGATGACCAACAAAGATGAAGTCCCTATCTACGAAAGGCATGACTGGGAGGCCACAAGGACT
CTCATTGAGTTCTTACTTCGTTTCAGTCAAGACAATGCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_015907

Insert Size:	1560 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_015907.2
RefSeq Size:	2100 bp
RefSeq ORF:	1560 bp
Locus ID:	51056
UniProt ID:	P28838
Cytogenetics:	4p15.32
Domains:	Peptidase_M17
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Arginine and proline metabolism, Glutathione metabolism, Metabolic pathways
MW:	56.2 kDa
Gene Summary:	Presumably involved in the processing and regular turnover of intracellular proteins. Catalyzes the removal of unsubstituted N-terminal amino acids from various peptides. [UniProtKB/Swiss-Prot Function]