

## Product datasheet for **SC315585**

### ACBD3 (NM\_022735) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ACBD3 (NM_022735) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACBD3
Synonyms:	GCP60; GOCAP1; GOLPH1; PAP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC315585 representing NM\_022735.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGCGGTGCTGAACGCAGAGCGACTCGAGGTGTCGTGACGGCCCTCACGCTCAGCCCGACCCG
GAGGAGCGGCCTGGGGCGGAGGGCGCCCGCTGCTGCCGCCACCGCTGCCACCGCCCTCGCCACCTGGA
TCCGGTCGCGGCCCGGGCGCCTCAGGGGAGCAGCCCGAGCCCGGGGAGGGCGGCTGGGGCGCGGCG
GAGGAGCGCGGCGGCTGGAGCAGCGCTGGGGTTTCGGCCTGGAGGAGTTGTACGGCCTGGCACTGCGC
TTCTTCAAAGAAAAAGATGGCAAAGCATTTTCATCCAATTATGAAGAAAAATTGAAGCTTGTGGCACTG
CATAAGCAAGTTCTTATGGGCCATATAATCCAGACACTTGTCTGAGGTTGGATTCTTTGATGTGTTG
GGGAATGACAGGAGGAGAGAATGGGCAGCCCTGGGAAACATGTCTAAAGAGGATGCCATGGTGGAGTTT
GTCAAGCTCTTAAATAGGTGTTGCCATCTCTTTCAACATATGTTGCGTCCCACAAAATAGAGAAGGAA
GAGCAAGAAAAAAAAGGAAGGAGGAAGAGGAGCGAAGGCGGCGTGAAGAGGAAGAAAGAGAACGTCTG
CAAAGGAGGAAGAGAAACGTAGGAGAGAAGAAGAGGAAAGGCTTCGACGGGAGGAAGAGGAAAGGAGA
CGGATAGAAGAAGAAAGGCTTCGGTTGGAGCAGCAAAAGCAGCAGATAATGGCAGCTTTAAACTCCAG
ACTGCCGTGCAGTTCAGCAGTATGCAGCCCAACAGTATCCAGGGAACACGAACAGCAGCAAATTCCTC
ATCCGCCAGTTGCAGGAGCAACACTATCAGCAGTACATGCAGCAGTTGTATCAAGTCCAGCTTGCACAG
CAACAGGCAGCATTACAGAAACAACAGGAAGTAGTAGTGGCTGGGTCTTCTTGCCTACATCAAAAA
GTGAATGCAACTGTACCAAGTAATATGATGTCAGTTAATGGACAGGCCAAAACACACACTGACAGCTCC
GAAAAAGAACTGGAACCAAGAGTGCAGAAGAAGCCCTGGAGAATGGACCAAAAGAATCTTCCAGTA
ATAGCAGCTCCATCCATGTGGACAGCAGCTCAGATCAAAGACTTCAAAGAGAAGATTGAGCAGGATGCA
GATTCGGTGATTACAGTGGGCCGAGGAGAAGTGGTCACTGTTGAGTACCCACCCATGAAGAAGGATCA
TATCTCTTTTGGGAATTTGCCACAGACAATTATGACATTGGGTTTGGGGTGATTTTTGAATGGACAGAC
TCTCAAACACTGCTGTGACGTCATGTGAGTCCAGCGATGACGACGAGGAGGAAGAAGAAAAC
ATCGGTTGTGAAGAGAAAGCCAAAAGAATGCCAAACAGCCTTTGCTGGATGAGATTGTGCCTGTGTAC
CGACGGGACTGTCATGAGGAGGTGTATGCTGGCAGCCATCAATATCCAGGGAGAGGAGTCTATCTCTC
AAGTTTGACAACCTACTCTTTGTGGCGTCAAATCAGTCTACTACAGAGTCTATTATACTAGATAA
```

```
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_022735
<b>Insert Size:</b>	1587 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_022735.3</a></u>
<b>RefSeq Size:</b>	3581 bp
<b>RefSeq ORF:</b>	1587 bp
<b>Locus ID:</b>	64746
<b>UniProt ID:</b>	<u><a href="#">Q9H3P7</a></u>
<b>Cytogenetics:</b>	1q42.12
<b>Domains:</b>	ACBP
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	60.6 kDa
<b>Gene Summary:</b>	The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is involved in the maintenance of Golgi structure and function through its interaction with the integral membrane protein giantin. It may also be involved in the hormonal regulation of steroid formation. [provided by RefSeq, Jul 2008]