

## Product datasheet for **SC315592**

### TRIM23 (NM\_033227) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TRIM23 (NM_033227) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIM23
Synonyms:	ARD1; ARFD1; RNF46
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL4</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >OriGene ORF sequence for NM\_033227 edited  
 ATGGCTACCCTGGTTGTAACAAGCTCGGAGCGGGAGTAGACAGTGGCCGGCAGGGCAGC  
 CCGGGGACAGCTGTAGTGAAGGTGCTAGAGTGTGGAGTTTGTGAAGATGTCTTTTCTTTG  
 CAAGGAGACAAAGTTCCCGTCTTTTGTGTTTGTGGCCATACCGTCTGTCACTGACTGTCTC  
 ACTCGCTACCTCTTTCATGGAAGAGCAATCCGTTGCCATTTGATCGACAAGTAACAGAC  
 CTAGGTGATTACAGGTGTCTGGGGATTGAAAAAATTTTGTCTTATTGGAGCTTTTGAA  
 CGACTGCAGAATGGCCATTGGTTCAGTATGGAGCTGCAGAAGAATCCATTGGGATATCT  
 GGAGAGAGCATCATTCGTTGTGATGAAGATGAAGCTCACCTTGCCCTGTATATTGCACT  
 GTGTGTGCAACTCATTTGTGCTCTGAGTGTCTCAAGTTACTCATTCTACAAAGACATTA  
 GCAAAGCACAGGCGAGTTCCTCTAGCTGATAAACCTCATGAGAAAATATGTGCTCTCAG  
 CACCAGGTGCATGCCATTGAGTTTGTGTTGCTTGGAAAGAAGTTGTCAAAGTACGCCACTC  
 ATGTGCTGTGCTGCAAAGAATATGGAAAACACCAGGGTCACAAGCATTAGTATTGGAA  
 CCAGAAGTAATCAGATCCGAGCATCAATTTAGATATGGCTCACTGCATACGGACCTTC  
 ACAGAGGAAATCTCAGATTATCCAGAAAATTAGTTGGAATTGTGCAGCACATTGAAGGA  
 GGAGAACAATCGTGGAAAGTGAATTGGAATGGCTCACACAGAACATGTACCAGGGACT  
 GCAGAGAATGCCCGGTCAATGATTCGAGCTTATTTTATGATCTACATGAAACTCTGTGT  
 CGTCAAGAAGAAATGGCTCTAAGTGTGTTGTGATGCTCATGTTCTGTAAGAAATGATTGG  
 CTCAGGCAGCAACAAGAAGATATGACTATTTTGTGTCAGAGGTTTCTGCAGCCTGCCTC  
 CACTGTGAAAAGACTTTGCAGCAGGATGATTGTAGAGTTGCTTGGCAAAACAGGAAATT  
 ACAAGGTTACTGGAAACATTGCAGAAACAGCAGCAGCAGTTTACAGAAGTTGCAGATCAC  
 ATTCAGTTGGATGCCAGCATCCCTGTCACTTTTACAAAGGATAATCGAGTTCACATTGGA  
 CCAAAAATGGAATTCGGGTCGTTACGTTAGGATTGGATGGTGGTGGAAAACTACTATC  
 TTGTTTAAAGTAAAACAGGATGAATTCATGCAGCCCATTCACAACAAATGGTTTTAACGTG  
 GAAACTGTAGAATATAAAAATCTAAAATCACTATTTGGGATGTAGGTGGAAAACACAAA  
 TTAAGACCATTGTGAAAACATTATTACCTCAACTCAAGCTGTTGTGTTTGTGTAGAT  
 AGCAGTCATAGAGACAGAATTAGTGAAGCACACAGCGAACTTGCAAAGTTGTTAACGGAA  
 AAAGAACTCCGAGATGCTCTGCTCCTGATTTTGTCTAACAAACAGGATGTTGCTGGAGCA  
 CTGTCACTAGAAAGAAATCACTGAACTACTCAGTCTCCATAAATTATGCTGTGGCCGTAGC  
 TGGTATATTCAGGGCTGTGATGCTCGAAGTGTGTTTTCAGATAATATGTGACCAATATACT  
 GGGAAAGAGGTAGTCACAGAGAAAGGTAA

- Restriction Sites:** Please inquire
- ACCN:** NM\_033227
- Insert Size:** 2500 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_033227.2](#), [NP\\_150230.1](#)

**RefSeq Size:** 2429 bp

**RefSeq ORF:** 1710 bp

**Locus ID:** 373

**UniProt ID:** [P36406](#)

**Cytogenetics:** 5q12.3

**Domains:** zf-B\_box, RING, RAB, BBC, SAR, ARF, arf

**Protein Families:** Druggable Genome

**Gene Summary:** The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein is also a member of the ADP ribosylation factor family of guanine nucleotide-binding family of proteins. Its carboxy terminus contains an ADP-ribosylation factor domain and a guanine nucleotide binding site, while the amino terminus contains a GTPase activating protein domain which acts on the guanine nucleotide binding site. The protein localizes to lysosomes and the Golgi apparatus. It plays a role in the formation of intracellular transport vesicles, their movement from one compartment to another, and phospholipase D activation. Three alternatively spliced transcript variants for this gene have been described. [provided by RefSeq, Jul 2008]  
Transcript Variant: This variant (beta) is alternatively spliced near the 3' end of the coding region. This introduces a frameshift and the encoded isoform beta has a unique C-terminus compared to isoforms alpha and gamma.