

## Product datasheet for **SC315602**

### Cylicin 1 (CYLC1) (NM\_021118) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cylicin 1 (CYLC1) (NM_021118) Human Untagged Clone
Tag:	Tag Free
Symbol:	CYLC1
Synonyms:	CYCL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315602 representing NM\_021118.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTTCCAAGTTGCTAAAAGTAAACATCAGAACATATGATAATTCATTCCAATCAGTGAATCA
AGCAGAAAATCATGGAATCAAAAACACTTTGCTTTGACATTTCCAAACCACTCCAGAGAGGTACAAAT
GATAAATCAAGACCTTTGAAATCACAAATAACAGTTACTAGACATGACAAAAGAAAAGTAGAAGAAGGC
CAGAAAACCAGCTCATAAATGGATAAGGCATTCTTTTCAGAAAAATTTGCAATGGCCACCCATTTACACA
GCTGCCAGGGAACAGACTCCATTAGACATCTTTATACTTCCAAAACCCATCTTAAAAAAGCAGAATAT
AAAAAGTCCAAGATGAAAAAGGAGGAACACCTTTGAAGAAAGATTCCAAGAAAAAGGAGGTTTCATAT
GCAACAAATCCAGAATCCAAGCAATAGTAGAAGAAAACTAAAAGACAAAATGAGGCAGATAAACT
CCCTTAAAATCATCACATGAAAATGAACAATCCAAGAAGTCAAAATCCAGTTCAGAACTAATCCAGAA
TCCAAAATTCTAAGACAGTCTCAAAAATTGTTCAAAAAGATAAGAAAGATTCAAAGAATTCGAAG
AAGACAAACACTGAATTCCTACATACAAAAGAACAATCCAAGAAAGATTGAAAGAGTCAAAGACTAGT
AATGATCCCATATCAGAGATTTGCTCAGAAAATAGTTTAAATGTTGATTTCTCATGTTAGTGGGACAG
TCTGATGATGAATCCATAAATTTGATGCATGGTTAAGGAATTACTCACAGAATAATTCAAAGAATTAT
TCTTTGAAGTATACAAAGTATACAAAAGAGGACACAAAAAGAATGCAAAGAAAAGCTCTGATGCTGAA
TCTGAAGACTCAAAGGATGCTAAGAAAGATTCAAAGAAAGTTAAGAAAAATGCAAGAAAGATGACAAG
AAAAAGGATGTAAGAAGGACACAGAGTCTACTGATGCTGAATCTGGAGACTCAAAGGATGAAAGGAAA
GATACAAAAGAGGATAAGAAAAAATTAAAGAAAGATGACAAGAAAAAGGACACAAAAGAGTACCCAGAG
TCTACTGATACTGAATCAGGAGATGCAAAGGATGCAAGAAATGATTCAAGAAATTTGAAGAAAGCTTCA
AAGAATGATGACAAGAAAAAGGATGCAAAGAAAATTACATTCTCTACTGATTCTGAATCTGAAGTGGAG
TCAAAGGAGAGTCAGAAAGATGAAAAAAGGATAAAAAAGATTCAAAGACAGATAATAAAAAAGTCTGTC
AAGAATGATGAAGAGTCTACTGATGCTGACTCTGAACCGAAGGGAGATTCAAAAAAGGGTAAAAAGGAT
GAAAAGAAGGGGAAGAAAGATTCAAAGAAAGATGACAAAAAGAAGGATGCAAAGAAAAATGCAGAATCT
ACTGAAATGGAATCTGATTTGGAGTTAAAGAAGGACAAGAAACACTCAAAGGAAAAGAAAGGTTCAAAG
AAAGATATCAAGAAGGATGCAAGAAAGGACACAGAGTCTACTGATGCTGAATTTGATGAATCTTCCAAG
ACAGGCTTTAAAACATCTACAAAATCAAAGTTCAGATACTGAATCTGAAGAGTCACTATATAACCT
GGGGCTAAGAAGAAAATTGATGAATCAGATGGCACATCTGCAAAATCAAAGATGGAAGGACTGGAATCA
AAGAGAGGATTCAGAATGTCATCAAAAAGACTACATTCAATGAAAAAGGGGAAAAGCAAGTACAGGT
AGAGTTCCATCAAGAGAAAAACCACACTCCCTGCTTGTGAGCCTTCTCTACCATCACCAAAGGTC
AGACGTCTTTGTTGGTGAAGATGCCTCCTCCACCTCCAAAACCAAGATATGCTCCTTTGCTGAAGCA
CCGTGGATTATAAGCTGCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_021118

**Insert Size:** 1956 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_021118.2](#)

**RefSeq Size:** 2143 bp

**RefSeq ORF:** 1956 bp

**Locus ID:** 1538

**UniProt ID:** [P35663](#)

**Cytogenetics:** Xq21.1

**MW:** 74.2 kDa

**Gene Summary:** This gene encodes a sperm head cytoskeletal protein. The encoded protein is associated with the calyx of spermatozoa and spermatids. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]  
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).