

Product datasheet for **SC315603**

TMEM16K (ANO10) (NM_018075) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM16K (ANO10) (NM_018075) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANO10
Synonyms:	SCAR10; TMEM16K
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315603 representing NM_018075.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAAGTGACCTTATCAGCTTTGGATACTTCTGAGAGTTCTTTCACACCTTTGGTGGTCATAGAAGCTT
GCTCAGGATGTCAAAGAAGAAACCAAGAATGGCTGAAAAACAGAATTATAGCTAAAAAAAAGATGGA
GGTGCCAGTTGTTGTTAGACCATTGTTAAATAAATATGAACAAGAACTAGAAAATCAGAAGCTTA
TATCTTGTGGTCCCAAGATTAGAATGTTACTAGGGCGAGAAGCAGTGGGATTGGTAAAAGAGTGC
AATGATAACACCATGAGAGCCTTACATACAGAACCAGACAGAAGCTTCAAAGGTTTTGATGATAACAAT
GATGATTTCTGACAATGGCAGAATGTCAATTCATTATCAAACATGAACTTAAAAATCTTAGAGCTAAA
GATGAAAAATGATCCCTGGTTACCCTCAGGCAAAGTTGTATCCAGGAAAATCATTGTTGAGAAGATTG
CTCACGCTGGCATCGTGATTACAGGTGTTCCACTGCATGACAGTGAAGCCCTGAAGAAGCTTGAGGAC
ACCTGGTACACTCGGTTTGCTTTGAAGTATCAGCCCATAGACAGTATTCTGGCTACTTTGGGGAAACA
ATTGCTCTGTACTTTGGATTTTTGGAGTATTTCACTTTTGCATTAATCCCATGGCTGTCATTGGGTTA
CCTTACTACTTGTGTTGGGAAAGACTATGACAAGTACGTGATCTTGCCTCGTTCAACCTCATCTGG
TCCACGGTGATTCTGAACTGTGGAAGCGTGGCTGTGCCAACATGACCTACAGGTGGGGGACACTGCTC
ATGAAGAGAAAAGTTTGGAGAGCCCCGGCCAGGATTTATGGTGTCTTGGGTATCAATTCATCACTGGG
AAGGAGGAGCCTCTGTACCCAGCTACAAGAGACAGTTGCCGATTTACCTGGTCTCCCTGCCATTCGTG
TGCTCTGCCTCTATTTCTCACTGTATGTCATGATGATTTACTTCGACATGGAGGTTTGGGCCTTGGGT
CTACATGAGAACAGCGGGTCTGAGTGGACCAGTGTCTGTTGTATGTGCCAGCATCATCTATGCCATT
GTGATTGAGATCATGAATCGTCTCTATCGATATGCTGCCGAGTTTTAACTTCATGGGAGAATCACAGA
TTGGAATCTGCCATCAGAACCATCTAATTCTGAAAGTTTTAGTGTCAACTTCCCTCAATTGCTTTGCC
TCACTCTTCTATATTGCCTTTGTCTTGAAGATATGAAGCTTTTGCGCCAGAGCTTGCCACTCTCCTA
ATTACCTCCCAGATCCTCAACCAATTATGGAATCTTTTCTTCTTATTGGCTCCAAAGGAAGCATGGT
GTGCGGGTGAAGAGGAAGGTGCAGGCTTTAAAGGCAGACATTGATGCTACATTATATGAACAAGTCATC
CTGAAAAAGAAATGGGAACTTATTTGGGCACCTTTGATGATTACTTGGAGTTATTCCTGCAGTTTGGT
TATGTGAGCCTTTTCTCCTGTGTTTACCATTAGCAGCTGCCTTTGCTGTGTTAAATAAATTCAGTAA
GTAATTCAGATGCCTTAAAAATGTGCAGGGTCTTCAAACGTCCATTCTCAGAACCTTCAGCCAATATT
GGTGTGTGGCAGTTGGCTTTTGAACGATGAGTGTATATCTGTGGTCACTAACTGTGCGCTGATTGGA
ATGTCACCACAAGTGAATGCAGTCTTCCAGAATCAAAGCAGACCTCATTGATTGTAGTAGCAGTG
GAGCACGCACTCCTGGCTTTAAAGTTTATACTTGCAATTTGCCATACCTGATAAGCCACGGCATATCCAG
ATGAAACTAGCCAGACTGGAATTTGAGTCTTTGGAGGCACTCAAGCAGCAGCAAATGAAGCTCGTGACC
GAGAACCTGAAGGAGGAACCAATGGAAGCGGGGAGGAGAAGGCAACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_018075

Insert Size: 1983 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018075.3](#)

RefSeq Size: 2743 bp

RefSeq ORF: 1983 bp

Locus ID: 55129

UniProt ID: [Q9NW15](#)

Cytogenetics: 3p22.1-p21.33

Domains: DUF590

Protein Families: Transmembrane

MW: 76.3 kDa

Gene Summary: The transmembrane protein encoded by this gene belongs to the anoctamin family of calcium-activated chloride channels, also known as the transmembrane 16 family. The encoded protein contains eight transmembrane domains with cytosolic N- and C-termini. Defects in this gene may cause autosomal recessive spinocerebellar ataxia-10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2016]
Transcript Variant: This variant (1) encodes isoform (1). Variants 1 and 11 encode the same isoform.