

Product datasheet for **SC315604**

MS4A14 (NM_001079692) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MS4A14 (NM_001079692) Human Untagged Clone
Tag:	Tag Free
Symbol:	MS4A14
Synonyms:	MS4A16; NYD-SP21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315604 representing NM_001079692.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGTCAACATCCCAGGACAGAAGGGCAACTCAGTCATCACTATAAAACCAAACGAACTGTATTG
ACTGCATTTCCCTACAGACCTCATAGCTCTCTGCTGGATTTTCTGAAGGGAGAGCCAAGAGTCTTGGGG
GCTACCCAGATCCTGCTTGTCTTAATCATTGTGGGCTTTGGAACATATTTGCACTTAATTACATCGGT
TTCTCCCAAAGACTTCCCTTGTGTCTCACAGGATATCCATTCTGGGAGCACTTATTGGTCAAGGT
GTCACGGGCATGAATGTTATCAGCTCCTTGGTTGCGATAACTGGGATTACTTTCACCATTCTCAGCTAC
AGACATCAAGACAAGTACTGCCAGATGCCATCCTTTGAAGAAATATGTGTTTTCAGTAGAACTCTTTTC
ATTGTTCTGTTTTCTTGCCTTCGGATGTTACTCAAAATAGTGAACAACCTGCCCCAGAAGAAAATGAT
CAATTACAATTTGTGCTTCAAGAAGAGTTTTCCAGTGATGATTCAACAACAAATGCACAATCTGTTATC
TTTGGAGGCTATGCTTTCTCAAGTTAACTCTCTAGGAGTCTTTAGTCTCCCAACCAGGTAATAAA
GGTAGAGAATTTGTGCCAGATGAACAAAAGCAAAGTATCCTTCCATCTCCAAATTTTCAGAGGAAGAA
ATTGAACCTTTGCCTCCACACTAGAGAAAAAGCCCTCAGAAAAATATGTCCATTAGCTAGACTCTACA
TTTAAACAAATGAAAGATGAAGATCTACAATCTGCTATTGTACAACCTTCTCAAATGCAAACCAAGCTT
CTGCAGGACCAAGCTGCGTCACTCCAAGTTTTTCCATCCATTCTGCACTAAAACCTCGAAGATATATCA
CCTGAAGACTTGCCATCCCAAGCTCTACAGTAGAAGGCTGTGAGAACAACCATGCCATCTAAGTCT
ACATCATCCCATGTCAACAGCTTCTAATCTGACAGCTAATGACCTGCCCCCTCAAGGCATACCTATCC
CAAGACACATCATCTCAAGATATGCTGTTTCATGACATGACATCCCAAGATATGCAATCCCTAGATATG
CTATCTCAAGACACACCATCCCACGCCATGCCACCTCAAGACATACCTTCCCAAGATATGCTATCCCAA
GCTCTATCAGCGCATGCCATATTACCTGAAGCCTCAACATCCCATTGTCAGTTCCCTGAAATACAA
CACCTACTTCAGCAGCCCCAGATCTTCAACCAGAAAACACTGAACCTCAAAACAGCAAAATTTTACAA
ATGTCATATCAAGATATTAGATCAGAAGTTATGGAAGAGACCAAGAATGGAATCTGAGGAGGAATC
CATAGAAGAAAATCCTCAAGACGGCATTCTTAAACCAGCAAACCAAGCCTTGCAATACTTAAGGAGA
CATTCTTTAGACGTGCAAGCCAAAGGCCAGAAATCCTCAAGAGGCATTCTTAGATCAGCAAAGCAAA
GGCTGGCAATCTCAAAGCAGAAATCCTTAGACCAGCAAATCAAAGACTGGCTATCCCCAAAGAGGCAC
TCCGTAGATAAGCAAGCTCAACTTAATCAAACCTAAAGAGCAACTCCAGATCAGCAAGCTGAAGATCAG
CAAGCCAAAGGGGAACAATACCCAGAAGGACAATCTAAAGATGGACAAGTTAAAGACCAGCAGACTGAT
AAGGAGCAAAACTCAAAGAAGCAAACCCAGGATCAGCAAACCTGAAGACCAGCCGCCCAAGAGAAGAAA
TCCCCGAAAGGACAATTCCAAAATGTTCAAGCCGAAGGACAGCAAGCTCAGGTGGAGAAAGTGCCAAAA
CTGTTATGCCAAGATTGAGAATCCCAAATACAGCAATACCAATTCTGGCAATTCCACAAAGGCAATCTC
CAGGCTGGACAACCCAGGACTGTCAATCTTTTGCCCAAGAATCCCTGACTGGATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001079692

Insert Size: 1989 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001079692.2</u>
RefSeq Size:	3302 bp
RefSeq ORF:	1989 bp
Locus ID:	84689
UniProt ID:	<u>Q96JA4</u>
Cytogenetics:	11q12.2
Protein Families:	Transmembrane
MW:	74.7 kDa
Gene Summary:	May be involved in signal transduction as a component of a multimeric receptor complex. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) has multiple coding region differences, compared to variant 4, which results in a shorter protein (isoform 2), compared to isoform 4. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.