

Product datasheet for **SC315612**

C2CD2 (NM_015500) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C2CD2 (NM_015500) Human Untagged Clone
Tag:	Tag Free
Symbol:	C2CD2
Synonyms:	C21orf25; C21orf258; TMEM24L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315612 representing NM_015500.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCATGGCCCGGCTGGGCTCGTGGCTCGGGGAGGCGCAGTGGCTCGCGCTGGTGTGCTCTTCGTC
CGGGCCCTGGCCACGGTAGCCCTGTACCTGGCGCAGTGGGCGCTGGCCAGGGCGGACCCAGCCCAG
CGGCGGGCGGTGGAGCCTGGAGAGGGGCCGCGCCCGGGTCCGACGCGCTGCTCTCCTGGATCCTGACG
CTGGGCGACTGGAGGAGCAGTGGCAGGCGGCCTGGGTGACCGCCCTGAACGAGGAGGCCGAGAGGAAA
GGGGGCCACCTTTCCTGTCTTTGAGGAGGACCCGCGCAGCAGGCACTGGAGCTGGTGGTGCAGGAG
GTCTCCAGCGTGTCTCAGGTGGCGGAGGAGAAGGTGGTGGTCTGTACGTGGTGGCCAGGCTATTGAG
TTCTTGGTCAGCGAGACGCTGCCTTGGGTGCTGGATGCCGGCTGTACGACATGCGGCTCTCCCTTTC
CATTTACAGCTGGAGTTCACATGAAGGAGAAGAGAGGACCTCCAGATTAGCTGGTCTTTCATCAGT
GTGCCGAAATGGCCGTTAATATCCAGCCAAAGCACTGGGGGAGGACCAGGTGGCTGAGACAAGTGGC
ATGTCTGACGTTCTCAAGGACATCTTGAAGCATTTGGCTGGTTCGCCTCTCCATCAGTGGTCTGATT
ACCAAGCCACGACTGTAAGGAAGCTCAGAACTTACAGTGTGCTGCATCTACTGCTCAGGAATCCTGT
CCTCTAAACCTCCAAGGCTCACGAGCTGAAGCTACTGGTGAGGAACATCCACGCTTGTGCTCAGC
GAGCCCGTGCCTCAGGCCACATTAATGCAGTGTGCGTGTGAGCTGAACGATCCTGTTGAGAGTTTC
TCCAGCACCTGACGAAAAACACTCCGGACCTCATGTGGGAAGAGGAATTCACCTTCGAGTTGAATGCC
AAGTCGAAGGAGTTACACCTGCAGATTTAGAGGCTGGGCGATCCTCAGAAGGTCTGCTGGCGACGGCG
ACAGTTCTCTGGACTTATTTAAGAAGCAGCCTTCTGGCCACAGAGCTTACGCTGACCGAGCGGTCT
GCCTGGCGCAGCTCGGTGCTGGGCTCGGTACGGCAGAGTCTCTTACATGGAACCTGGTGAATGAAA
TCCTGGCCATCCCTCCCTGTTCTGCTGCAAAAATAGAAAAGGACCGCACGGTGTGATGCCCTGTGGG
ACTGTGGTCACTACTGCTGCTGTGAAGACCAAGCCTCGCGTGCACGTGGGGAGGGCGTCCCGGCTG
AGCTCTGATTCTCCGGTGAAGACTCCCATCAAGGTGAAGGTGATCGAGAAGGACATCTGTCCAGGCC
ATCGCCTGCCGACGCGCCCGTCCAGCAAAACTCTCTTCTTCCAGACACAGAATTGTTGGTGTGAAT
GGTTCGGATCCAGTGGCTGAAGTGGCCATTCGACAGCTCAGTGAATCTTCAAAGCTGAACTCAAGTCG
CCACGGAAGAAAAGCACTATTATCATATCAGGGATCTCCAAGACCTCACTATCTCAGGACCACGACGCT
GCCCTGATGCAGGGCTACACGGCTCTGTGGACAGCACCCACCAGGAGGACGCCCATCCCATCCGGAG
AGGGCGGCAGCCTCTGCCCGCCAGAGGAAGCCGAGTCAAGCCAGGCATCCCTTGCCCCAAGCCCCAG
GAGGACGAGCTAGACTCCTGGACTTGGAGAAGGAGCCACAGGCCGCGCATGGAGCAGCCAGGTCTGTG
CTGACCCCGACGGTGTGAGCTGTGAGAGGCTCCATGAGTGTCTTGGAGCCGGGCACTGCCAAAAG
CATAAAGGAGGAATTCTAAGGAAAGGTGCAAAGCTGTTCTTCCGCGGGCGGCATCAACAGAAAGACCCA
GGCATGAGTCAGTCACACAATGACCTTGTGTTCTTGGAGCAGCCAGAGGGTTCCCGGAGGAAAGGCATC
ACCCTCACCAGGATCCTGAACAAGAAGTGTCTCTCAGGCACAGAAAACAAGAACCATGAACGGTGGC
CCCGTGGAGCCCTGCACGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_015500

Insert Size: 2091 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015500.1](#)

RefSeq Size: 6349 bp

RefSeq ORF: 2091 bp

Locus ID: 25966

UniProt ID: [Q9Y426](#)

Cytogenetics: 21q22.3

Protein Families: Transmembrane

MW: 75.5 kDa