

Product datasheet for **SC315613**

RNF214 (NM_001077239) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RNF214 (NM_001077239) Human Untagged Clone
Tag:	Tag Free
Symbol:	RNF214
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001077239, the custom clone sequence may differ by one or more nucleotides

```

ATGGCAGCGTCTGAGGTTGCTGGTGTGGCCAAATGCCCCAGTCTCCGGAATCTTCT
AGTTTATGTGCTTCCAAATCAGACGAAGGTCTCCAGATGGTCTAAGCACAAAGACTCT
GCACAGAAGCAGAAGAACTCGCCTCTGTTGAGTGAAGTAGCCAAACAATAACCAAGGAG
AATAACAGAAATGTCCATTTGGAGCACTCAGAGCAGAATCCTGGTTCATCAGCAGGTGAC
ACCTCAGCAGCGCACACAGGTGGTTTTAGGAGAAAACCTTGATAGCCACAGCCCTTTGTCTT
TCTGGCAGTGGGTCTCAGTCTGATTTGAAGGATGTGGCCAGCACAGCAGGAGAGGGGG
GACACAAGCCTTCGGGAGAGCCTCCATCCAGTCACTCGGTCTCTTAAGGCAGGGTGCCAT
ACTAAGCAGCTTGCCCTCCAGGAATTGCTCTGAAGAGAAAATCCCCACAAACCTCCATCCTA
AAGGAAGGTAACAGGGACACAAGCTTGGATTTCCGACCTGTAGTGTCTCCAGCAAATGGG
GTTGAAGGAGTCCGAGTGGATCAGGATGATGATCAAGATAGCTCTTCCCTGAAGCTTTCT
CAGAACATTGCTGTACAGACTGACTTTAAGACAGCTGATTCAGAGGTAACACAGATCAA
GATATTGAAAAGAATTTGGATAAAAATGATGACAGAGAGAACCCCTGTTGAAAGAGCGTTAC
CAGGAGGTCCTGGACAAACAGAGGCAAGTGGAGAATCAGCTCCAAGTGAATTAAGCAG
CTTCAGCAAAGGAGAGAAGAGGAAATGAAGAATACCAGGAGATATTAAGGCTATTCAG
GATGTGACAATAAAGCGGGAAGAAAACAAAGAAGAAGATAGAGAAAAGAGAAGGAGTTT
TTGAGAAGGAGCAGGATCTGAAAGCTGAAATTGAGAAGCTTTGTGAGAAGGGCAGAAGA
GAGGTGTGGGAAATGGAAGTGGATAGACTCAAGAAATCAGGATGGCGAAATAAATAGGAAC
ATTATGGAAGAGACTGAACGGGCTGGAAGGCAGAGATCTTATCACTAGAGAGCCGGAAA
GAGTTACTGGTACTGAACTAGAAGAAGCAGAAAAAGAGGCAGAATTGCACCTTACTTAC
CTCAAGTCAACTCCCCAACACTGGAGACAGTTCGTTCCAACAGGAGTGGGAGACGAGA
CTGAATGGAGTTCGGATAATGAAAAGAATGTTCTGTACCAATTTAATAGTCATATCCAG
TTAGTGAGGAACGAGCCAAGCTGAGCAGCCTTCCCAAATCCCTACTCCCACTTTACCT
CCACCCCATCAGAGACAGACTTCATGCTTCAGGTGTTTCAACCCAGTCCCTCTCTGGCT
CCTCGGATGCCCTTCTCCATTGGGCAGGTCACAATGCCCATGGTTATGCCAGTGAGAT
CCCCGCTCCTTGTCTTTCCCAATCCTGAACCCTGCCCTTTCCAGCCAGCCAGCCTTCC
TCACCCCTTCTGGCTCCCATGGCAGAAATAGCCCTGGCTTGGGTTCCCTTGTGACCCCC
CACGGTCCACACATGCCCTGCCGCTCCATCCCACCTCCCCAGGCTTGGGCGGTGTT
AAGGCTTCTGCTGAAACTCCCCGGCCCAACCAGTAGACAAACTGGAGAAGATCCTGGAG
AAGCTGCTGACCCGGTTCACACAGTGAATAAGGCCAGATGACCAACATTCTTCAGCAG
ATCAAGACAGCAGTACCACCATGGCAGGCTGACCATGGAGGAACTTATCCAGTTGGTT
GCTGCACGACTGGCAGAACATGAGCGGGTGGCAGCAAGTACTCAGCCACTTGGTCGCATC
CGGGCCTTGTTCCTGCTCCACTGGCCCAAATCAGTACCCCAATGTTCTTGCCTTCTGCC
CAAGTTTCATATCCTGGAAGGTCTTACATGCTCCAGCCACCTGTAAGCTATGCTAATG
TGCCAGAAACTCGTCCAGCCAGTGAGCTGCATCCAATGGCGTGTACCCATGTATTGCAC
AAGGAGTGTATCAAATTTGGGCCAGACCAACACAAATGACACTTGTCCCTTTGTCCA
ACTCTTAAA
    
```

Restriction Sites: Please inquire

ACCN: NM_001077239

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001077239.1</u> , <u>NP_001070707.1</u>
RefSeq Size:	2695 bp
RefSeq ORF:	2112 bp
Locus ID:	257160
UniProt ID:	<u>Q8ND24</u>
Cytogenetics:	11q23.3