

## Product datasheet for **SC315614**

### EXOC3L4 (NM\_001077594) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	EXOC3L4 (NM_001077594) Human Untagged Clone
Tag:	Tag Free
Symbol:	EXOC3L4
Synonyms:	C14orf73
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC315614 representing NM\_001077594.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGCCATACCACAGACAGACTCTGGGCCGGAGCTGCAGAGTCCAAGGAGGCTGAGGAGCCACAG
ACTCCAGCTCAGGGCTCCCGCGAACAAGCAGCAGGAAAGAGCCCAATGCCACCACAAGGATGGCACA
AGGCTGGGCTGGGCTCCCTGAGGCAGGCTTCTCCCGGCCAGCCAGCGGGCTTTGACCCAGGCTCC
AAGGAAGATACGGGCTGTTCCGGCGAAGCTCCTGCTCCCTGTTCCGGTCTTCCGGCAAGCCCTGAAT
GACGGCCAGCTACCGGCCATTCCAGGCCACTCCTGAGGTGCCCTCGGGGTCATGAATGGTGTGACG
CAGCAGGCATCCACTGGGCAGCGTCTGAGGAAGTAAACCCGAGGCAGAAGGCAAAATCCGTGGCCGAC
CTCATTACTGAACGGCAACTGCTGGCGGCTTTCGAACAGCTTCTGCGCTGGAGACGCTGCTGGTGGCC
GAGAAGGCCTCGCGCACCTTTGAGCAGGACCTACGGCCTTCGCGCGGCGCCTATGGACGTGTGCCTG
CTTACGACGGGCTGGCAGCCGAGATCGGCGCCATCGTGCAGGAGACGCTGGACAGCGACGGTGTGGAC
GCGGCCGCGCTGGCCGAGCTGGCCGCGTGGTGCAGCGGAGGAGGAAGCCACCCTTCTCCCCCGAC
GACGGCGACTTCTGCGCAGCCGCGCCCTGGCGCCAGCACTGGGAGGAGGCGGTGCGGCGAAGCGCT
CAGGAGCGCTGCGGCGGCCGGGCGGGTGGGCCCTTCGGGAGGCGGAGGGCGCGTCCGGTTGGCC
CAGCTTCTGGCCGAGCTGGTGGCTTGGTTCGCCGCGACCTGCAGAAGGTGCGGCAGGAGGTGCAGCCC
GCGTATGCGGCGGCCGGCTTCCAGCGTGGGAGGTCTATCTGCGTGCCTTCCACAGCGCCGTGGCCAG
CGCTCCAGGAGCTCGCGCGCAGCAGCCGCGGCTGCGAGCAGCTCTATCTGCTGGACTGGGCCGCC
AACGTCTACGGCAGTCTGACTTCTGGGCGCCCCGGGCTGGCGCTGCCCGCAGCCGCTGCCTCCG
CTCCTGGCGCCGACGTGTGGCCCGACTGGAGAGCGACTACACCAGTCTCTGGAGGCCAAGATCGCA
AGCTGTTGACAGCATTTGCAGCTGGAGCAGAGTCACTGGGCGGCCCGAGGTCCTCCGAGGTGCTG
CAGGGCCTCTACCAGGCGCCGCTGCTCATGGACGTCCATATGCTCGTGGCCGAGCACGTGAAGGCGGCC
GGCGCCATCTCCGCGGAGCTGGAGGCCACCACCCTGCGAATCTGCACGCGGGCGCTCGGCCCTTCTGTG
CCCAGGTTTAAAAGGCTTTTCTGGCGTCGGAGGCGGTGAGCGAGCCGACCTGGGCGCCTACATCAAC
GCCTGCGAGGAGCTCAGGACCAGTCTTCTCTCCAGGTTCCAGGAACCCAAGAGGAGCTGGAGAAGCCC
CTGGTGACGGCCACCTGCAGTTCAGAAAGCACTTGCTTCAGGGCTTGCAGCGTGAGTTGCAGCCGCTC
TTCAGGGTTGTGTGACCCAGGACTGGCTGACGAGGACTGGCTGCATCCCTCATGGACAAGGTGGTG
ACCTTCGCCGGTATCTCCAGCGTGTGGCCGGCCGCGGGCACAGGAGACTCTGCAGGAGGTGCACCGG
TTCGTGGTCCGCGAGTACCTGGCGGGGCGCTGAGGCCACGGGAGCGGTTCCGGGCGATGGAGCGCATG
CATGGCTCCCAGAAGATGAGCCTGGATGCCAGGCCATCAGCGACACCTTCCAGGGCCTGGGTCCGAG
GCCACATGGTTGGACCAAGCCATCCAGTGCCTGGTGGATCCTGGGCGAGACCTACAAAGATGACATC
CAGCGGCACCTGGAGACTCTTATCCGGAGTACCCCGACATCAGGCGGGACCACATACTGGCCATTCTG
GCGCTGCGCCGACTGGGCCGCCAGCGGAACCAGCATCTTTCAGCACACTCAAGACCTGCTGAGAGCT
GCGGCCGGGGCGGGTGGCGAGGCCCTCGGGCCGCGTCTTTCGAGGAGATCAAGGTGCCCAGT
GCCATGGCTGTGCTGATCACCTGCGTTAG
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001077594

**Insert Size:** 2169 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001077594.1](#)

**RefSeq Size:** 2591 bp

**RefSeq ORF:** 2169 bp

**Locus ID:** 91828

**UniProt ID:** [Q17RC7](#)

Cytogenetics: 14q32.32

MW: 79.9 kDa