

Product datasheet for **SC315617**

HNRNPUL2 (NM_001079559) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HNRNPUL2 (NM_001079559) Human Untagged Clone
Tag:	Tag Free
Symbol:	HNRNPUL2
Synonyms:	HNRNPUL2; SAF-A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315617 representing NM_001079559.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGTGAAGCGGCTGAAAGTGACCGAGCTGCGGTGCGAGCTGCAGCGCGGGGCTGGACTCGCGC
GGCCTCAAGGTGGATCTGGCGCAGCGGCTGCAGGAGGCGCTGGACGCCGAGATGCTCGAGGACGAGGCC
GGCGGCGCGGGGCCGGGCCGGCGGGGCTGCAAGGCGGAGCCTCGGCCGTGTGGCCGCGTCGGGCGGC
GGCCCGGGCGGGGACGAGGAGGAGGACGAAGAGGAGGAGGAGGAGGACGAGGAGGCGCTGCTTGAGGAC
GAGGACGAGGAGCCACCCCTGCTCAAGCCTTGGGTGAGCGCGCAGCCGCCCGGAGCCCCGGAG
GCGGCAGCCATGGAGGCCGCGGCCGAGCCAGATGCTTCCGAGAAGCGCGGAGGCCACGGCCGGGTCA
GGCGGGGTAAATGGTGCGAAGAGCAGGGCCTCGGCAAGAGGGAGGAAGACGAACCCGAGGAGCGGAGC
GGGACGAGACGCCGGGATCCGAGGTGCCGGTGACAAGGCCCGCAGGAACAGGAGATGACCAGGAT
AGTGAAAAGTCAAACCAGCAGGCTCAGATGGTGAGCGCGGGGGTAAAGAGACAGCGGATGAGAAG
GATGAACATGGCCGAGCTTACTATGAATCCGAGAGGAGGCTTACCACAGCCGCTCAAAGTCTCCACTG
CCTCCTGAAGAAGAGGCAAAAGATGAGGAGGAGGATCAAACCTCTGTGAACCTGGACACGTATACCTCG
GATCTGCATTTTCAAGTGAGCAAGACCGCTATGGAGGGCAGCCACTTTTCTCAGAGAAGTCCCCACC
CTTTGGTCTGGGCAAGGAGTACTTACGGAGTGACAAAGGAAAAGTCTGCTTTGAGGCAAGGTAACC
CAGAATCTCCCAATGAAAGAAGGCTGCACAGAGGTCTCTCTCCTTCGAGTTGGGTGGTCTGTTGATTTT
TCCCGTCCACAGCTTGGTGAAGATGAATTTCTTACGGTTTCGATGGACGAGGACTCAAGGCAGAAAAT
GGACAATTTGAGGAATTTGGCCAGACTTTTGGGAGAAATGATGTTATTGGCTGCTTGCTAATTTGAG
ACTGAAGAAGTAGAATTTCTTCTCCAAGAATGGAGAAGACCTAGGTGTGGCATTCTGGATCAGCAAG
GATTCCTTGGCAGACCGGGCCCTTCTACCCCATGTCCTCTGCAAAAATTTGTTGTAGAATTAATCTTC
GGTCAGAAGGAGGAGCCCTTCTTCCACCACCAAGAAGATTTGTGTTTATTTCATGCTGTGCTGTTGAG
GAGCGTGTACGCACTGCAGTCCCTCCCAAGACCATAGAGGAATGTGAGGTGATTCTGATGGTGGGACTA
CCCGGATCTGAAAGACCCAGTGGGCACTGAAATATGCAAAAGAAAACCTGAGAAAAGATACAATGTC
CTGGGAGCTGAGACTGTGCTCAATCAAATGAGGATGAAGGGTCTCGAGGAGCCAGAGATGGACCCAAA
AGCCGAGACCTTTTAGTTCAGCAAGCCTCCAGTGCCTTAGTAAGCTGGTCCAGATTGCTTCCCGGACA
AAGAGGAATTTATTCTTGATCAGTGAATGTGTACAATTCTGGCCAACGGCGGAAGCTATTGCTGTTT
AAGACCTTCTCTCGAAAGTGGTGGTGGTGTCCCTAATGAGGAAGATTGGAAGAAGAGGCTGGAGTTG
AGGAAGGAAGTAGAGGGAGATGATGTGCCTGAATCTATAATGCTGGAGATGAAAGCCAATCTCTTTTG
CCTGAAAAATGCGACTATATGGATGAGGTGACATATGGGAGCTGGAGAAGGAGGAAGCTCAGCCCAT
GTCACTAAGTACAAGGAGGAGGCAAGGAAGCTTCTGCCCCCTCCGAGAAGCGGACAAATCGCCGAAAC
AACCAGAAACAAGCGTAACCGGCAGAACCGAAGCGGGGCCAAGGCTATGTGGGCGGGCAGCGCCGAGGC
TACGACAACCGGGCTACGGGCAGCAGTACTGGGGCAGCCTGGAACAGAGGGGGTTACCGTAATTTT
TATGATCGATACAGGGGAGACTATGATCGATTTTACGGGCAGATTATGAGTACAACAGATACAGAGAC
TATTACAGACAATACAATCGGGATTGGCAGAGTTACTACTACCACCACCCCAAGGACAGAGACCGGATAC
TACAGGAATTACTACGGGTACCAAGGGTATCGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001079559

Insert Size: 2244 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001079559.2</u>
RefSeq Size:	5156 bp
RefSeq ORF:	2244 bp
Locus ID:	221092
UniProt ID:	<u>Q1KMD3</u>
Cytogenetics:	11q12.3
Protein Families:	Druggable Genome
MW:	85.1 kDa