

Product datasheet for SC315668

BEX4 (NM_001080425) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BEX4 (NM_001080425) Human Untagged Clone
Tag:	Tag Free
Symbol:	BEX4
Synonyms:	BEXL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315668 representing NM_001080425. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGAGTCCAAGAGGAAGTACGCGCAACAATCTCAACGGGGAAAAATGCCCAACAAGAAACGAAGGA
 GGGGAGCAGGCCCCACGCAGAAATGAAGAAGAATCCCGCATTTGGGAGGGGGTGAAGGCCAGAAGCCT
 GGAGGAAATATCAGGCGGGGCGAGTTAGGCGACTTGCCCTAATTTTCGATGGGCCATACCTAATAGG
 CATATTGAGCACAATGAAGCGAGAGATGATGTAGAAAGGTTTGTAGGGCAGATGATGGAATCAAGAGA
 AAGACTAGGGAACAGCAGATGAGGCACTATATGCGCTTCCAACTCCTGAACCTGACAACCATTATGAC
 TTTTGCCTCATACCTTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001080425
Insert Size:	363 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080425.3</u>
RefSeq Size:	1282 bp
RefSeq ORF:	363 bp
Locus ID:	56271
UniProt ID:	<u>Q9NWD9</u>
Cytogenetics:	Xq22.1
MW:	14.1 kDa
Gene Summary:	<p>This gene is a member of the brain expressed X-linked gene family. The proteins encoded by some of the other members of this family act as transcription elongation factors which allow RNA polymerase II to escape pausing during elongation. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Aug 2011]</p> <p>Transcript Variant: This variant (2) contains an alternate non-coding exon in the 5' UTR. Variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>