

Product datasheet for **SC315743**

GSG1 (NM_001080555) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	GSG1 (NM_001080555) Human Untagged Clone
Tag:	Tag Free
Symbol:	GSG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315743 representing NM_001080555. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCGATCCCTCTCAACTGACTCAAAATGTTTGCCTCACCCAGGAGATGGAGCTCTCGAAGGCCTTC
TCTGGCCAGCGGACACTCCTATCTGCCATCCTCAGCATGTATCACTCAGCTTCTCCACAACATCCCTG
CTCAGCAACTACTGGTTTGTGGGCACACAGAAGGTGCCAAGCCCTGTGCGAGAAAGGTCTGGCAGCC
AAGTGCTTTGACATGCCAGTGTCCCTGGATGGAGATACCAACACATCCACCCAGGAGGTGGTACAATAC
AACTGGGAGACTGGGGATGACCGGTTCTCCTCCGGAGCTCCGGAGTGGCATGTGGCTATCCTGTGAG
GAAACTGTGGAAGAACCAGCACTGCTCCATCCCAAGTCCCTGGAAACAATTTAGAGCCCTTCGGTCCAGT
GGTACAGCGGCAGCAAAAGGGGAGAGGTGCCGAAGTTTCATTGAACTTACACCACCAAGAGAGAGAA
ATCCTATGGTTATCCCTGGGAACGCAGATCACCTACATCGGACTTCAATTCATCAGCTTCTCCTGCTA
CTAACAGACTTGCTACTCACTGGAACCCCTGCCTGTGGGCTCAAAGTGGAGCCTTTGCTGCTGTTCC
TCTGTCTGTCCAGGTCTCCTGGGATGGTGGCCACATGATGTATTACAAGTCTTCCAAGCGACTGTC
AACTGGGTCCAGAAGACTGGAGACCACATGTTTGAATTATGGCTGGGCCTTCTACATGGCCTGGCTC
TCCTTACCTGCTGCATGGCGTCGGCTGTACCACCTTCAACACGTACACCAGGATGGTGTGGAGTTC
AAGTGCAAGCATAGTAAGAGCTTCAAGGAAAACCCGAAGTGCCTACCACATCACCATCAGTGTTCCT
CGGCGGCTGTCAAGTGCAGCCCCACCGTGGGTCTTTGACCAGTACCACCAGTATCATAATCAGCCC
ATCCACTCTGTCTGAGGGAGTCGACTTCTACTCCGAGCTGCGGAACAAGGGATTTCAAAGAGGGGCC
AGCCAGGAGCTGAAAGAAGCAGTTAGGTCATCTGTAGAGGAAGAGCAGTGTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001080555

Insert Size: 1089 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080555.2](#)

RefSeq Size: 2637 bp

RefSeq ORF: 1089 bp

Locus ID:	83445
Cytogenetics:	12p13.1
Protein Families:	Transmembrane
MW:	40.7 kDa
Gene Summary:	<p>May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and 5' coding region, includes an alternate in-frame exon in the central coding region, and uses an alternate splice site that causes a frameshift in the 3' coding region, compared to variant 1. The encoded isoform (4) has distinct N- and C-termini and is longer than isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>