

## Product datasheet for **SC315777**

### ZNF275 (NM\_001080485) Human Untagged Clone

#### Product data:

|                           |                                                                                            |
|---------------------------|--------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                        |
| Product Name:             | ZNF275 (NM_001080485) Human Untagged Clone                                                 |
| Tag:                      | Tag Free                                                                                   |
| Symbol:                   | ZNF275                                                                                     |
| Mammalian Cell Selection: | Neomycin                                                                                   |
| Vector:                   | pCMV6-Entry (PS100001)                                                                     |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                       |
| Fully Sequenced ORF:      | >SC315777 representing NM_001080485.<br>Blue=Insert sequence Red=Cloning site Green=Tag(s) |

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG  
GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**  
ATGATGAGTCATCCGTGTGTCTCTTTTGGGCGTTCTGTTTTAAATCCTGCCTTGGTCCCTCACCTG  
GCACAAGGACAAGTGCTACTGGTGTGACACCCATCGCCAACACTGATCTGCTAAGTACTCTGAAAGC  
ACCTCCGCGACCCGACACCAGATGAAGGGGAAGATGCCAGCCACAGGAGATGGCGTCCACAAGCTTC  
CCAAGGGCCAGTGGTCCCAGTCTGAATTCAGACAGCACGGGGACTCTGACGGGAAGAGAGGGAGCCCA  
CAGAATCTGCCATAGAATCATTTTGTCTGTAAGAGTGTGGGACACCTTTCCGCTTAAAGTCTGT  
CTTGTCAGCACCAGAGAGTCCACAGTGAGGAGAAGGGCTGGGAATGTGGCGACTGCGGGAAGGTCTTT  
AGGGGGGTGGCGAGTTTAATGAGCACAGGAAAAGCCACGTAGCTGCGGAGCCCCAGCCCGGCCCCAGT  
AGGGCCCTGGAGAATGCCGCGGAGAAGAGGGAGCAGATGGAGAGGGAGGCAAAGCCCTTCGAGTGCAG  
GAGTGCAGAAAACGGTTTAAGAAGAATGCAGGCCTGAGTCAGCATCTGAGGGTCCACAGCAGAGAGAAG  
CCCTTTGATTGCGAGGAATGCGGGCGGTCTGTTCAAAGTCAACACCCACCTTCCGACATCAGAACTT  
CACACTTCGAAAAGCCTTTCGCTGCAAGGCGTGCAGCAGGGATTTCTGGATCGCCAGGAGCTTCTC  
AAGCACCAGCGCATGCACACTGGCCACCTGCCCTTCGACTGCGACGACTGCGGCAAGTCTTCCGAGGG  
GTCAACGGGTGGCCGAGCACCAGCGCATCCACAGTGGGGCCAAGCCATACGGGTGTCCCACTGCGGC  
AAGCTCTTCCGAAGGAGCTCGGAGCTACCAAGCACCGCGGATCCACACGGGCGAGAAGCCCTACGCC  
TGGCGCCTTCGTGGGCCCTC**TGA**  
**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC

|                    |              |
|--------------------|--------------|
| Restriction Sites: | SgfI-MluI    |
| ACCN:              | NM_001080485 |
| Insert Size:       | 990 bp       |



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001080485.3</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 6320 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 990 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 10838                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Cytogenetics:</b>          | Xq28                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>MW:</b>                    | 37 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Gene Summary:</b>          | This gene encodes a zinc finger protein that appears to be conserved in eutheria. Its function has not yet been established. [provided by RefSeq, Jul 2010]                                                                                                                                                                                                                                                                                                                                                             |