

Product datasheet for **SC315778**

CPEB1 (NM_001079534) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CPEB1 (NM_001079534) Human Untagged Clone
Tag:	Tag Free
Symbol:	CPEB1
Synonyms:	CPE-BP1; CPEB; CPEB-1; h-CPEB; hCPEB-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC315778 representing NM_001079534.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTTTCCCAACCTCTGCGCAAGAATCTTCCCGTGGCCTCCCAGATGCAAATGACTTGTGCCTTGGC
CTGCAGTCCCTCAGTCTGACAGGCTGGGACCGACCCCTGGAGCACCCAGGACTCAGATTCTCAGCCCAG
AGCAGCACACACTCGGTACTGACATGCTCCATAACCCACTGGGAAATGCTCTAGGAAAACCCCTTG
AGTTCCTGCCTCTGGATCCCCTTGGGTCTGACTTGGTGGACAAGTTCCAGCACCTCAGTTAGAGGA
TCACGCCTGGACACCCGGCCATCCTGGACTCTCGATCTAGCAGCCCTCTGACTCAGACACCAAGTGGC
TTCAGCTCTGGATCAGATCATCTCTCAGATTTGATTTCAAGCCTTCGATTTCTCCACCTCTGCCCTTC
CTGTCTCTGTCAAGGGTGGTCCCAGAGACCTTTAAAGATGGGGTAGGGTCTCGGATGGACCAAGAG
CAAGCTGCTCTTGTGCAGTCACTCCCTCCCAACCAAGTCTTCAAGAGATGGCCAGGAGCTTCTGTG
TGGCCATCTGGACCTCTCGAAGCTCCAAAGACCCCTTCAGCATAGAGAGAGAGGCCAGGCTGCAC
CGACAAGCTGCAGCTGTGAATGAAGCCACTGTACCTGGAGTGGCCAGCTTCTCCCGGAACTATAAG
AACCCATCTACTCTTGAAGGTGTTTCTAGGAGTGTTCTTGGGATATTACAGAAGCTGGATTAGTT
AACACCTCCGTGTTTTGGCTCTTTGAGTGTGGAGTGGCCTGGTAAGGATGGCAAGCATCCCCGGTGT
CCTCCAAAGGTATGTGTATCTGGTCTTGAAGTACGAGAAAGTCTGTCCGATCCTTGTTCAGGCTTGC
TCTCATGACCCGCTGAGCCAGATGGCCTGAGTGAATATTATTTCAAGATGTCCAGCCGAAGGATGCGC
TGCAAGGAGGTGCAGGTGATCCCTGGGTATTAGCCGACAGTAACCTTGTCCGGAGCCATCTCAGAGG
CTTGACCCAGCAGGACGGTGTGTGGTGTCTGCATGGAATGCTAAATGCTGAGGCCCTGGCAGCC
ATCTTGAACGACCTATTTGGTGGAGTGGTGTATGCCGGGATTGACACAGATAAGCACAAGTATCCATT
GGTCTGGTCTGTGACTTTCAATAACCAACGGAGTTACCTGAAAGCAGTCAGCGCTGCTTTGTGGAG
ATCAAAACCAACAAGTTCACAAAGAAGTTCAGATTGACCCCTACCTAGAAGATTCTCTGTGTCATATC
TGCAAGTCTCAGCCTGGTCTTTCTTCTGTGAGATCAGGTCTGCTTCAAATACTTCTGCCGGAGCTGC
TGGCACTGGCGGCACAGCATGGAGGGCCTGCGCCACCACAGCCCTGATGCGGAACCAGAAGAACCGA
GATTCCAGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001079534

Insert Size:	1461 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001079534.1
RefSeq Size:	3095 bp
RefSeq ORF:	1461 bp
Locus ID:	64506
UniProt ID:	Q9BZB8
Cytogenetics:	15q25.2
Protein Pathways:	Dorso-ventral axis formation, Oocyte meiosis, Progesterone-mediated oocyte maturation
MW:	53.6 kDa
Gene Summary:	<p>This gene encodes a member of the cytoplasmic polyadenylation element binding protein family. This highly conserved protein binds to a specific RNA sequence, called the cytoplasmic polyadenylation element, found in the 3' untranslated region of some mRNAs. The encoded protein functions in both the cytoplasm and the nucleus. It is involved in the regulation of mRNA translation, as well as processing of the 3' untranslated region, and may play a role in cell proliferation and tumorigenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (4) contains a distinct 5' UTR and lacks an in-frame portion of the 5' coding region, compared to variant 1. The resulting isoform (3) has a shorter N-terminus, compared to isoform 1. Variants 3, 4, and 5 encode the same isoform.</p>