

Product datasheet for **SC315788**

LRRC14B (NM_001080478) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC14B (NM_001080478) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC14B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315788 representing NM_001080478.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACACAATGAGGTCACTCCGCTTCATTTCTGCAGAAGCTCTGGTGTCCACCCCAAGTGGCCCGG
CAGAGCCTGGACAGCGTGGCCCAACCTCTACCACTCCTGTTCAAAGCCAGCTACCTGCTGGAGCAG
GCGGAGGTGACGCGCGCGGTGCTGGGGCGCTGGCCCTGGAGGAGTTCCGGCTGGGAGCGCTGCTGGGT
CCTGGTGGCGACCAACCCCAAGGACCTGCGCGACAGAACCTGCAGGGCTGCCTGGAGGCGCTGGTGGCG
GGCCTCGCGGACCACGTGCTGCAGGACCGAGCGCGCGCGGTGCGGGTGGCTGACCTCACGGGCATC
CGAGATGTGCAGGTGCAGCGGTGCCGTGCGGGAGGGCGCTGGGAGGTGGGGCCGACCCAGCTGCTG
GCCAGGACCTGCTGTGAGCTGCAGGCAGAGCCCTCGAGCCGGGCGCCCGTCGAGGTCTCGCCGAC
CTCTTCGTCAGTGGGCAACTTCGAGGCGGTGGTGCAGGCTCTGAGGCCAGCGGGCCCGGCCCTCTG
CGGGTGCAGTCCCTCTGTTCCGGGCGGACAGCCTGAGCCCAAGCAGCTCCTGCAGTGTGCTGCTG
GCTGGCCCGGTGCCCTGCGCAAGCTGGAGGTGGTGCACAACGTGCGGCTGCATGCGGGCCACGTGCAG
CAGCTTCTGGCCAGGTGGGCTTCCCCGGCTGGCCTCGCTCACCTGCCACCAAGGCCTTTGATGCA
CCCCCACCTACGCTCCACTCCGACGGCGAGGACCCCTCCTCGCTCCATCGCCGGGAGCTCAGC
AAGATGGCGCAGCTCACTGAGCTCAGTGTGGCCTTCTCCACGCTGACCGGAAGATCCCGACGCTGCTT
GGCCCCCTACAGACCCCGCTGCGAGTACTGGACCTGGCCAAGTGTGCCCTGAACACACGGACATGGCC
TTCTTGGCAGACTGTGCCACGCTGCCACCTGGAGGTGCTGGACCTCAGTGGACACAACCTGGTCAGC
CTGTACCCCTCGACCTTCTTCAGGCTGCTCAGCAGGCTTCCCGACGCTGAGGATCCTGACACTGGAG
GAGTGTGGCATCGTAGACAGCCACGTTGGCATGCTGATCCTGGGCTGAGCCCTGCCACCGGCTGCGC
CAGCTCAAGTTCTCGGGAACCGCTGTCGGCCGCGCCCTCCGGCGCTCTTCAACGCACTCTGTGAG
CTCCCCGAGCTGCGCTGCATTGAGTTCCCGGTGCCCAAGGACTGCTACCCGAGGGTGCCGCTACCCA
CAGGACGAGCTGGCCATGTCCAAGTTCAACCAGCAGAAATACGACGAGATCGCCGAGGAGCTGCGTGCC
GTGCTGCTGCGGGCTGACCGAGAGGACATCCAAGTCTCCACACCTCTCTTTGGAAGTTTTGACCCAGAC
ATTCAAGAAACAAGCAATGAGCTTGGTGTCTTCTGCTGCAAGCTTTCAAACGCTCTAGAAAACCTTC
TCCAGAGCACTCAAACAATAGAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001080478

Insert Size: 1545 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080478.1</u>
RefSeq Size:	1573 bp
RefSeq ORF:	1545 bp
Locus ID:	389257
UniProt ID:	<u>A6NHZ5</u>
Cytogenetics:	5p15.33
MW:	56.8 kDa
Gene Summary:	The protein encoded by this gene is a leucine-rich repeat containing protein that is a member of the PRAME family. [provided by RefSeq, Apr 2017]