

Product datasheet for **SC315797**

ARHGAP9 (NM_001080156) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP9 (NM_001080156) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP9
Synonyms:	10C; RGL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315797 representing NM_001080156.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCAGAGCCCCCTGTGTACTGTAACCTGGTGGACCTTCGCCGCTGTCCTCGGTCCCCACCCCCAGGC
CCTGCATGCCCCCTGCTGCAGAGGCTGGATGCCCTGGGAGCAGCACCTGGACCCCAACTCTGGACGCTGC
TTCTACATAAATTCAGTACTGGCTGCAAGTCTTGGGAAGCCCCCGCCGAGTCGCAGCGAGACGAAC
CCTGGCTCCATGGAGGGACACAGACCCTGAAGAGGAACAATGATGTCCTGCAACCTCAGGCAAAGGGC
TTCAGATCTGACACAGGGACCCAGAACCCTTGACCCACAGGGTTCAGTACGCTCAGCCAACGCACC
TCGCAGCTTGACCCTCCAGCCTTGACGGCCCTCGACCTCGCCGAGCTCCTGGACGACCCCATGAG
GTGAAAAGTCCGGTCTGCTCAACATGACCAAGATTGCCAAGGGGGGCGCAAGCTCAGGAAGAAGTGG
GGCCCGTCTTGGGTGGTGTAAACGGGTAAACAGCCTGGTGTCTACCGAGAGCCACCGCCGACAGCGCCC
TCCTCAGGCTGGGACCAGCGGGTAGCCGCCCCAAAGTAGCGTGGACCTGCGCGGGGCGGCCCTGGCG
CACGGCCGCCACCTGTCCAGCCGCGCAACGTCTGCACATCCGCACGATCCCTGCCACGAGTTCCTG
CTGCAGTCGGACCACGAGACAGAGCTGCGAGCCTGGCACCAGCGCCTGCGGACTGTCATCGAGCGGCTG
GATCGGGAGAACCCCTGGAGCTGCGTCTGTGGGCTCTGGACCCGCGGAGCTGAGCGCCGGGAGGAC
GAAGAAGAGGAGTCGGAGCTGGTGTCCAAGCCGCTGCTGCGCCTCAGCAGCCGCCGAGCTCCATTCGG
GGCCCCGAAGGCACCGAGCAGAACCCTGCGCAACAACTAAAGCGGCTCATCGGAAGAGACCGCCC
TTACAAAGCCTGCAGGAGCGGGTCTGCTCCGAGACCAGGTGTTCCGGTCCAGTTGGAATCACTCTGC
CAGCGGAAGGAGACACGGTGCCAGCTTTTTGCGGCTCTGCATTGCTGCTGTGGATAAAAGAGGTCTA
GATGTGGATGGCATTATCGGGTGAAGCGGAACCTGGCAGTGGTCCAGAAGCTTCGCTTTCTGGTGGAC
AGAGAGCGTGCAGTCACTCCGATGGGAGTATGTGTTCCAGAAGAGGAGGACAAAGAAGTCCGGTTA
GATTTGGACAGTACTGAGTGGATGACATTCATGTGGTACCGGAGCCCTGAAGCTTTTTCTCCGGGAG
CTGCCACGCTCTGGTCCACCCTGCTGCTGCCCATTTCCGTGCTGCCCTTGCAGTCTCCGAATCA
GAGCAGTCCCTCTCAGATAACAAGATTAATAGGCTCAATGCCAAAGCCCAACCATGACTCTACGG
TACCTCTGGAGCATTTATGCAGGGTATAGCACACTCAGATAAGAATCGCATGACACCCCAACCTG
GGAATTGTTTGGACCAACCCTGTTTCGGCCAGAGCAGGAGACATCTGACCCAGCAGCCCATGCTCTC
TACCCAGGCAGCTGGTCCAGTGTGCTACCAACTTACCAGCCTTTCCTCCGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001080156
Insert Size:	1644 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080156.2</u>
RefSeq Size:	2355 bp
RefSeq ORF:	1644 bp
Locus ID:	64333
UniProt ID:	<u>Q9BRR9</u>
Cytogenetics:	12q13.3
MW:	61.2 kDa
Gene Summary:	<p>This gene encodes a member of the Rho-GAP family of GTPase activating proteins. The protein has substantial GAP activity towards several Rho-family GTPases in vitro, converting them to an inactive GDP-bound state. It is implicated in regulating adhesion of hematopoietic cells to the extracellular matrix. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region, initiates translation at a downstream start codon, compared to variant 1. The resulting isoform (3) has a shorter N-terminus compared to isoform 1. Variants 3 and 11 both encode the same isoform (3).</p>