

Product datasheet for **SC315800**

LRRIQ4 (NM_001080460) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRIQ4 (NM_001080460) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRIQ4
Synonyms:	LRRC64
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315800 representing NM_001080460.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCAAAGACATAAAATCAGTAGAACATTCACCTAAAATTCATCAGAGAAATGATCCACAGCACGTC
AATGATAGAACATTTTTTCATTGATGCCTAATCAGAGCTTGACTGCCATTCTTTGGAGATCTTCACA
TTCACAGAATTAGAAGAAGTGCATTTGGAGAATAACCAGATTGAAGAAATTCGCCAGGAGATTCAGCGT
TTAAAGAACATCAGGGTCTCTACCTGGATAAGAACAACCTGAGGAGCCTGTGCCCGCGCTGGGGCTG
CTGAGCAGCCTGGAGAGCCTGGACCTGAGCTACAACCCCATCTTCTCCTCCTCCCTTGCTGTTGTCAGC
TTCCTCCACGCCCTGCGCGAGCTCCGGCTACCAGACCGACCTGAAGGAAATTCCTCGTCATCTTT
AAAAACCTCCACCATCTCGAGCTGCTCGGACTGACCGAAACACCTGAAATGTCTGCCCAAGGAAATA
GTGAACCAGACCAAGCTGAGGGAGATCTACCTGAAGCGAAACAGTTTGAAGTTTTCCCGAGGAGCTC
TGTGTTCTACACCCTGAAATCATTGACCTGGACGAGAAACAAATAGGTGCCATCCAGAAGAGATC
GGACACCTGACGGGGCTGAGAAGTTCTATATGGCTTCTAACACCTTCCCGTTCTGCCCGCTCCTTG
TGCCAGTGTAGCCAACTGTCGGTCTCGATTTATCCACAACTCCTCCACTCCATCCCGAAGAGCTTC
GCCGAGCTCAGGAAGATGACGAAATCGGGCTGAGCGGGAACCGCCTGGAGAAGGTGCCACGCCTCATT
TGCAGGTGGACCTCGCTGCACCTGCTCTACCTGGGAAACACCGCCTGCACAGGCTGCGGGGCTCCTTC
AGGTGCCTGGTCAACTTGCCTTCTGGACCTAAGCCAGAACCATCTGCACCACTGCCCGCTGCAGATC
TGTGCACTGAAGAACCTTGAAGTCTGGGACTGGATGACAATAAAATAGGACAGTTACCTTCAGAATTG
GGCTCACTTTCAAACCTGAAGTACTTGGACTAACAGGAAATGAGTTCCTTTCTTTCCGGAGGAAGTC
CTTTCTTAGCGTCTTTAGAGAAATTATACATTGGGCAAGACCAGGGATTCAAACCTACCTATGTGCCA
GAACACATTAGGAACTGCAGAGTCTCAAAGACTATATATAGAGAACAATCATCTGGAGTACCTGCC
GTATCCTTGGGGTCAATGCCTAACCTAGAAGTTCTTGATTGCCGACAAATTTGCTTAAGCAACTTCCA
GATGCCATTTGCCAAGCACAAGCTTTGAAAGAATTACGGCTGGAGGACAATTGCTCACCATCTTCCA
GAGAATTTGGATTCCCTAGTGAATCTTAAGTTCTGACTGATGGACAATCCCATGGAAGAACCCCA
AAAGAAGTGTGTCTGAAGGCAATGAGGCCATATGAAAATACCTCAAGGAAAACAGAAACAGGAATATA
ATGGCAACAAAGATTGAGCATGGTGGCGTGGAACAATGGTACAGAGAGGATTTGGGAAATTCGGTGAA
CTACTAAAACCACAAAAGAAAGGAAAGACCTCTCCAAAAGATAAGAAAGGAAAGAAGGATGTAAGGA
AAACCAGGAAAGGAAAAAGAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
  
```

- Restriction Sites:** SgfI-MluI
- ACCN:** NM_001080460
- Insert Size:** 1683 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080460.1](#)

RefSeq Size: 1825 bp

RefSeq ORF: 1683 bp

Locus ID: 344657

UniProt ID: [A6NIV6](#)

Cytogenetics: 3q26.2

MW: 64 kDa