

Product datasheet for SC315801

ZNF878 (NM_001080404) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF878 (NM_001080404) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF878
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001080404, the custom clone sequence may differ by one or more nucleotides

```

ATGCGCTCGGGGTTCTGTCCTATCTGGGTTGTACCTTTCAGAGGATGAAGGGGCAGGG
ACGGCCACCGAGTGAGGAAGACAACAGTGTCCCTTCTAGCAGTTCGGAGCAAGTTTGGC
CGAAGGCCAAATATTGCAAGTAAAGGATTCGGTGGCCTTTGAGGATGTGGCTGTGAACCTC
ACCCAGGAGGAGTGGGCTTTGCTGGATCCTTCCAGAAAAATCTCTACAGGGAAGTGATG
CAGGAAACCTTGAGGAACCTGACCTCCATAGGAAAAAATGGAACAACCACTACATTGAA
GATGAGCACCAAAATCCTAGGAGAAACCTAAGAAGACTTATAGGGGAGAGACTCTCTGAA
AGTAAAGAAAGTCATCAGCATGGAGAAGTTTGGACACAGGTTCCAGATGACACACTGAAG
AAGAAAACTCCTGGAGTACAATCATATGAAAGCAGTGTGTGTGGAGAAAATCGGCATAGGT
CTTTCATCCCTTAATAGGCACCTCAGAGCCTTTAGTTATTCCAGTAGCCTTGCAATACAT
GGAAGAACTCACACTGGGAAAAAGCCTTATGAATGTAAGGAATGTGGAAAGCATTTCAGG
TTTCCAGTTCTGTTCTGAGACATGAAAGAATCCACTCTGCAAAAAAACCTATGAATGT
AAGCAGTGTGGGAAAGCATTTCAGTTTCCAGTTCTGTTCTGAGACATGAAAGAATCCAC
TCTGCAAAAAAACCTATGAATGTAAGCAGTGTGGGAAAGCATTATCTTATCTTGAAGC
TTTCAGACACACATGAGAATGCACACTGGAGAGAGACCTATAAATGTAAACATATGTGGG
AAAGCCTTTTTTCTCCAGTTCGTTAAAAAGACACGAGAAAAGTCACACTGGAGAGAAA
CGCTATAAATGCAAGCAATGTGATAAAGCCTTCAATTGTCCAGTTCCTTTCAATATCAT
GAAAGGACTCACAGTGGAGAGAAACCTATGAGTGTACACAATGTAGGAAAGCCTTCAGA
TCTGTCAAGTACCTGCGAGTACATGAAAGAAAACACACTGGAGAGAAACCTATGAGTGT
AAGCTATGTGGTAAGGGATTTATTTCTTCCACTTCCTTTTCGCTATCATGAAAAGACTCAC
ACTGGAGAGAAACCTTATGAATGTAAGCAATGTGTGAAAGCCTTCAGTTTGTCAAGGAT
CTTCGAATACATGAAAGGACACACACTGGAGAGAAACCTTTGAATGTAAACAATGTGGG
AAAACCTTCACCTCTTCCAATTCCTTTCACTATCATGAAAGGACTCACACTGGAGAGAAA
CCCTATGAGTGTAAAGCAATGTGGGAAAGCCTTCAGATCTGCCTCAGTCTTCAAAAGCAC
ATACGAACTCACACAGGAGAGAGAAACCTATGGATGTAAGCAATGTGGTAAAGTCTTTAGA
GTTGCCTCACAACCTTAAATGCATGAAAGGACTCACACAGGAGAGAAACCTATGAGTGT
AAGCAATGTGGAAAGCCTTCATTTCTTCTAATTCTATTGCTATCATAAAGGACTCAC
ACTGGAGAGAAACCTATAAATGTAAACAGTGTGGGAAAGCCTTCATTTCTTCCAACCTCT
TTTCTCTACCATGAAAGGATTCATACTGGAGAGAAACCTATGAGTGTAAAGCAATGTGGT
AAAGCCTTTAGATCTGCCTCAATCCTTCAAAAGCATGTAAGGACTCACGCTGGC
  
```

Restriction Sites: Please inquire


[View online »](#)

ACCN:	NM_001080404
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080404.1, NP_001073873.1</u>
RefSeq Size:	1737 bp
RefSeq ORF:	1737 bp
Locus ID:	729747
UniProt ID:	<u>C9JN71</u>
Cytogenetics:	19p13.2
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]