

Product datasheet for **SC315802**

TTC39A (NM_001080494) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39A (NM_001080494) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39A
Synonyms:	C1orf34; DEME-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315802 representing NM_001080494.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCCAGAAGGGCCACAAAGACTCCCTTTACCCTTGTGGAGGGACTCCTGAGAGCAGCCTCCATGAG
GCCCTGGACCAGTGCATGACCGCCCTGGACCTCTTCCTACCAACCAGTTCTCAGAAGCACTCAGCTAC
CTCAAGCCCAGAACCAAGGAAAGCATGTACCCTCACTGACATATGCCACCATCCTGGAGATGCAGGCC
ATGATGACCTTTGACCCTCAGGACATCCTGCTTCCGGCAACATGATGAAGGAGGCACAGATGCTGTGT
CAGAGGCACCGGAGGAAGTCTTCTGTAACAGATTCCTTCAGCAGCCTGGTGAACCGCCCCACGCTGGGC
CAATTCCTGAAGAGGAAATCCACGCTGAGGTCTGCTATGCAGAGTGCCTGCTGCAGCGAGCAGCCCTG
ACCTTCTGCAGGACGAGAACATGGTGAGCTTCATCAAAGGCGGCATCAAAGTTCGAAACAGCTACCAG
ACCTACAAGGAGTGGACAGCCTTGTTTCAGTCTCACAATACTGCAAGGGTGAGAACCACCCGACTTT
GAAGGAGGAGTGAAGCTTGGTGTAGGGCCTTCAACCTGACACTGTCCATGCTTCTCTACTAGGATCTG
AGGCTGTTGGAGTTTGTGGGTTTTTCAGGAAACAAGACTATGGGCTGCTGCAGCTGGAGGAGGGAGCG
TCAGGGCACAGCTTCCGCTCTGTGCTCTGTGTCATGCTCCTGCTGTGCTACCACACCTTCCCTCACCTTC
GTGCTCGTACTGGAAACGTCAACATCGAGGAGGCCGAGAAGCTCTTGAAGCCCTACCTGAACCGGTAC
CCTAAGGGTGCCATCTTCTGTTCTTTGACGGGAGGATTGAAGTCATTAAGGCAACATTGATGCAGCC
ATCCGGCGTTTCGAGGAGTGTGTGAGGCCCAGCAGCACTGGAAGCAGTTCACCACATGTGCTACTGG
GAGCTGATGTGGTGTTCACCTACAAGGCCAGTGAAGATGTCCTACTTCTACGCCGACCTGCTCAGC
AAGGAGAACTGCTGGTCCAAGGCCACCTACATTTACATGAAGGCCGCCTACCTCAGCATGTTTGGGAG
GAGGACCACAAGCCGTTTCGGGGACGACGAAGTGAATTAATTCGAGCTGTGCCAGGCCTGAAGCTCAAG
ATTGCTGGGAAATCTCTACCCACAGAGAAGTTTGCCATCCGGAAGTCCCGCGCTACTTCTCCTCCAAC
CCTATCTCGCTGCCAGTGCCTGCTCTGGAAATGATGTACATCTGGAACGGCTACGCCGTGATTGGGAAG
CAGCCGAAACTCACGGATGGGATACTTGAGATTATCACTAAGGCTGAAGAGATGCTGGAGAAAGGCCCA
GAGAACGAGTACTCAGTGGATGACGAGTGTGGTGAATTTGTTGAAAGCCTGTGTCTGAAATACCTG
GGCCGTGTCCAGGAGGCCGAGGAGAAATTTAGGAGCATCTCTGCCAATGAAAAGAAGATTAATATGAC
CACTACTTGATCCCAAACGCCCTGCTGGAGCTGGCCCTGCTGCTTATGGAGCAAGACAGAAACGAAGAG
GCCATCAAATTTTGAATCTGCCAAGCAAACTACAAGAATTACTCCATGGAGTCAAGGACACACTTT
CGAATCCAGGCAGCCACACTCAAGCCAAGTCTTCCCTAGAGAACAGCAGCAGATCCATGGTCTCATCA
GTGTCCTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001080494

Insert Size: 1737 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080494.3](#)

RefSeq Size: 2743 bp

RefSeq ORF: 1737 bp

Locus ID: 22996

UniProt ID: [Q5SRH9](#)

Cytogenetics: 1p32.3

MW: 66.2 kDa