

Product datasheet for **SC315803**

ASPG (NM_001080464) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ASPG (NM_001080464) Human Untagged Clone
Tag:	Tag Free
Symbol:	ASPG
Synonyms:	C14orf76; GPA/WT; LYSOLP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315803 representing NM_001080464.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGCGCGCGTGGGGCCGAGCGGAGGCTGCTGGCCGTCTACACCGCGGCACCATTGGCATGCGG
AGTGAGCTCGGCGTGTGTGCCCGGACGGGCTGGCTGCCATCCTGAGGACACTGCCCATGTTCCAT
GACGAGGAGCACGCCGAGCCCGGGCTCTCTGAGGACACCCTGGTGTACCCCGGCCAGCCGAAC
CAGAGGATTCTCTACACCGTGTGGAGTGCCAGCCCTCTTCGACTCCAGTGACATGACCATCGCTGAG
TGGGTTTGCTTGCCAGACCATCAAGAGGCACTACGAGCAGTACCACGGCTTTGTGGTCATCCACGGC
ACCGACACCATGGCCTTTGCTGCCTCGATGCTGCTTTCATGCTGGAGAACCCTGCAGAAGACTGTCATC
CTCACTGGGGCCAGGTGCCCATCCATGCCCTGTGGAGCGACGGCGGTGAGAACCCTGCTGGGGGCACTG
CTCATGGCTGGCAGTATGTGATCCAGAGGTCTGCCTTTTCTCCAGAATCAGTGTTCGGGGCAAC
CGGGCAACCAAGGTAGACGCTCGGAGGTTGCGAGCTTTCTGCTCCCGAACCTGCTGCCTCTGGCCACA
GTGGGTGCTGACATCACAATCAACAGGGAGCTGGTGCGGAAGGTGGACGGGAAGGCTGGGCTGGTGGT
CACAGCAGCATGGAGCAGGACGTGGGCCTGCTGCGCCTTACCCTGGGATCCCTGCCGCCCTGGTTCGG
GCCTTCTTGACGCTCCCTGAAGGGCGTGGTCATGGAGACCTTCGGTTCAGGGAACGGACCCACCAAG
CCCGACCTGCTGCAGGAGCTGCGGGTGGCCACCGAGCGCGGCTGGTTCATCGTCAACTGTACCCACTGC
CTCCAGGGGGCTGTGACCACAGACTATGACAGTGGCATGGCCATGGCGGGAGCCGGCGTCATCTCAGGC
TTCGACATGACATCGGAGGCCGCCCTGGCCAAGCTATCGTATGTGCTGGGCCAGCCAGGGCTGAGCCTG
GATGTCAGGAAGGAGCTGCTGACCAAGGACCTTCGGGGGAGATGACGCCACCCTCGGTGGAAGAGCGC
CGGCCCTCACTGCAGGGCAACACGCTGGGCGGTGGGTCTCCTGGCTCCTCAGTCTGAGCGGCAGCCAG
GAGGCAGATGCCCTGCGGAATGCCCTGGTGCCAGCCTGGCCTGTGCTGCTGCCACGCCGCTGACGTG
GAGGCGCTGCAGGCGCTTGTGGAGCTGGGCAATGACCTGGGCTGGTGGACTTTAACGGCCAAACCCCA
CTGCACGCGGCCGCCCGGGAGGCCACACAGAGCAGTACCATGCTGCTGCAGAGAGGTGTGGACGTG
AACACCCGGGACACGGATGGCTTCAGCCGCTGCTGCTGGCCGTGCGGGGAGGCATCCGGGTGTCATT
GGGTGCTGCGGAAGCCGGGGCTCCCTGTCCACCCAGGAGCTGGAGGAAGCAGGGACGGAGCTGTGC
AGGCTGGCATACAGGGCCGACCTCGAAGGCTGCAGGTGTGGTGGCAGGCAGGGGCTGACCTGGGGCAG
CCGGGCTATGACGGGCACAGCGCCCTGCACGTCGAGAGGCAGCCGGGAACCTGGCAGTGGTGGCCTTT
CTACAGAGCCTGGAGGTGCGGTGGTGGCCAGGCCCCATGCCAGAAGTGTGCTGCTGGTGTCTAA
ACGGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001080464

Insert Size: 1722 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080464.2</u>
RefSeq Size:	1972 bp
RefSeq ORF:	1722 bp
Locus ID:	374569
UniProt ID:	<u>Q86U10</u>
Cytogenetics:	14q32.33
MW:	60.9 kDa
Gene Summary:	Exhibits lysophospholipase, transacylase, PAF acetylhydrolase and asparaginase activities. [UniProtKB/Swiss-Prot Function]