

Product datasheet for **SC315809**

Huntingtin Associated Protein 1 (HAP1) (NM_001079870) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Huntingtin Associated Protein 1 (HAP1) (NM_001079870) Human Untagged Clone
Tag:	Tag Free
Symbol:	HAP1
Synonyms:	HAP2; hHLP1; HIP5; HLP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315809 representing NM_001079870.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCGCCCGAAGAGGTTGGGCCGGTGTGCGCGGGGAGCCGGCTCGGACCCGGGGACCCAGCAGCACTC
ACCTGTGCACCTTCGCCCTCAGCCAGTCCCGCTCCGGAGCCCTCTGCGCAGCCGACGGCAGGGGCACT
GGACAGAGAGTAGGATCCCGAGCCACCTCTGGATCCCAGTTCCTCTCGGAAGCCCGACCGGAGCTCGC
CCGGCCTCGGAGGCTGGAGCCAAGGCAGGAGCCCGGCCCGCTCCGCATTCTCGGCCATCCAAGGGGAT
GTCCGGTCTATGCCCGACAATTCGGACGCGCCGTGGACCCGCTTCGATTCCAAGGGCCGTTGGTTCC
CGGGCCACTGGCCGGGGACTGGAAGGCAGCGGGCATCTGGAAGACGCCAGCCGCTACGTTGGCCGG
CGACCCGGGGTGTCCGGCCCTGAGCGCGCCGCTTTATTCGGGAGCTGGAGGAAGCACTGTGTCTAAC
CTACCTCCGCCAGTCAAAAAGATCACCCAGGAAGACGTCAAAGTATGTTATATTTGCTGGAGGAGGTG
TGTACTGCCTTCCTGATCCAGCTTCTCCACCTGTCTGGGAGAGCGTTACCTATGGGATGGTCTGCAG
AGAGAGAGGGACCTGAACACTGCAGCTCGCATCGGCCAGTCCCTGGTGAACAGAACAGTGTGTTGATG
GAGGAGAACAGCAAGCTGGAAGCCCTGCTGGGCTCAGCCAAGGAGGAGATTTTATACCTCAGACACCAG
GTGAACCTGCGGGATGAGCTCCTCCAGCTCTACTCAGATTCTGATGAGGAGGATGAGGATGAAGAAGAG
GAGGAGGAAGAAAAGGAGGCAGAAGAGGAACAGGAAGAAGAAGCAGAGGAAGACCTGCAGTGTGCT
CATCCCTGTGATGCCCTAAGCTGATTTTCGAGGAGGCATTGCTGCACCAGCACCCTGCCACAGCTG
GAAGCCTTGACAGGAGAAGCTGAGGCTGCTGGAGGAGGAGAATCATCAGCTGAGAGAAGAGGCCTCTAA
CTCGACACTCTTGAGGATGAGGAACAGATGCTCATTCTGGAGTGTGTGGAGCAGTTTTCGGAGGCCAGC
CAACAGATGGCTGAGCTGTCCGAGGTGCTGGTGTCTCAGGCTGGAAAATGAAACGGCAGCAGCAGGAG
GTCGCTCGGCTGCAGGCCAGGTGCTGAAGCTGCAGCAGCGCTGCCGGATGGAGACTTTCCTGGTTTC
CAGGAGACGCTGGCTGAGGAGCTCAGAAGCTCTTAAGGAGGATGATCTCAGACCCTGTGTATTTTATG
GAGAGGAATTATGAGATGCCAGAGGGACACATCCAGCCTAAGGTATGATTTTCGCTACAGTGAGGAT
CGAGAGCAGGTGCGGGGTTTGAGGCTGAGGAAGGTTGATGCTGGCAGCGGATATCATGCGGGGGAA
GATTTACGCCTGCGGAGGAGTTCGTGCCCCAGGAGGAGCTGGGGCTGCCAAGAAGGTGCCGGCTGAG
GAAGGGGTGATGGAAGAGGCAGAGCTGGTGTGAGAGGAGACCGAGGGCTGGGAGGAGGTGGAAGTGGAG
CTGGATGAGGCAACGCGGATGAACGTGGTACATCAGCCCTGGAGGCCAGCGGCTTGGGCCCTTACAC
CTGGACATGAATTATGTCCTCCAGCAGCTGGCCAAGTGGCAAGATGCCATTACAGGCGGCAGCTGAGG
TGAAGATGCTCCAGAAAGGTGAGTGCCCCACGGGGCCCTCCCTGCCGCCAGCCGACAAGCTGCAGA
TCGTCGTGCCGATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001079870

Insert Size: 1809 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001079870.1](#)

RefSeq Size: 3893 bp

RefSeq ORF: 1809 bp

Locus ID: 9001

UniProt ID: [P54257](#)

Cytogenetics: 17q21.2

Protein Pathways: Huntington's disease

MW: 67.6 kDa

Gene Summary: Huntington's disease (HD), a neurodegenerative disorder characterized by loss of striatal neurons, is caused by an expansion of a polyglutamine tract in the HD protein huntingtin. This gene encodes a protein that interacts with huntingtin, with two cytoskeletal proteins (dynactin and pericentriolar autoantigen protein 1), and with a hepatocyte growth factor-regulated tyrosine kinase substrate. The interactions with cytoskeletal proteins and a kinase substrate suggest a role for this protein in vesicular trafficking or organelle transport. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (3) uses an alternate in-frame splice site the 5' coding region, and lacks an alternate in-frame exon in the 3' coding region, compared to variant 2. The encoded isoform (3) is shorter than isoform 2. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The extent of this transcript is supported by transcript alignments.